



I'm here for all your outdoor journeys.

Just like my brother, I'm always there when you want to cover some real distance, giving you the best possible support. With my outdoor outfit, I feel most at home in the fresh air and I'm the perfect companion for walks in nature. My large wheels roll easily over almost any surface, and my drum and parking brakes give you extra safety along the way.



Optional anti-tip mechanism
protects you
from tipping back.

Frame colours Size 1



Frame colours Size 2



Frame colours Size 3+4





data	new			
	Size 1	Size 2	Size 3	Size 4
Handle height	34 - 58 cm	51 - 75 cm	68 - 94 cm	87 - 113 cm
Height of forearm supports	46 - 64 cm	63 - 87 cm	79 - 103 cm	98 - 124 cm
handle spacing	22 - 36 cm	30 - 44 cm	38 - 52 cm	43 - 57 cm
Seat height	33 cm	45 cm	60 cm	70 cm
Overall width	56 cm	66 cm	73 cm	80 - cm
Total length	66 cm	75 cm	82 cm	92 cm
Frame height (folded)	25 cm	32 cm	34 cm	37 cm
Transport dimensions (L x W x H)	66 x 56 x 25 cm	75 x 66 x 32 cm	82 x 73 x 34 cm	92 x 80 x 37 cm
Max. load	30 kg	50 kg	70 kg	100 kg
Weight	7,1 kg	8,3 kg	9,8 kg	10,2 kg
HMV-No.	10.46.02.3052			

Did you know that I am one of the few walking trainers which has an outdoor approval? With larger and wider wheels and a drum and parking brake, I am well-equipped for big adventures outdoors. Brilliant!



You can find more information on accessories in the price list.

