

GMFCS
Level III-V

kg
max. 65




(-)5°-25°

(-)5°-20°

30-42 cm

46-58 cm


28-38 cm



I am your strong therapy chair.

While you actively participate in everyday life, I support you in every way. My super stable frame is designed for active kids and teens like you. My seat unit can easily be tilted forwards and backwards with a lever. Even to the point where you can sit back and relax. For a secure base, I provide up to 3 weights on my frame. So you can always trust me – together we are strong.







Headrest
adjustable in three dimensions,
with reinforced holder,
upholstered incl. cover.



**Headrest with
angle-adjustable
side guide**
is flexible, practical
and individually
adaptable.

Gas pressure pistons
for height and tilt in space adjustment
with locking mechanisms to reduce
movement in the base.





Data	mädita hd.
Seat depth	30 - 42 cm
Seat width	28 - 38 cm
Back height with short back	42 - 52 cm
Back height with long back	53 - 63 cm
Seat height from floor	47,5 - 58 cm
Lower leg length	30 - 48 cm
Total width	53 cm
Total length	77 cm
Tilt-in-space	(-)5° - 20°
Back angle	(-)5° - 25°
Max. load	65 kg
Weight	29 kg
HMV - Nr.	26.11.05.0021

**Super stable with a small footprint:
four weights of 8 kg each can be
attached to my base frame for stabilisation.**



You can find more information on fabrics
on pages 176 - 179.

You can find more information on accessories & harnesses
in the price List.



Scan this QR code
and watch videos about
the product directly on our
homepage.

