

GMFCS
Levels I-V

kg
max. 120



cm
42-98

cm
18-23

I am your tried-and-tested tricycle with E-drive.

I am here to help when you start to get tired. Whether on short or long trips into town or the in countryside. Because I have a silent direct drive, you can enjoy your days out without being disturbed. In addition, my optional start-up assistance supports you, so that you get the tour started. For experienced drivers, my drive is also available in an "advanced version". This gives your a more powerful drive. Even my "beginner version" can be controlled in the Eco, standard or power mode using the display.







DirectPower E-Bike drive
is designed brushless and without gears, minimizes maintenance and wear and tear and eliminates motor noise.



LCD-Display

shows you information on the support level, speed and battery status – and is easy to read, even in bright / sunny surroundings.



Button or twist handle operation

allow you to set the mode, and the start-up assistance if necessary, including during travel.



Frame colours 16"

- Window grey / Pure white
- Telemagenta (pink) / Pure white
- Turquoise blue / Pure white

Frame colours 20" & 24" / 20"

- Window grey / Pure white
- Pastel blue / Pure white
- Kiwi sorbet green / Pure white

Frame colours 24", 26" & 26" (XL)

- Window grey / Pure white
- Ivory / Pure white
- Black / Pure white



No matter where the journey takes you, you can always rely on me.



The battery can also be removed for charging. The charge level is shown on the display.



data	16"	20"	24 / 20"	24"	26"	26" (XL)
Pedal to top edge of saddle*	42 – 62 cm	49 – 71 cm	53 – 75 cm	55 – 77 cm	65 – 87 cm	72 – 94 cm
Pedal to top edge of saddle**	45 – 60 cm	52 – 67 cm	56 – 71 cm	58 – 81 cm	68 – 91 cm	75 – 98 cm
Saddle support tube to front tuber	39 cm	45 cm	54 cm	52 cm	58 cm	
Wheelbase	80 – cm	91 cm	99 cm	104 cm	112 cm	
Access point height	18 cm	20 cm		22 cm	23 cm	
Overall length	120 cm	140 cm	153 cm	165 cm	178 cm	
Overall width	65 cm	69 cm	75 cm			
Max. load	60 kg	80 kg	90 kg	100 kg	120 kg	
Weight	24,5 kg	27,kg	28,5kg	29,5kg	31,5 kg	31,7 kg
HMV-No.	22.51.04.0005					

* Measurement with standard seat post / ** Measurement with T-seat post



Further information on accessories and harnesses
can be found in the price list.

