



activity chair.

Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The dynamic therapy chair.

- Seat and back unit incl. pad
- Back angle adjustment via gas pressure spring from (-)5° to 20°
- Back height and seat depth adjustment
- Incl. 2-point pelvic harness
- Tilt-in-space adjustment via gas pressure spring
 - for the standard base (-) 15° to 15°
 - for the Hi/Lo base (-) 15° to 25°
- Frame colour: Effect black



HMV-No.: 26.11.05.0046
(Standard base)

HMV-No.: 26.11.05.0047
(Hi/Lo base)

Basic model



activity chair.

Seat unit (standard)

Tilt-in-space adjustment via gas pressure spring (-) 15° to 20°

Price / piece

Size 0 (72 00 200)	1.028 €
Size 1 (72 01 200)	1.050 €
Size 2 (72 02 200)	1.200 €



activity chair.

Seat unit (dynamic)

Back angle adjustment via gas pressure spring from (-)5° to 20°; a gas pressure spring enables self-generated movement to be accommodated - the spring option offers the user 10° of movement.

Price / piece

Size 0 (72 00 300)	1.173 €
Size 1 (72 01 300)	1.200 €
Size 2 (72 02 300)	1.387 €

Dimensions				
		Size 0	Size 1	Size 2
Seat depth		20 - 30 cm	28 - 41 cm	38 - 51 cm
Seat width		30 cm	36 cm	43 cm
Seat width (with hip support pads)		18 - 23 cm	22 - 29 cm	28 - 36 cm
Chest width		17 - 29 cm	17 - 30 cm	24 - 36 cm
Back height		32 - 39 cm	39 - 48 cm	48 - 61 cm
Height headrest - seat		37 - 53 cm	44 - 62 cm	50 - 75 cm
Height armrest - seat		13 - 19 cm	18 - 27 cm	18 - 27 cm
Lower leg length when used with footrest		18 - 30 cm	20 - 38 cm	29 - 54 cm
Seat height above floor	with short / long feet	24 - 32 cm / 47 - 55 cm	32 - 42 cm / 48 - 58 cm	42 - 52 cm / 51 - 61 cm
	with short / long feet and dual castors	34 - 42 cm / 47 - 55 cm	41 - 51 cm / 48 - 58 cm	42 - 52 cm / 51 - 61 cm
	with Hi/Lo base	25,5 - 63,5 cm	30,5 - 58,5 cm	40,5 - 61 cm
Frame width	with short / long feet	53 / 58 cm	59 / 62 cm	68 / 68 cm
	with short / long feet and dual castors	57 / 60 cm	62 / 65 cm	71 / 71 cm
	with Hi/Lo base	66 cm	68,5 cm	73,5 cm
Tilt-in-space	with standard base / with Hi/Lo base	(-) 15° - 15° / (-) 15° - 25°		
Back angle		(-) 5° - 20°		
Total length/frame depth		54-73 cm	57-73 cm	77-90 cm
Total width		53-60 cm	59-65 cm	68-70,5 cm
Total height incl. headrest		61-108 cm	76-120 cm	92-136 cm
Total height without headrest		56-94 cm	71-106 cm	90-122 cm
Recommended body height		81-122 cm	102-142 cm	127-188 cm
max. load		34 kg	68 kg	113 kg
Weight of standard base frame		12,2 kg	15,9 kg	24,5 kg
Weight of Hi/Lo frame		25,4 Kg	39 Kg	44,5 kg

The dimensions given apply to both seat units with all bases. If certain dimensions differ between the base frames, these are listed specifically.

Base frame			
Standard base activity chair. incl. height-adjustable, non-slip feet; seat angle adjustment by gas pressure spring from (-) 15° to 15°		Dynamic base activity chair. incl. height-adjustable, non-slip feet; seat angle adjustment by gas pressure spring from (-) 15° to 15°; a gas pressure spring enables self-generated movement to be accommodated - the spring option offers the user room to move.	
			
Price / piece		Price / piece	
Size 0 (72 00 400)	730 €	Size 0 (72 00 500)	842 €
Size 1 (72 01 400)	743 €	Size 1 (72 01 500)	871 €
Size 2 (72 02 400)	871 €	Size 2 (72 02 500)	1.016 €
		Hi/Lo base activity chair. incl. 4 double castors (100 mm, braked, 1x with directional lock for size 0) & 1 directional assist wheel (for sizes 1 & 2); height adjustment (by gas pressure spring for size 0 and hydraulically for sizes 1+2); seat angle adjustment by gas pressure spring from (-) 25° to 15°	
			
		Price / piece	
		Size 0 (72 00 600)	2.270 €
		Size 1 (72 01 700)	3.187 €
		Size 2 (72 02 700)	3.381 €

Pad colour			
Blue		Red	
			
		Price	
		for size 0 (72 05 201)	incl.
		for size 1 (72 00 201)	incl.
		for size 2 (72 02 201)	incl.
		Price	
		for size 0 (72 05 301)	incl.
		for size 1 (72 00 301)	incl.
		for size 2 (72 02 301)	incl.

Feet			
	Short feet – standard configuration, if nothing else stated		Long feet – for all sizes
	Price / set		Price / set
	for size 0 (72 04 004)	incl.	72 01 004
	for size 1 (72 00 004)	incl.	incl.
	for size 2 (72 05 004)	incl.	

Feet			
	4 Dual castors (75 mm) short version, with brakes, 1 x with direction lock		4 Dual castors (75 mm) long version, with brakes, 1 x with direction lock, for all sizes
	Upcharge / set		Upcharge / set
	for size 0 (72 06 004)	393 €	72 03 004 393 €
	for size 1+2 (72 02 004)	393 €	

Accessories			
Headrest (Standard) depth, angle and height adjustable 	Headrest (Shell-shaped) depth, angle and height adjustable 		Headrest (Lateral guide) depth, angle and height adjustable, wing pads can be adjusted from almost flat to 90° 
	Price / piece		Price / piece
	72 00 010	352 €	72 00 011 378 € 72 00 018 485 €

Accessories					
Lateral pads – height and width adjustable 			Lateral pads – height and width adjustable, with chest strap 		Push handles – recommended for all activity chair . with dual castors and models with Hi/Lo base 
	Depth	Height		Depth	Height
for size 0-2	13 cm	10 cm	for size 0+1	13 cm	10 cm
for size 1-2	16,5 cm	13 cm	for size 2	16,5 cm	13 cm
Price / pair			Price / set		Price / pair
for size 0-2 (72 01 005)			for size 0 (72 00 006)	393 €	72 00 012 250 €
for size 1-2 (72 02 005)			for size 1 (72 01 006)	393 €	
			for size 2 (72 02 006)	421 €	

Accessories					
Armrests – height adjustable 			Hip support pads – padded and angle adjustable (only for use together with the armrests) 		Therapy tray – only for use together with the armrests 
	Width	Length		Height	Length
for size 0	6 cm	22 cm	for size 0	8 cm	14 cm
for size 1+2	7,5 cm	30,5 cm	for size 1+2	10 cm	23,5 cm
Price / pair			Price / pair		Price / piece
for size 0 (72 05 007)			for size 0 (72 05 008)	393 €	for size 0 (72 00 009) 462 €
for size 1+2 (72 00 007)			for size 1+2 (72 00 008)	393 €	for size 1 (72 01 009) 462 €
					for size 2 (72 02 009) 462 €

Accessories

Backrest extension		Spinal support and tube edge		Abduction wedge depth adjustable	

Accessories

Knee lateral guide 	Knee support pad with key lock, depth adjustable 
Price / pair 72 00 014 309 €	Price / piece for size 0 (72 00015) 348 € for size 1 (72 01015) 348 € for size 2 (72 02015) 378 €

Harness

Chest strap 	Chest-shoulder strap  normal dünn	Thigh straps 	Groin harness 
Price / piece for size 0 (72 00020) 153 € for size 1 (72 01020) 153 € for size 2 (72 02020) 153 €	Price / piece normal for size 0 (72 00021) 227 € normal for size 1 (72 01021) 227 € normal for size 2 (72 02021) 227 € thin for size 1 (72 04021) 227 € thin for size 2 (72 03021) 227 €	Price / pair for size 0 (72 00022) 153 € for size 1 (72 01022) 153 € for size 2 (72 02022) 153 €	Price / piece for size 0 (72 00023) 268 € for size 1 (72 01023) 268 € for size 2 (72 02023) 268 €

Accessories

Footplate – height adjustable, can be folded up and swung away for transfer, recommended for models with Hi/Lo base. 	Ankle positioner – with socket lock, only for use together with the footplate 
Price / piece for size 0 (72 00029) 702 € for size 1 (72 01029) 702 € for size 2 (72 02029) 702 €	Price / pair for size 0+1 (72 01026) 153 € for size 2 (72 02026) 153 €

Accessories

Sandals with foot straps only for use together with the footplate  <table border="1"> <thead> <tr> <th></th><th>Price / pair</th></tr> </thead> <tbody> <tr> <td>for size 0 (72 00027)</td><td>153 €</td></tr> <tr> <td>for size 1 (72 01027)</td><td>153 €</td></tr> <tr> <td>for size 2 (72 02027)</td><td>153 €</td></tr> </tbody> </table>		Price / pair	for size 0 (72 00027)	153 €	for size 1 (72 01027)	153 €	for size 2 (72 02027)	153 €	Wedges – only for use together with the footplate and sandals  <table border="1"> <thead> <tr> <th></th><th>Price / pair</th></tr> </thead> <tbody> <tr> <td>for size 0 (72 00028)</td><td>122 €</td></tr> <tr> <td>for size 1 (72 01028)</td><td>122 €</td></tr> <tr> <td>for size 2 (72 02028)</td><td>122 €</td></tr> </tbody> </table>		Price / pair	for size 0 (72 00028)	122 €	for size 1 (72 01028)	122 €	for size 2 (72 02028)	122 €
	Price / pair																
for size 0 (72 00027)	153 €																
for size 1 (72 01027)	153 €																
for size 2 (72 02027)	153 €																
	Price / pair																
for size 0 (72 00028)	122 €																
for size 1 (72 01028)	122 €																
for size 2 (72 02028)	122 €																