



**Invoice address:**

**Dispatch address:**

See invoice address

**Order**

**Cost estimate**

Customer No.: \_\_\_\_\_

Orderer: \_\_\_\_\_

Date: \_\_\_\_\_

Signature: \_\_\_\_\_

**Commission:**

**The versatile bathing and shower cradle.**


- Black, powder-coated frame with removable, washable cover
- Incl. back angle, knee angle and head angle adjustment
- Incl. chest strap

**HMV-No.: 04.40.04.0006**



**Dimensions**

|   |                                                              | Size 1     | Size 2 | Size 3 | Size 4 |
|---|--------------------------------------------------------------|------------|--------|--------|--------|
| A | Seat width                                                   | 30 cm      | 35 cm  |        | 40 cm  |
| B | Seat depth                                                   | 30 cm      | 34 cm  | 38 cm  | 44 cm  |
| C | Back height                                                  | 45 cm      | 61 cm  | 71 cm  | 82 cm  |
| D | Standing space length                                        | 69 cm      | 74 cm  | 88 cm  | 93 cm  |
| E | Lower leg length                                             | 19 cm      | 20 cm  | 28 cm  |        |
| F | Total width                                                  | 36 cm      | 41 cm  |        | 46 cm  |
|   | Back angle                                                   | 90° - 180° |        |        |        |
|   | Knee angle                                                   | 0° - 90°   |        |        |        |
|   | Headrest angle                                               | 0° - 90°   |        |        |        |
|   | Seat height                                                  | 0* / 5 cm  |        |        |        |
|   | Max. load                                                    | 35 kg      | 45 kg  | 55 kg  | 75 kg  |
|   | Weight                                                       | 4 kg       | 4,1kg  | 4,3kg  | 4,8kg  |
|   | * in combination with the cover with pelvic mould depression |            |        |        |        |



### Basic model

| basti. size 1          | basti. size 2            | basti. size 3            | basti. size 4<br>back angle adjustment via<br>stainless steel ratchet joint |
|------------------------|--------------------------|--------------------------|-----------------------------------------------------------------------------|
| Price / piece          | Price / piece            | Price / piece            | Price / piece                                                               |
| <b>56 01 000</b> 989 € | <b>56 02 000</b> 1.051 € | <b>56 03 000</b> 1.107 € | <b>56 04 000</b> 1.255 €                                                    |

### Accessories

#### Cover with hollow for pelvis and adjustable tube edge



#### Head guide pad – adjustable with velcro



| Upcharge / piece                |       | Price / piece                   |       |
|---------------------------------|-------|---------------------------------|-------|
| for size 1 ( <b>56 01 001</b> ) | 160 € | for size 1 ( <b>56 01 002</b> ) | 261 € |
| for size 2 ( <b>56 02 001</b> ) | 160 € | for size 2 ( <b>56 02 002</b> ) | 261 € |
| for size 3 ( <b>56 03 001</b> ) | 191 € | for size 3 ( <b>56 03 002</b> ) | 261 € |
| for size 4 ( <b>56 04 001</b> ) | 191 € | for size 4 ( <b>56 04 002</b> ) | 261 € |

### Accessories

#### Additional retainer strap – adjustable with velcro for the chest or pelvis



#### Shower cradle subframe basti. – mobile with 4 castors with brakes



| Price / piece                       |       | Price / piece    |       |
|-------------------------------------|-------|------------------|-------|
| for size 1 ( <b>56 01 003</b> )     | 126 € | <b>56 00 200</b> | 868 € |
| for size 2 + 3 ( <b>56 03 003</b> ) | 126 € |                  |       |
| for size 4 ( <b>56 04 003</b> )     | 126 € |                  |       |

### Dimensions - subframe

|          |                                                              |             |
|----------|--------------------------------------------------------------|-------------|
| <b>A</b> | Length                                                       | 92 cm       |
| <b>B</b> | Width                                                        | 57 cm       |
| <b>C</b> | Seat height with subframe                                    | 61* / 66 cm |
|          | * in combination with the cover with pelvic mould depression |             |

