



120

boris.

Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:


The comfortable toilet and shower chair.

- Anatomically-shaped seat pan made of PE, colour steel blue
- Subframe made of stainless steel
- Seat angle adjustment via gas pressure spring:
130 N (size 0), 200 N (sizes 1-3) and 2x 130 N (size 4)
- Mechanical height adjustment
- 4 dual castors (100 mm) with locking devices

HMV-No.: 33.40.04.1014



Presettings



Presettings on the toilet and shower chair at Schuchmann.
Please provide us with the following body measurements:

I Back height _____

II Pelvic width
(Seat width rear) _____

III Seat depth _____


IV Lower leg length _____

52 00 080

Piece

incl.

Dimensions



		Size 0	Size 1	Size 2	Size 3	Size 4
A	Seat width (rear)	26 cm	31 cm	36 cm	41 cm	46 cm
B	Seat depth	31 cm	36 cm	41 cm	46 cm	48 cm
C	Seat shell width (front)	31 cm	36 cm	41 cm	46 cm	51 cm
D	Back height	49 cm	51 cm	56 cm	61 cm	65 cm
E	Tilt-in-space	(-) 5° - 30°				
F	Can be positioned over a toilet	47 - 57 cm	45 - 55 cm	42 - 52 cm		
G	Total width	56 cm	61 cm	63 cm	71,5 cm	74 cm
H	Total length	65 cm			70 cm	
I	Max. lower leg length	38 cm		48 cm		
J	Front edge of hygiene opening to back edge of subframe	47 cm	49 cm	51 cm	54 cm	57 cm
K	Front opening width	10 cm				
L	Rear opening width	12 cm		13,5 cm		
M*	Rear opening width	-	-	18,5 cm		
	Max. load	60 kg	80 kg	100 kg		150 kg
	Weight	28 kg	30 kg	33 kg	35 kg	38 kg

* only applies to sizes 2b, 3b and 4b

Basic model

boris. size 0	boris. size 1	boris. size 2	boris. size 3	boris. size 4
Price / piece	Price / piece	Price / piece	Price / piece	Price / piece
52 00 100 2.750 €	52 01 100 2.781 €	Size 2a (52 02 100) 2.781 € Size 2b (52 02 200) 2.859 €	Size 3a (52 03 100) 2.859 € Size 3b (52 03 200) 2.941 €	Size 4a (52 04 100) 3.204 € Size 4b (52 04 200) 3.292 €

Accessories

PU-Headrest – blue, adjustable in height and depth with stainless steel bracket		Tray – with padded arm rest segment and holder	
Price / piece		Price / piece	
52 00 001 281 €		for size 0 (52 00 012) 452 € for size 1 (52 01 002) 452 € for size 2 (52 02 002) 452 € for size 3 (52 03 002) 452 € for size 4 (52 04 002) 452 €	

Accessories

PU insert pad – two-piece for opening form a


Price / piece

for size 0 (52 00 101)	397 €
for size 1 (52 01 101)	397 €
for size 2a (52 02 101)	397 €
for size 3a (52 03 101)	397 €
for size 4a (52 04 101)	397 €

PU insert pad – two-piece for opening form b


Price / piece


for size 2b (52 02 201)	446 €
for size 3b (52 03 201)	446 €
for size 4b (52 04 201)	446 €

Accessories

Toilet bucket – incl. lid and bracket as assembly set

Abductor wedge – with integrated splash guard (Note: a subsequent assembly at size 4 has to take place in our factory)


	Width	Height from seat
for size 0 + 1	7 cm	14 cm
for size 2 - 4	8 cm	17,5 cm

Footplate – height, knee and foot angle adjustment, also folding


	Adjustable lower leg length
for size 0 + 1	30 – 38 cm
for size 2 - 4	35 – 48 cm

PU-Calf plate – height adjustable (only for use in combination with the footplate)


Height	15 cm
Width	38 cm


Piece / set	Price / piece	Price / piece	Price / piece
for size 0+1+3 (52 00 102)	203 €	for size 0+1 (52 01 103)	218 €
for size 2+4 (52 02 102)	203 €	for size 2+4 (52 02 103)	218 €
		for size 3 (52 03 103)	218 €
		for size 0+1 (52 01 004)	390 €
		for size 2-4 (52 02 004)	390 €
		52 01 009	246 €

Accessories

Positioning jacket – made of spacer fabric incl. strap holder for shoulder straps, colour: blue


Pelvic harness – made of spacer fabric, colour: blue


Strap holder

Strap holder – for shoulder straps when using other harnesses apart from the adjacent positioning jacket


Piece / set	Price / piece	Art.-No.	Usage	max. shoulder width	Price / piece
Size 1 (52 01 010)	301 €	52 01 005	with headrest	30 cm	127 €
Size 2 (52 02 010)	301 €	52 01 008	without headrest	30 cm	127 €
Size 3 (52 03 010)	301 €				

Accessories

Ankle cuffs – with velcro fastener and idler pulley; colour: black


	Ankle circumference	Strap width
Size 1	12 – 19 cm	2,5 cm
Size 2	15 – 23 cm	2,5 cm
	Price / pair	
Size 1 (52 01 018)	236 €	
Size 2 (52 02 018)	236 €	