



dynamic pacer.

Order Cost estimate

Invoice address:

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Dispatch address:

See invoice address

Commission:

The dynamic gait trainer.

- Modular set with various options to combine different bases and upper frames (except size 4)
- Upper frame to accommodate accessories, continuous height adjustment via gas spring support (for size 1-4)



HMV-No.: 10.46.02.3005
(Size 0)

HMV-No.: 10.46.02.3044
(with standard base)

HMV-No.: 10.46.02.3046
(with utility base)



Frame colours				
Colour selection for size 0	Colour selection for size 1	Colour selection for size 2	Colour selection for size 3	Colour selection for size 4
Piece	Piece	Piece	Piece	Piece
Lime (92 00 153) incl. White (92 00 155) incl.	Lime (92 00 153) incl. White (92 00 155) incl.	Red (92 00 150) incl. White (92 00 155) incl.	Blue (92 00 152) incl. Black (92 00 156) incl.	Black (92 00 156) incl.

Basic model pacer mini. (Size 0)



pacer mini. (size 0)

one-piece base frame with gripping frame, 4 swivel castors with grinding brake, recoil stop, parking brake and directional fixed spring

	Price / piece
Size 0 without a distance meter (3300 200)	1.449 €
Size 0 without a distance meter (3300 201)	1.775 €

Basic model - Upper frames (sizes 1-4)

Standard upper frames



Dynamic Upper frames – Additionally, with integrated and separately adjustable vertical and horizontal dynamic for a physiological natural gait.



Price / piece

Price / piece

Size 1 **(33 01 200)**
Size 2 **(33 02 200)**
Size 3 **(33 03 200)**
Size 4 **(33 04 200)**

1.165 €
1.633 €
1.947 €
2.718 €

Size 1 **(33 01 202)**
Size 2 **(33 02 202)**
Size 3 **(33 03 202)**
Size 4 **(33 04 202)**

2.016 €
2.718 €
3.059 €
3.857 €

Basic model - Bases (Size 1-4)

Standard base

incl. four 4" (100 mm) castors (size 1) or four 5½" (140 mm) castors (sizes 2–4) with parking brake, friction brake, backwards stop and direction lock



Utility base

incl. two 8" (200 mm) front wheels with direction lock and two 11 ½" (290 mm) rear wheels with parking brake, friction brake and backwards stop, for indoor and outdoor use (for outdoor use the function and parking handbrake must be ordered)



Accessories

Functional and parking handbrake

can only be used on the Utility base frame and in conjunction with the upholstered armrests or the underarm supports (can only be retrofitted with a frame hole)



Price / piece

Price / piece



Piece / pair

Size 1 without odometer **(33 01 204)** 752 €
Size 1 with odometer **(33 01 205)** 1.052 €
Size 2 without odometer **(33 02 204)** 996 €
Size 2 with odometer **(33 02 205)** 1.328 €
Size 3 without odometer **(33 03 204)** 1.333 €
Size 3 with odometer **(33 03 205)** 1.671 €
Size 4 without odometer **(33 04 204)** 2.357 €
Size 4 with odometer **(33 04 205)** 2.639 €

Size 2 without odometer **(33 02 206)** 1.494 €
Size 2 with odometer **(33 02 207)** 1.826 €
Size 3 without odometer **(33 03 206)** 1.831 €
Size 3 with odometer **(33 03 207)** 2.168 €

for size 2 **(33 02 011)** 556 €
for size 3 **(33 03 011)** 611 €

Dimensions



			Size 0	Size 1	Size 2	Size 3	Size 4
A	Elbow height		39 – 52 cm	47 – 70 cm	61 – 89 cm	81 – 119 cm	86 – 124 cm
B	Total width	Standard base	52 cm	58 cm	66 cm	71 cm	80 cm
		Utility base	–	–	76 cm	81 cm	–
C	Total length	Standard base	57 cm	66 cm	76 cm	91 cm	101 cm
		Utility base	–	–	91 cm	103 cm	–
D	Frame height	Standard-/ Utility base	32 cm	41 – 53 cm	55 – 69 cm	72 – 98 cm	79 – 104 cm
E	Base height (without upper frame)	Standard-/ Utility base	–	28 cm	37 cm	41 cm	47 cm
	Weight	Standard base	4,6 kg	5,2 kg	6,8 kg	7 kg	10,2 kg
		Utility base	–	–	8,4 kg	8,7 kg	–
		Standard upper frames	–	2,7 kg	3,4 kg	4,4 kg	4,9 kg
		Dynamic Upper frames	–	3,4 kg	5,2 kg	6,1 kg	6,8 kg
	Dynamic upper frame movement	Vertical	–	3,5 cm	8 cm	8 cm	8 cm
		Horizontal	–	4 cm	5 cm	5 cm	5 cm
	Max. load		23 kg	34 kg	68 kg	91 kg	113 kg

Accessories

Multi-position saddle (MPS)

The multi-position saddle (MPS) is the ultimate positioning accessory for gait. It quickly fine tunes pelvis positioning for anterior or posterior use.

Adjustable in five directions to accommodate each client's posture, the MPS easily attaches to both the dynamic and standard upper frames.



Seat pad Dimensions	Size 1+2	Size 3
Front width (a)	7,5 cm	11,5 cm
Middle width (b)	4,5 cm	6,5 cm
Rear width (c)	9 cm	13,5 cm
Length (d)	24 cm	27,5 cm

Saddle depth – adjusts to optimally position the pelvis in relation to the upper body to encourage forward motion during ambulation.

Adjustment range: 13 cm



Saddle angle – is adjustable to position the pelvis in anterior or posterior tilt for the most effective gait pattern.

Forward: 7.5° and 15°
Back: 7.5° and 15°



Saddle height – adjusts to finetune the position of the pelvis during gait.

Adjustment range: 22 cm



Hip corral – depth and height adjustments provide cushioning and support for the pelvis at the level of ASIS.

Adjustment range depth: 8 cm
Adjustment range height: 8 cm



Price / piece

Size 1 for size 1	(33 01012)
Size 2 for size 2-4	(33 02012)
Size 3 for size 2-4	(33 03012)

1.685 €
1.743 €
1.893 €

Accessories

MPS Padding cover to reduce frictional resistance


Price / piece

for MPS Size 1+2	(33 01024)
for MPS Size 3	(33 03024)

136 €
136 €

Accessories

Pelvic support – incl. bracket and a pair of handholds




	Size 0	Size 1	Size 2
Key user dimension*	15 – 20 cm	15 – 25 cm	20 – 36 cm
a	17 cm	22 cm	30 cm
b	30 cm	41 cm	51 cm
c	18 cm	23 cm	30 cm

* half of thigh circumference

Piece / set

Size 0 (33 00015)	762 €
Size 1 (33 01015)	830 €
Size 2 (33 02015)	902 €

Hip positioner – adjustable, incl. bracket and a pair of handholds



	Size 1	Size 2
Stride length (a)	6 cm	9 cm
Width (b)	3 cm	4,5 cm
Length (c)	23 cm	28 cm

Piece / set

Size 1 (33 01003)	665 €
Size 2 (33 02003)	739 €

Pad – for hip positioner



Price / piece

for size 1 (33 01004)	241 €
for size 2 (33 02004)	241 €

Accessories

Chest prompt – incl. bracket

	Range of circumference	Height
Size 1	36 – 71 cm	11 cm
Size 2	56 – 102 cm	17 cm
Size 3	71 – 127 cm	20 cm
		Piece / set
Size 1 (33 01002)		885 €
Size 2 (33 02002)		958 €
Size 3 (33 03002)		958 €



Chest unit – with back strap and holder

	Height	Range of circumference
for size 0	10 cm	36 – 71 cm
for size 1	10 cm	36 – 71 cm
for size 2	12,5 cm	56 – 102 cm
for size 3	15 cm	71 – 127 cm
		Piece / set
for size 0 (33 00020)		553 €
for size 1 (33 01020)		553 €
for size 2 (33 02020)		603 €
for size 3 (33 03020)		632 €



Accessories

Grip rings – incl. bracket



Piece / pair

33 00 008

351 €

Handles – with bracket



	Height	Width
Size 1	12,5cm	9cm
Size 2	18cm	9cm
Piece / pair		
Size 1 (33 01 019)	362 €	
Size 2 (33 02 019)	384 €	

Accessories

Arm prompts – incl. handles and bracket



Padded surface

Size 1

11x18 cm

Adjustment range (elbow to fist)

Size 1

13x22 cm

20 – 32 cm

27 – 41 cm

Piece / pair

Size 1 (**33 01 001**)

758 €

Size 2 (**33 02 001**)

758 €

Arm platforms – incl. bracket (recommended for sizes 2–4)


without handles

Padded surface

with handles

34 x 14 cm

34 x 14 cm

Piece / pair

without handles (**33 00 010**)

711 €

with handles (**33 01 010**)

842 €

Accessories

Additional handholds for all sizes



Piece / pair

33 00 002

215 €

Therapy tray – angle adjustable and transparent cover



Inner dimensions

22x28 cm

Max. load

5 kg

for size 0 (**33 00 001**) 724 €
for size 1–4 (**33 00 013**) 724 €

Thigh prompts – incl. bracket



Length

Size 1

13 – 20 cm

Size 2

13 – 28 cm

Max. leg circumference

38 cm

51 cm

Piece / pair

Size 1 (**33 01 005**)


661 €

Size 2 (**33 02 005**)

785 €

Accessories

Ankle prompts for standard base



Piece / set

for size 0 (**33 00 016**)

357 €

for size 1 (**33 01 016**)

357 €

for size 2 (**33 02 016**)

357 €

for size 3+4 (**33 03 016**)

403 €

Ankle prompts for utility base



Piece / set

for size 2 (**33 02 017**) 357 €

for size 3 (**33 03 017**) 403 €

Attendant guide bar



Price / piece

for size 0 (**33 00 009**) 290 €

for size 1 (**33 01 022**) 290 €

for size 2+3 (**33 02 022**) 290 €

Accessories tote

35 cm x 40 cm x 23 cm



Price / piece

33 00 014

164 €