



follow-me momo.

Order Cost estimate

Invoice address:

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Dispatch address:

See invoice address

Commission:

The clever parent-child tandem coupling for momo tricycles.

- To suit all **Schuchmann tricycles** with free wheeling and / or switchover to free wheeling
- For all 26" – 29" adult's bicycles with and without rear suspension with quick release wheel fixation
- Wheels with full axle or wheel hub gears can be prepared using optionally available adaptors
- The child's seat, luggage carrier and bags can still be fully used on the leading bicycle
- It is possible to attach and remove the tricycle within seconds
- Incl. special axle nuts



Dimensions

	Size 1	Size 2	Size 3
Required size of the leading bicycle	26" – 29"		
Max. tricycle size	12" (to be used also for 16" + 20")	16/12" + 16" + 20"	24/20" + 24" + 26"
Max. traction weight	40 kg	60 kg	80 kg
Weight	5,1kg		8,7kg

Basic model

follow-me momo. size 1 for 12"		follow-me momo. size 2 for 16/12" + 16" + 20"		follow-me momo. size 3 for 24/20" + 24" + 26"	
Price / piece		Price / piece		Price / piece	
39 00 100		39 01 100		39 02 100	
546 €		546 €		805 €	


Use with third-party bikes and two-wheelers with supports (e.g. momo therapy bike.) is not permitted!

Adapter for rear axle of towing bike

Hub or chain gears and quick release with a max. frame width of 150 mm	Hub or chain gears and quick release with a max. frame width of 150 mm to 175 mm	Hub or chain gears with full axle and nuts	Hub gears and clickbox	Hub or chain gears and thru axle
				
<i>i</i> One of these four adapters is included. Each additional adapter will be charged at 95 € / set.				
Piece / set	Piece / set	Piece / set	Upcharge / Set	Upcharge / Set
39 01 101 incl.	39 02 101 incl.	39 01 102 incl.	39 01 103 72 €	39 01 104 104 €

i When ordering, we need pictures of both sides of the axle!

Dimension query thru axle



Axle width (X): _____ Thread size (Y): M12 x 1,0 mm 1,5 mm 1,75 mm Thread length (Z): _____
 Manufacturer: Syntace DT-Swiss Shimano Sram Maxle _____

Accessories

Suspension bracket – for luggage carrier with „Top-Light“, incl. carabiner hooks and assembly instructions (incl. for size 3)



Preparing the follow-me momo. for installation.

Adjustment of the **follow-me momo.** to the wheel size when ordering a new **momo tricycle.** (at Schuchmann).

Price / piece		Net price (not discountable)	
39 01 109 29 €		39 01 110 42 €	

Note

Adult bike with rear engine

If the adult bike on which the **follow-me momo.** is to be mounted is equipped with a rear engine, then an adapter other than described above is needed. Please contact us in such case and we will work with you to develop an adapter solution.

i If the cable of a rear motor is routed through the axle, compatibility with the **follow-me momo.** is not given.

momo motion. or momo tricycle. with electric drive in combination with the follow-me momo.

Generally we do not recommend to pair a **momo motion.** or **momo tricycle.**, which is equipped with an electric drive, via **follow-me momo.** with an adult bike. This can lead to dangerous situations due to the support that comes from the electric drive.

Please contact us if you still require this combination. We will gladly try to find a solution together with you.