



franka.

47

Invoice address:

Dispatch address:

See invoice address

**Order
Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The universal subframe.

- Tilt adjustment of the seat pan via gas pressure spring 500 N
- Height adjustment via gas pressure spring (800 N for size 0; 1200 N for size 1; 1600 N for size 2)
- 4 castors (75 mm), 2 with locking devices
- incl. mounting head for trapezoidal and third-party adapters and footrest insertion


HMV-No.: 26.99.01.0022



Dimensions

		Size 0	Size 1	Size 2
A*	Height of seat pan support	18 - 55 cm (20 - 57 cm**)		
	Height of seat pan support with electric height adjustment	26 - 55 cm (28 - 57 cm**)		
	Height of seat pan support with hydraulic height adjustment	20 - 56 cm (22 - 58 cm**)		
A	Height of seat pan support with combi subframe	-	25 - 64 cm (27 - 66 cm**)	
B	Total length	58 cm	65 cm	75 cm
	Total length with combi subframe	71 cm	76 cm	86 cm
C	Total width	52 cm	54 cm	64 cm
	Total width with combi subframe	62 cm	66 cm	76 cm
D	Frame height	27 cm		
E	Tilt-in-space	(-) 5° - 35°		
	Max. load	30 kg	50 kg	75 kg
	Weight	8 kg	10 kg	12 kg



* with 75 mm castors ** incl. trapezoidal adapter



Basic model					
franka. size 0		franka. size 1		franka. size 2	
Price / piece		Price / piece		Price / piece	
2100000	1.319 €	2101000	1.440 €	2102000	1.505 €

Frame colours (all available from stock)					
Traffic white		Telegrey 4		Anthracite grey	
					
Piece		Piece		Piece	
9200101	incl.	9200011	incl.	9200017	incl.


Push bar / handle		Accessories			
Push bar – angle adjustable, mounted on the support head 		Push handle – height and angle adjustment (mounted on the base frame) 		Anti tip – adjustable in height 	
Price / piece		Price / piece		Price / piece	
for support head size 0 (2100013) 262 €		for size 0+1 (2100001) 356 €		for stand. & electric height adjustment (2100009) 281 €	
for support head size 1 (2101013) 262 €		for size 2 (2100201) 356 €		for hydr. height adjustment (2100109) 281 €	
				Wheel base extension – pluggable, mounted in the front, length = 10 cm 	
				Price / set	
				2100209 191 €	

Seat angle adjustment			
Trigger lever – mounted directly on the gas pressure spring (standard configuration, if nothing else stated) 		Trigger lever – mounted on the push bar/handle 	
Price / piece		Upcharge / piece	
2100002 incl.		2100014 38 €	
		Trigger lever – with fixing plate for assembly to the back of the seat pan 	
		Upcharge / piece	
		2100015 62 €	
		Safety bolt for Trigger lever, to prevent unintentional release of the gas pressure spring (Not possible with release lever „directly on the gas spring“..) 	
		Price / piece	
		2101014 82 €	


Accessories			
Conversion kit for combi subframe with 12 1/2" wheels, screw-in axles, drum brake, steering wheels 6" and Push handle 		Hydraulic height adjustment retrofitting not possible! (not possible together with the conversion kit for combi subframe or electric height adjustment) 	
Price / set		Upcharge / set	
for size 0 (2100101) 933 €		for size 1 (2101050) 589 €	
for size 1 (2101001) 933 €		for size 2 (2102050) 589 €	
for size 2 (2102001) 933 €			
		Electric height adjustment by a 24-volt electric motor; incl. battery and hand control unit 	
		Upcharge / set	
		2100107 1.034 €	
		Additional gas dampers (2 pcs) to support the lift function (not possible in combination with the hydraulic or electric height adjustment!) 2x 300 N Recommended from 30 kg load 2x 500 N Recommended from 30 kg load	
		Price / set	
		2x 300 N (2100005) 133 €	
		2x 500 N (2100006) 133 €	

Adapter

Support head size 0 – for trapezoidal adapter size 0 and external adapter




Trapezoidal adapter lower part size 0 – for assembly on the subframe



	Length	Front width	Rear width
Size 0	20 cm	21 cm	14,5 cm

Trapezoidal adapter top part size 0 – for assembly on the seat pan



	Length	Front width	Rear width
Size 0	19 cm	16,5 cm	9,8 cm

Price / piece

2102003

incl.

Price / piece

2100016

217 €


Price / piece

2100017


143 €

Adapter

Support head size 1 – for trapezoidal adapter size 1 and external adapter



Trapezoidal adapter lower part size 1 – for assembly on the subframe



	Length	Front width	Rear width
Size 1	26 cm	26 cm	17,5 cm

Trapezoidal adapter top part size 1 – for assembly on the seat pan



	Length	Front width	Rear width
Size 1	25 cm	21,5 cm	12,7 cm

Price / piece

2100003

incl.

Price / piece

2101016

217 €

Price / piece

2101017

143 €

Footplate



Footplate standard – without fixed position, foldable, with non-slip surface and 4 slots each for fixing straps


Footplate with lock – foldable, with non-slip surface and 4 slots each for fixing straps

Footplate with angle adjustment – foldable, with non-slip surface and 4 slots each for fixing straps

*The minimum lower leg length is achieved by repositioning the screws

Calf plate

Calf plate for Footplate covered with black flexi-fabric



	Size 1	Size 2	Size 3
Width	27 cm	32 cm	38 cm
Height	6,5 cm	6,5 cm	6,5 cm

Price / piece

Size 1 (**2101023**)
Size 2 (**2102023**)
Size 3 (**2103023**)

377 €
377 €
377 €

Price / piece

Size 1 (**2101024**)
Size 2 (**2102024**)
Size 3 (**2103024**)

435 €
435 €
435 €

Price / piece

Size 1 (**2101025**)
Size 2 (**2102025**)
Size 3 (**2103025**)

483 €
483 €
483 €

Price / piece

for size 1 (**2101026**)
for size 2 (**2102026**)
for size 3 (**2103026**)

271 €
271 €
271 €

Padded cover for footplate

Padded cover for footplate – covered with black Flexi cover, height: 3 cm (not possible in combination with foot straps)

Central parking brake for rear castors incl. 4 dual castors (100 mm), operated by using the brake lever which is mounted to the seat system

Central parking brake

Central parking brake for rear castors incl. 4 dual castors (100 mm), **front with locking devices**, operated by using the brake lever which is mounted to the seat system

Price / piece

for size 1 (**2101027**)
for size 2 (**2102027**)
for size 3 (**2103027**)

133 €
133 €
133 €

Piece / set

for size 0 (**2100011**)
for size 1 (**2101011**)
for size 2 (**2102011**)

467 €
467 €
467 €

Piece / set

for size 0 (**2100111**)
for size 1 (**2101111**)
for size 2 (**2102111**)

504 €
504 €
504 €

Castors



4 Castors – 2 with locking devices (125 mm not possible for size 0)

4 Castors – all with locking devices (125 mm not possible for size 0)

Upcharge / Set

75 mm (**2104075**)
100 mm (**2104100**)
125 mm (**2104125**)

incl.
64 €
101 €

Upcharge / Set

75 mm (**2106075**)
100 mm (**2106100**)
125 mm (**2106125**)

42 €
101 €
138 €