



Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The mobile therapy chair which grows with the child.

- Incl. armrests, footplate with slots for Harness, push handle, standard seat pad with cover
- Adjustable seat height with gas pressure spring support
- 4 castors (75 mm), 2 with locking devices
- Adjustment options:
 - Tilt-in-space (via gas pressure spring 500N)
 - Back angle
 - Seat depth and seat width
 - Angle, height and width adjustable armrests
 - Height adjustable backrest
 - Folding footplate, knee angle, foot angle and height adjustable

HMV-No.: 26.11.05.0021




Basic model

madita. size 0		madita. size 1		madita. size 1b	madita. size 2
Price / piece		Price / piece		Price / piece	Price / piece
73 00 000	1.735 €	73 01 000	1.735 €	73 01 200 1.860 €	73 02 000 1.985 €
with low subframe (73 00 100)		with low subframe (73 01 100)			
	1.735 €		1.735 €		

Frame colours

Anthracit metallic		Telegrey 4 (white grey)	
☒		☐	
Price		Price	
92 00 014	incl.	92 00 011	incl.

Presettings



Presettings on the therapy chair at Schuchmann.
Please provide us with the following body measurements:

A Axilla height	_____	D Pelvic width	_____
B Chest width	_____	E Seat depth	_____
C Back height	_____	F Lower leg length	_____


7300080

Price

incl.

Dimensions

		Size 0	Size 1	Size 1b	Size 2
A	Seat depth*	15–25 cm	21–30 cm	26–36 cm	30–42 cm
B	Seat width	16–26 cm	20–30 cm	20–30 cm	28–38 cm
C	Back height with short back	35–43 cm	35–43 cm	42–52 cm	42–52 cm
C	Back height with long back	42–52 cm	42–52 cm	53–63 cm	53–63 cm
D	Seat height	46–58 cm	46–58 cm	46–58 cm	46–58 cm
D	Seat height with low subframe	35–40 cm	35–40 cm	–	–
E	Lower leg length	17–25 cm	20–32 cm	28–40 cm	29–49 cm
F	Total width	45 cm	45 cm	53 cm	53 cm
G	Length of base frame	66 cm	66 cm	77 cm	77 cm
	Tilt-in-space	(–)5° – 37°	(–)5° – 37°	(–)5° – 37°	(–)5° – 37°
	Back angle	(–)5° – 25°	(–)5° – 25°	(–)5° – 25°	(–)5° – 25°
	Max. load	30 kg	30 kg	45 kg	50 kg
	Weight	16 kg	17 kg	18 kg	19 kg



* reduced by 1 cm in combination with adjustable pelvic guide

Covers

Dark Grey / Pink
(No. 92 + 114)

Dark Grey / Azur
(No. 92 + 110)

Dark Grey / Deep orange
(No. 92 + 107)


Price

9100222

incl.

Price

9100225

incl.

Price

9100213

incl.

Covers

**Black/
Dark Grey**
(No. 105 + 92)



Price

9100205

incl.

Name embroidery

Name embroidery on the back pad
please fill out the form on
page 155!
An embroidery of
Flexi fabrics
is not possible.



Minor position

Without seat and back upholstery and cover
for the installation of custom cushions

Price reduction

7300059

Net –104 €

Coated covers

Flexi Black / Light grey
(No. 43 + 50)


Upcharge

9100226

250 €

Flexi Blue
(No. 44)


Upcharge


9100044

250 €


Configuration options - back (standard design or with integrated upper arm guide)

Back – incl. pad and cover		for size 0	for size 1	for size 1b	for size 2
Short back: Long back:		35 – 43 cm 42 – 52 cm	35 – 43 cm 42 – 52 cm	42 – 52 cm 53 – 63 cm	42 – 52 cm 53 – 63 cm
		Price / piece		Price / piece	
	Short standard design:	73 00 102	185 €	73 01 202	185 €
	Long standard design:	73 00 103	185 €	73 01 003	185 €
	Short with upper arm guide:	73 00 129	295 €	73 01 229	295 €
	Long with upper arm guide:	73 00 030	295 €	73 01 030	295 €

Configuration options - back

Back unit – fixed (standard configuration = without push handles, if nothing else stated)		Back unit – removable		Adjustable limitation of the back angle – to determine the individual max. back angle
				
Price / piece		Price / piece		Price / set
for size 0+1 with push handles	(73 01 004) incl.	for size 0+1 with push handles	(73 01 006) 370 €	73 01 111 62 €
for size 2 with push handles	(73 02 004) incl.	for size 2 with push handles	(73 02 006) 370 €	
for size 0+1 without push handles	(73 01 005) incl.	for size 0+1 without push handles	(73 01 007) 370 €	
for size 2 without push handles	(73 02 005) incl.	for size 2 without push handles	(73 02 007) 370 €	
for size 0+1 with push loop	(73 01 053) incl.			

Headrests

Headrest (standard) padded, incl. cover, height and angle adjustable	Headrest (shell-shaped) padded, incl. cover, height and angle adjustable	Headrest (lateral guide) padded, incl. cover, height and angle adjustable	Headrest (lateral guide, angle adjustable) padded, incl. cover, height, depth and angle adjustable side guides are each angle adjustable
			
Price / piece	Price / piece	Price / piece	
for size 0 (73 00 108) 254 €	for size 0+1 (73 01 110) 360 €	for size 0+1 (73 01 009) 340 €	
for size 1 (73 01 008) 254 €	for size 2 (73 02 110) 360 €	for size 2 (73 02 009) 340 €	
for size 2 (73 02 008) 254 €			
			</

Headrests
Seat angle adjustment

Universal adapter for tailored headrests in the version 15 x 15 mm	Trigger lever – mounted directly on the gas spring (standard configuration, if nothing else stated)	Trigger lever – mounted on the push bar/handle (not possible with removable back unit)
Price / piece	Price / piece	Price / piece
73 01 031 85 €	73 00 085 incl.	73 01 085 incl.

Configuration possibilities - lateral pads (Step 1: Select the bracket)

Chest width – mounted on madita. size 0	Chest width – mounted on madita. size 1 + 1b	Chest width – mounted on madita. size 2	Bracket	Upcharge / pair
13 – 20 cm	16 – 23 cm	26 – 32 cm	Fixed (73 00 012)	incl.
			Flip away (73 00 013)	259 €
20 – 27 cm	23 – 30 cm	32 – 38 cm	Fixed (73 01 012)	incl.
			Flip away (73 01 013)	259 €
___ – ___ cm	___ – ___ cm	___ – ___ cm	Fixed and tailored (73 09 012)	70 €
			Flip away and tailored (73 09 013)	330 €


Split footplates

Split footplates – individually adjustable, adjustable in height, width and angle, foldable, fixed to the seat unit, with non-slip surface and 4 slots (each) for attaching straps



Upcharge / pair

for size 0	(73 00 052)	282 €
for size 1	(73 01 052)	282 €
for size 1 b	(73 01 152)	282 €
for size 2	(73 02 052)	282 €



Split footplates with foot cups – individually adjustable, adjustable in height, width and angle, foldable, fixed to the seat unit

	Min. front width	Max. front width	Min. rear width	Max. rear width	Length
Size 1	8.7 cm	10.4 cm	5.7 cm	7.4 cm	17.4 cm
Size 2	9.5 cm	12 cm	6.7 cm	9.2 cm	20.1 cm
Size 3	11.5 cm	14 cm	8.1 cm	10.6 cm	23.8 cm

Upcharge / pair

Size 1	(73 01 064)	282 €
Size 2	(73 02 064)	282 €
Size 3	(73 03 064)	282 €

Accessories

Therapy tray – made of wood with edge



	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

Price / piece

for size 0	(73 00 024)	249 €
for size 1	(73 01 024)	249 €
for size 2	(73 02 024)	249 €

Therapy tray – transparent, with edge



	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

Price / piece

for size 0	(73 00 037)	340 €
for size 1	(73 01 037)	340 €
for size 2	(73 02 037)	340 €

4 Castors – 2 with locking devices



75 mm	(73 00 058)	incl.
100 mm	(73 01 058)	58 €
125 mm	(73 02 058)	92 €

4 Castors – all with locking devices



75 mm	(73 00 060)	35 €
100 mm	(73 01 060)	92 €
125 mm	(73 02 060)	127 €

Accessories

Wheel base extension
inserts at rear, length = 10 cm (not possible together with central parking brake)

Central parking brake for rear castors
incl. 4 dual castors (100 mm), operated by using the brake lever which is mounted to the seat system



Harness assembly kit for the seat
consisting of 2 brackets with clip buckle and 2 clip buckles (drilled)

Harness assembly kit for the back
consisting of 2 threaded plates with clip buckles and 2 clip buckles (drilled)

Price / pair

73 01 028 168 €

Parking lever mounted left	(73 00 048)	425 €
Parking lever mounted right	(73 00 548)	425 €
Parking lever mounted left; castors front with brakes	(73 00 148)	459 €
Parking lever mounted right; castors front with brakes	(73 00 648)	459 €

Upcharge / set

Upcharge / set

20 mm	(73 00 025)	44 €
25 mm	(73 01 025)	44 €

Upcharge / set

20 mm	(73 00 026)	44 €
25 mm	(73 01 026)	44 €

Accessories

Incontinence coating
in the seat area



Price / piece

73 00 020 196 €

Padded cover for footplate – covered with black Flexi cover, height: 3 cm (not possible in combination with foot straps)



Price / piece

for size 0	(73 00 038)	110 €
for size 1	(73 01 038)	110 €
for size 2	(73 02 038)	110 €

Padded cover for footplate
covered with black Flexi cover, height: 3 cm (not possible in combination with foot straps)



Price / pair

for size 0	(73 00 039)	142 €
for size 1	(73 01 039)	142 €
for size 2	(73 02 039)	142 €

Table upholstery – covered with black Flexi cover, height: 1 cm



Price / piece

for size 0	(73 00 054)	124 €
for size 1	(73 01 054)	124 €
for size 2	(73 02 054)	124 €


Harness

The **madita.**-harness is assembled with the harness assembly sets (see page 12). Further fixing aids and configuration options can be found starting on page 156.

Harness

Chest-shoulder strap – made of elastic, breathable and moisture repellent material, three straps with clips, sewn strap to the left shoulder; colour: black

	Size 0	Size 1	Size 2	Size 3
A	13 cm	13 cm	18 cm	18 cm
B	17.5 cm	19 cm	20 cm	21 cm
C	15 cm	15 cm	19 cm	28 cm
D	3.5 cm	5 cm	5 cm	5 cm
E	25.5 cm	28 cm	34 cm	44 cm
Strap width	2 cm	2 cm	2.5 cm	2.5 cm




Price / piece

Size 0 (73 00 032)	177 €
Size 1 (73 01 032)	177 €
Size 2 (73 02 032)	177 €
Size 3 (73 03 032)	177 €

Positioning jacket – with zipper, made of elastic, breathable material, shoulder straps sewn at back, pelvis straps with clips; colour: black

	Size 0	Size 1	Size 2
A	11 cm	11 cm	18 cm
B	21 cm	21 cm	22 cm
C	12 cm	12 cm	15 cm
D	2.5 cm	3.5 cm	3.5 cm
E	25.5 cm	28 cm	34 cm
Strap width	2 cm	2 cm	2.5 cm




Price / piece

Size 0 (73 00 056)	206 €
Size 1 (73 01 056)	206 €
Size 2 (73 02 056)	206 €

Harness

T-shape groin harness – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black

	Size 0	Size 1	Size 2
A	23 cm	25 cm	28 cm
B	15 cm	17 cm	20 cm
C	33 cm	35 cm	38 cm
D	24 cm	26 cm	28 cm
Strap width	2 cm	2 cm	2.5 cm




Price / piece

Size 0 (73 00 033)	148 €
Size 1 (73 01 033)	148 €
Size 2 (73 02 033)	148 €

Groin harness (standard) – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black

	Size 1	Size 2	Size 3
A	24 cm	29 cm	33 cm
B	14 cm	17 cm	22 cm
C	35 cm	43 cm	49 cm
Strap width	2 cm	2.5 cm	2.5 cm



Price / piece

Size 1 (73 01 034)	148 €
Size 2 (73 02 034)	148 €
Size 3 (73 03 034)	148 €

Harness

Pelvic harness – with pad and adjustable clip on both sides; colour: black




	Size 0	Size 1	Size 2
Pad width	4 cm	4 cm	5 cm
Pad length per side	16 cm	18 cm	28 cm
Strap length per side	30 cm	30 cm	40 cm
Strap width	2 cm	2.5 cm	2.5 cm

Price / piece

Size 0 (73 00 035)	89 €
Size 1 (73 01 035)	89 €
Size 2 (73 02 035)	89 €

4-point pelvic harness – with pad and adjustable clip on both sides; colour: black




	Size 0	Size 1	Size 2
Pad width	4 cm	4 cm	5 cm
Pad length per side	16 cm	18 cm	28 cm
Strap length per side	30 cm	30 cm	40 cm
Strap width	2 cm	2.5 cm	2.5 cm

Price / piece

Size 0 (73 00 036)	155 €
Size 1 (73 01 036)	155 €
Size 2 (73 02 036)	155 €

Harness

Ankle cuffs – with velcro fastener and idler pulley colour: black (not possible in combination with padded cover for footplate/split footplate)




	Ankle circumference	Strap width
Size 1	12 – 19 cm	25 mm
Size 2	15 – 23 cm	25 mm
Size 3	19 – 27 cm	40 mm

Price / pair

Size 1 (73 01 061)	212 €
Size 2 (73 02 061)	212 €
Size 3 (73 03 061)	212 €

Foot positioner – with straps, incl. velcro fastener, colour: black (not possible in combination with padded cover for footplate/split footplate)



Velcro strip to open

	A	Strap width
Size 1	6.5 cm	2.5 cm
Size 2	8 cm	2.5 cm
Size 3	9.5 cm	2.5 cm

Price / pair

Size 1 (78 01 034)	196 €
Size 2 (78 02 034)	196 €
Size 3 (78 03 034)	196 €