



madita-fun.

Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The modern and integrative therapy chair.

- Incl. armrests, footplate with slots for harness, push handle, standard seat pad with cover
- Height adjustment via gas pressure spring 1200 N
- 4 castors (75 mm), 2 with locking devices
- Adjustment options:
 - Tilt-in-space (via gas pressure spring 500 N)
 - Back angle
 - Seat depth and seat width
 - Angle, height and width adjustable armrests
 - Height adjustable backrest
 - Folding footplate, knee angle, foot angle and height adjustable

HMV-No.: 26.11.05.0022




Basic model

madita-fun. size 0	madita-fun. size 1	madita-fun. size 1b	madita-fun. size 2
Price / piece	Price / piece	Price / piece	Price / piece
74 00 100 2.445 €	74 01 000 2.445 €	74 01 100 2.635 €	74 02 000 2.750 €

Frame colours

Anthracit metallic	Telegrey 4 (white grey)
Price	Price
92 00 014 incl.	92 00 011 incl.

Presettings



Presettings on the therapy chair at Schuchmann.
Please provide us with the following body measurements:


A Axilla height	_____	D Pelvic width	_____
B Chest width	_____	E Seat depth	_____
C Back height	_____	F Lower leg length	_____

74 00 080

Price
incl.

Dimensions

		Size 0	Size 1	Size 1b	Size 2
A	Seat depth*	15 - 25 cm	21 - 30 cm	26 - 36 cm	30 - 42 cm
B	Seat width	16 - 26 cm	20 - 30 cm	20 - 30 cm	28 - 38 cm
C	Back height with short Back	35 - 43 cm	35 - 43 cm	42 - 52 cm	42 - 52 cm
C	Back height with long Back	42 - 52 cm	42 - 52 cm	53 - 63 cm	53 - 63 cm
D	Seat height	25 - 62 cm	25 - 62 cm	25 - 62 cm	25 - 62 cm
	Seat height with electric height adjustment	32 - 62 cm	32 - 62 cm	32 - 62 cm	32 - 62 cm
	Seat height with hydraulic height adjustment	27 - 63 cm	27 - 63 cm	27 - 63 cm	27 - 63 cm
E	Lower leg length	17 - 25 cm	20 - 32 cm	28 - 40 cm	29 - 49 cm
F	Total width	54 cm	54 cm	64 cm	64 cm
G	Length of base frame	61 cm	61 cm	75 cm	75 cm
	Total length with combi subframe	76 cm	76 cm	86 cm	86 cm
	Tilt-in-space	(-)8° - 37°	(-)8° - 37°	(-)8° - 37°	(-)8° - 37°
	Back angle	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°
	Max. load	30 kg	30 kg	45 kg	50 kg
	Weight	19 kg	20 kg	22 kg	24 kg



* reduced by 1 cm in combination with adjustable pelvic guide

Covers

Dark Grey / Pink (No. 92 + 114)		Dark Grey / Azur (No. 92 + 110)		Dark Grey / Deep orange (No. 92 + 107)	
	Price		Price		Price
9100222	incl.	9100225	incl.	9100213	incl.

Covers

Name embroidery

Minor position

Black/ Dark Grey (No. 105 + 92)		Name embroidery on the back pad please fill out the form on page 155! An embroidery of Flexi fabrics is not possible.		Without seat and back upholstery and cover for the installation of custom cushions
	Price			Price reduction
9100205	incl.			7300059 Net -104 €

Coated covers

Flexi Black/ Light grey (No. 43 + 50)		Flexi Blue (No. 44)	
	Upcharge		Upcharge
9100226	250 €	9100044	250 €


Configuration options - back (Standard design or with integrated upper arm guide)

Back – incl. pad and cover		for size 0	for size 1	for size 1b	for size 2
Short back: Long back:		35 – 43 cm 42 – 52 cm	35 – 43 cm 42 – 52 cm	42 – 52 cm 53 – 63 cm	42 – 52 cm 53 – 63 cm
		Price / piece	Price / piece	Price / piece	Price / piece
	Short standard design:	73 00 102 185 €	73 01 202 185 €	73 01 102 185 €	73 02 002 208 €
	Long standard design:	73 00 103 185 €	73 01 003 185 €	73 01 103 185 €	73 02 003 208 €
	Short with upper arm guide:	73 00 129 295 €	73 01 229 295 €	73 01 129 295 €	73 02 029 314 €
	Long with upper arm guide:	73 00 030 295 €	73 01 030 295 €	73 01 130 295 €	73 02 030 314 €

Configuration options - back

Back unit – fixed (standard configuration = with push handles, if nothing else stated)			Back unit – removable			Adjustable limitation of the back angle – to determine the individual max. back angle	
							
Price / piece			Price / piece				
for size 0+1 with push handles	(73 01 004)	incl.	for size 0+1 with push handles	(73 01 006)	370 €	73 01 111	62 €
for size 2 with push handles	(73 02 004)	incl.	for size 2 with push handles	(73 02 006)	370 €		
for size 0+1 without push handles	(73 01 005)	incl.	for size 0+1 without push handles	(73 01 007)	370 €		
for size 2 without push handles	(73 02 005)	incl.	for size 2 without push handles	(73 02 007)	370 €		
for size 0+1 with push loop	(73 01 053)	incl.					

Headrests

Headrest (standard) padded, incl. cover, height and angle adjustable	Headrest (shell-shaped) padded, incl. cover, height and angle adjustable	Headrest (lateral guide) padded, incl. cover, height and angle adjustable	Headrest (lateral guide, angle adjustable) padded, incl. cover, height, depth and angle adjustable, side guides are each angle adjustable																																
																																			
Price / piece	Price / piece	Price / piece	<table><tr><td></td><td>A</td><td>B</td><td>C</td></tr><tr><td>Size 1</td><td>13.5cm</td><td>10.5cm</td><td>10cm</td></tr><tr><td>Size 2</td><td>15cm</td><td>10.5cm</td><td>13.5cm</td></tr><tr><td>Size 3</td><td>17cm</td><td>12.5cm</td><td>16.5cm</td></tr><tr><td colspan="4">Price / piece</td></tr><tr><td>Size 1</td><td>(73 01 063)</td><td colspan="2">588 €</td></tr><tr><td>Size 2</td><td>(73 02 063)</td><td colspan="2">588 €</td></tr><tr><td>Size 3</td><td>(73 03 063)</td><td colspan="2">588 €</td></tr></table>		A	B	C	Size 1	13.5cm	10.5cm	10cm	Size 2	15cm	10.5cm	13.5cm	Size 3	17cm	12.5cm	16.5cm	Price / piece				Size 1	(73 01 063)	588 €		Size 2	(73 02 063)	588 €		Size 3	(73 03 063)	588 €	
	A	B	C																																
Size 1	13.5cm	10.5cm	10cm																																
Size 2	15cm	10.5cm	13.5cm																																
Size 3	17cm	12.5cm	16.5cm																																
Price / piece																																			
Size 1	(73 01 063)	588 €																																	
Size 2	(73 02 063)	588 €																																	
Size 3	(73 03 063)	588 €																																	
for size 0 (73 00 108) 254 € for size 1 (73 01 008) 254 € for size 2 (73 02 008) 254 €	for size 0+1 (73 01 110) 360 € for size 2 (73 02 110) 360 €	for size 0+1 (73 01 009) 340 € for size 2 (73 02 009) 340 €	A = Width contact surface B = Height headrest C = Width lateral guide																																

A = Width contact surface
B = Height headrest
C = Width lateral guide

Headrests
Seat angle adjustment

Universal adapter for tailored headrests in the version 15 x 15 mm	Trigger lever – mounted directly on the gas spring (standard configuration, if nothing else stated)	Trigger lever – mounted on the push bar/handle (not possible with removable back unit)
Price / piece	Price / piece	Price / piece
73 01 031 85 €	74 00 085 incl.	74 01 085 incl.

Configuration possibilities - lateral pads (Step 1: Select the bracket)

Chest width - mounted on madita-fun. size 0	Chest width - mounted on madita-fun. size 1 + 1b	Chest width - mounted on madita-fun. size 2	Bracket	Upcharge / pair
13 - 20 cm	16 - 23 cm	26 - 32 cm	Fixed (7300012)	incl.
			Flip away (7300013)	259 €
20 - 27 cm	23 - 30 cm	32 - 38 cm	Fixed (7301012)	incl.
			Flip away (7301013)	259 €
___ - ___ cm	___ - ___ cm	___ - ___ cm	Fixed and tailored (7309012)	70 €
			Flip away and tailored (7309013)	330 €

(Step 2: Select the lateral pads)

	Lateral pads size 0 padded, incl. cover (standard configuration for madita-fun. size 0, if nothing else stated) Height: 7 cm Depth: 9 cm	Lateral pads size 1 padded, incl. cover (standard configuration for madita-fun. size 0, if nothing else stated) Height: 8 cm Depth: 13 cm	Lateral pads size 2 padded, incl. cover (standard configuration for madita-fun. size 0, if nothing else stated) Height: 11 cm Depth: 16 cm	Tailored lateral pads padded, incl. cover Height: ___ cm Depth: ___ cm
	Price / pair	Price / pair	Price / pair	Price / pair
	7300011 106 €	7301011 106 €	7302011 106 €	7309011 208 €




(Step 2: Select the lateral pads)

	Lateral pads size 0, long (trunk guidance) padded, incl. cover Height: 13 cm Depth: 7 cm Min. adjustable underarm height in combination with adjustable pelvic guide: 22 cm	Lateral pads size 1, long (trunk guidance) padded, incl. cover Height: 15 cm Depth: 8 cm Min. adjustable underarm height in combination with adjustable pelvic guide: 24 cm	Lateral pads size 2, long (trunk guidance) padded, incl. cover Height: 17 cm Depth: 9 cm Min. adjustable underarm height in combination with adjustable pelvic guide: 26 cm	Tailored lateral pads, long (trunk guidance) padded, incl. cover Height: ___ cm Depth: ___ cm
i Lateral pads (long) with flip away bracket can be folded 20° to the side.				
	Price / pair	Price / pair	Price / pair	Price / pair
	7300055 186 €	7301055 186 €	7302055 186 €	7309055 288 €

Thigh guides

	Thigh guides for size 0 padded, incl. cover, adjustable angle, depth and width 12 cm = standard configuration for size 0	Thigh guides for size 1 + 1 b padded, incl. cover, adjustable angle, depth and width 14 cm = standard configuration for sizes 1 + 1 b	Thigh guides for size 2 padded, incl. cover, adjustable angle, depth and width 19 cm = standard configuration for size 2
	Price / pair	Price / pair	Price / pair
	12 cm (7300114) 192 €	14 cm (7301114) 192 €	19 cm (7302114) 197 €
	14 cm (7300115) 192 €	19 cm (7301115) 192 €	24 cm (7302115) 197 €

Seat option

Standard seat pad – incl. cover (standard configuration, if nothing else stated) 	Seat pad – incl. cover with anatomic leg guidance 	Adjustable pelvic guidance – padded, incl. cover  
Price / piece	Upcharge / piece	Price / piece
for size 0 (7300018) incl.	for size 0 (7300019) 217 €	for size 0 – seat width A = 15 - 24 cm (7300021) 384 €
for size 1 (7301018) incl.	for size 1 (7301019) 217 €	for size 1 + 1 b – seat width A = 18 - 28 cm (7301021) 384 €
for size 1 b (7301118) incl.	for size 1 b (7301119) 217 €	for size 2 – seat width A = 26 - 36 cm (7302021) 407 €
for size 2 (7302018) incl.	for size 2 (7302019) 217 €	


Seat option

Abductable leg guidance – with leg length compensation, padded, incl. cover (already includes thigh guides)

	for size 0	for size 1+1b	for size 2
A	0° – 20°	0° – 20°	0° – 20°
B	7 – 10 cm	9 – 12 cm	10 – 14 cm
C	3 cm	3 cm	3 cm

Price / piece

for size 0	(7300022)	520 €
for size 1	(7301022)	520 €
for size 1 b	(7301122)	520 €
for size 2	(7302022)	540 €



A = Abduction angle
B = Thigh width
C = Leg length compensation

Accessories

Knee abduction wedge – pluggable, using quick-release axle, padded, incl. cover

	Width	Depth
Size 0	7 cm	10 cm
Size 1	8 cm	14 cm
Size 2	10 cm	17 cm



Price / piece

Size 0	(7300023)	225 €
Size 1	(7301023)	225 €
Size 2	(7302023)	225 €


Split footplates

Split footplates – individually adjustable, adjustable in height, width and angle, foldable, fixed to the seat unit, with non-slip surface and 4 slots (each) for attaching straps



Upcharge / pair

for size 0	(7300052)	282 €
for size 1	(7301052)	282 €
for size 1 b	(7301152)	282 €
for size 2	(7302052)	282 €



Split footplates with foot cups – individually adjustable, adjustable in height, width and angle, foldable, fixed to the seat unit, with magnetic lock

	Min. front width	Max. front width	Min. rear width	Max. rear width	Length
Size 1	8.7 cm	10.4 cm	5.7 cm	7.4 cm	17.4 cm
Size 2	9.5 cm	12 cm	6.7 cm	9.2 cm	20.1 cm
Size 3	11.5 cm	14 cm	8.1 cm	10.6 cm	23.8 cm

Upcharge / pair

Size 1	(7301064)	282 €
Size 2	(7302064)	282 €
Size 3	(7303064)	282 €

Accessories

Therapy tray – made of wood with edge



	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

Price / piece

for size 0	(7300024)	249 €
for size 1	(7301024)	249 €
for size 2	(7302024)	249 €

Therapy tray – transparent, with edge



	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

Price / piece

for size 0	(7300037)	340 €
for size 1	(7301037)	340 €
for size 2	(7302037)	340 €

4 Castors – 2 with locking devices



4 Castors – all with locking devices



Upcharge / set

75 mm	(7300058)	incl.
100 mm	(7301058)	58 €
125 mm	(7302058)	92 €

Upcharge / set

75 mm	(7300060)	35 €
100 mm	(7301060)	92 €
125 mm	(7302060)	127 €

Accessories

Central parking brake for rear castors – incl. 4 dual castors (100 mm), operated by using the brake lever which is mounted to the seat system



Conversion kit for combi subframe

complete with 12 1/2" wheels, screw-in axles, drum brake, steering wheels 6" and push bar



Price / set

for size 0+1	Parking lever mounted left	(74 01 051)	424 €
for size 0+1	Parking lever mounted right	(74 01 551)	424 €
for size 1b+2	Parking lever mounted left	(74 02 051)	424 €
for size 1b+2	Parking lever mounted right	(74 02 551)	424 €
for size 0+1	Parking lever mounted left – castors front with brakes	(74 01 151)	459 €
for size 0+1	Parking lever mounted right – castors front with brakes	(74 01 651)	459 €
for size 1b+2	Parking lever mounted left – castors front with brakes	(74 02 151)	459 €
for size 1b+2	Parking lever mounted right – castors front with brakes	(74 02 651)	459 €

Price / set

for size 0+1	(74 01 001)	853 €
for size 1b+2	(74 02 001)	853 €

Accessories

Hydraulic height adjustment

Retrofitting not possible! Not possible together with the conversion kit for combi subframe!



Electric height adjustment

by a 24 Volt electric motor; incl. battery and hand control



Additional gas dampers (2 pcs)

to support the lift function

(not possible in combination with the hydraulic or electric height adjustment!)

2x 300 N	recommended from 30 kg load
2x 500 N	recommended from 40 kg load

Anti tip – adjustable in height



Upcharge / set

Price / set

Price / set

Price / piece

for size 0+1	(74 01 050)	537 €	74 00 102	939 €	2x 300 N	(74 00 003)	119 €	for standard & electric height adjustment	(74 00 004)	256 €
for size 1b+2	(74 02 050)	537 €			2x 500 N	(74 01 003)	119 €	for hydr. height adjustment	(74 01 004)	256 €

Accessories

Wheel base extension

insertable, inserts at front, length = 10cm



Incontinence coating

in the seat area



Harness assembly kit for the seat

consisting of 2 brackets with clip buckle and 2 clip buckles (drilled)

Harness assembly kit for the back

consisting of 2 threaded plates with clip buckles and 2 clip buckles (drilled)

Price / set

Price / piece

Price / set

Price / set

74 00 005	176 €	73 00 020	196 €	20 mm	(73 00 025)	44 €	20 mm	(73 00 026)	44 €
				25 mm	(73 01 025)	44 €	25 mm	(73 01 026)	44 €

Padding

Padded cover for footplate

covered with black Flexi cover, height: 3 cm (not possible in combination with foot straps)



Padded cover for split footplate

covered with black Flexi cover, height: 3 cm (not possible in combination with foot straps)



Padded tray cover

covered with black Flexi cover, height: 1cm



Price / piece

Price / piece

Price / piece

for size 0	(73 00 038)	110 €	for size 0	(73 00 039)	142 €	for size 0	(73 00 054)	124 €
for size 1	(73 01 038)	110 €	for size 1	(73 01 039)	142 €	for size 1	(73 01 054)	124 €
for size 2	(73 02 038)	110 €	for size 2	(73 02 039)	142 €	for size 2	(73 02 054)	124 €


Harness

The **madita-fun.**-harness is assembled with the harness assembly sets (see page 19). Further fixing aids and configuration options can be found from page 156.

Harness

Chest-shoulder strap – made of elastic, breathable and moisture repellent material, three straps with clips, sewn strap to the left shoulder; colour: black

	Size 0	Size 1	Size 2	Size 3
A	13 cm	13 cm	18 cm	18 cm
B	17.5 cm	19 cm	20 cm	21 cm
C	15 cm	15 cm	19 cm	28 cm
D	3.5 cm	5 cm	5 cm	5 cm
E	25.5 cm	28 cm	34 cm	44 cm
Strap width	2 cm	2 cm	2.5 cm	2.5 cm




Price / piece

Size 0 (73 00 032)	177 €
Size 1 (73 01 032)	177 €
Size 2 (73 02 032)	177 €
Size 3 (73 03 032)	177 €

Positioning jacket – with zipper, made of elastic, breathable material, shoulder straps sewn at back, pelvis straps with clips; colour: black

	Size 0	Size 1	Size 2
A	11 cm	11 cm	18 cm
B	21 cm	21 cm	22 cm
C	12 cm	12 cm	15 cm
D	2.5 cm	3.5 cm	3.5 cm
E	25.5 cm	28 cm	34 cm
Strap width	2 cm	2 cm	2.5 cm




Price / piece

Size 0 (73 00 056)	206 €
Size 1 (73 01 056)	206 €
Size 2 (73 02 056)	206 €

Harness

T-shape groin harness – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black

	Size 0	Size 1	Size 2
A	23 cm	25 cm	28 cm
B	15 cm	17 cm	20 cm
C	33 cm	35 cm	38 cm
D	24 cm	26 cm	28 cm
Strap width	2 cm	2 cm	2.5 cm




Price / piece

Size 0 (73 00 033)	148 €
Size 1 (73 01 033)	148 €
Size 2 (73 02 033)	148 €

Groin harness (standard) – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black

	Size 1	Size 2	Size 3
A	24 cm	29 cm	33 cm
B	14 cm	17 cm	22 cm
C	35 cm	43 cm	49 cm
Strap width	2 cm	2.5 cm	2.5 cm



Price / piece

Size 1 (73 01 034)	148 €
Size 2 (73 02 034)	148 €
Size 3 (73 03 034)	148 €

Harness

Pelvic harness – with pad and adjustable clip on both sides; colour: black



	Size 0	Size 1	Size 2
Pad width	4 cm	4 cm	5 cm
Pad length per side	16 cm	18 cm	28 cm
Strap length per side	30 cm	30 cm	40 cm
Strap width	2 cm	2.5 cm	2.5 cm

Price / piece

Size 0 (73 00 035)	89 €
Size 1 (73 01 035)	89 €
Size 2 (73 02 035)	89 €

4-point pelvic harness – with pad and adjustable clip on both sides; colour: black




	Size 0	Size 1	Size 2
Pad width	4 cm	4 cm	5 cm
Pad length per side	16 cm	18 cm	28 cm
Strap length per side	30 cm	30 cm	40 cm
Strap width	2 cm	2.5 cm	2.5 cm

Price / piece

Size 0 (73 00 036)	155 €
Size 1 (73 01 036)	155 €
Size 2 (73 02 036)	155 €

Harness

Ankle cuffs – with velcro fastener and idler pulley, colour: black (not possible in combination with padded cover for footplate/split footplate)




	Ankle circumference	Strap width
Size 1	12 – 19 cm	25 mm
Size 2	15 – 23 cm	25 mm
Size 3	19 – 27 cm	40 mm

Price / pair

Size 1 (73 01 061)	212 €
Size 2 (73 02 061)	212 €
Size 3 (73 03 061)	212 €

Foot positioner – with straps, incl. velcro fastener, colour: black (not possible in combination with padded cover for footplate/split footplate)



Velcro strip to open

	A	Strap width
Size 1	6.5 cm	2.5 cm
Size 2	8 cm	2.5 cm
Size 3	9.5 cm	2.5 cm

Price / pair

Size 1 (78 01 034)	196 €
Size 2 (78 02 034)	196 €
Size 3 (78 03 034)	196 €