



madita-fun mini.

Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The therapy chair for early aid.



- Incl. armrests, footplate with slots for harness, upholstery for footplate, ergonomically shaped push handle, standard seat pad incl. cover, 4 Castors (75 mm) – 2 with locking devices
- Height adjustable per gas pressure spring (800 N)
- Frame colour: Telegrey 4 (white grey)
- Adjustment options:
 - Tilt-in-space (by gas pressure spring 250 N, (-)5° - 37° / 0° - 42°) the trigger lever is mounted on the push handle
 - Back angle
 - Seat depth and seat width
 - Armrest, height, width and angle adjustable
 - Height adjustable backrest
 - Footplate, knee angle and height adjustable

HMV-No.: 26.11.05.0022




Dimensions



A	Seat depth	12 - 19 cm
B	Seat width	14 - 22 cm
C	Height short back	28 - 33 cm
C	Height long back	30 - 35 cm
D	Seat height	25 - 62 cm
E	Lower leg length	12 - 22 cm
F	Total width	52 cm
G	Length of base frame	58 cm
	Tilt-in-space	(-)5° - 37° / 0° - 42°
	Back angle	(-)5° - 25°
	Max. load	20 kg
	Weight	16 kg







Basic model		Tilt-in-space presettings	
madita-fun mini.		Please select a default setting ((-5° - 37° = standard configuration, if nothing else stated)	
	Price / piece	Price / piece	Price / piece
3100 000	2.674 €	(-5° - 37° (3100 001)	incl. 0° - 42° (3100 101)

Covers			
Slate grey (No. 124) 	Soft grey (No. 127) 	Flexi blue (No. 44) 	Cover according to customer request
Price	Price	Upcharge	Upcharge
9100 124	incl. 9100 127	incl. 9100 044	272 € 9100 999
			278 €


Headrests		
Headrest (standard) padded, height, depth and angle adjustable, incl. cover 	Headrest (Shell-shaped) padded, height, depth and angle adjustable, incl. cover 	Headrest (Lateral guide, angle adjustable) padded, incl. cover, height, depth and angle adjustable, side guides are each angle adjustable  A = 13,5 cm B = 10,5 cm C = 10,0 cm
Price / piece	Price / piece	Price / piece
3100 021	287 € 3100 124	405 € 3100 026
		651 €

Configuration options - back																																			
	Back incl. pad and cover (back (short) is standard unless otherwise specified)		Back with integrated upper arm guidance incl. pad and cover with adjustable hollow form in the lower trunk area through upholstery wedges (height approx. 3 cm)																																
	<table><tr><td></td><th>Height</th><th>Width</th></tr><tr><td>Back (short)</td><td>28 - 33 cm</td><td>24 cm</td></tr><tr><td>Back (long)</td><td>30 - 35 cm</td><td>24 cm</td></tr><tr><td colspan="3">Price / piece</td></tr><tr><td>Back (short)</td><td>(31 00 002)</td><td>205 €</td></tr><tr><td>Back (long)</td><td>(31 00 003)</td><td>205 €</td></tr></table>		Height	Width	Back (short)	28 - 33 cm	24 cm	Back (long)	30 - 35 cm	24 cm	Price / piece			Back (short)	(31 00 002)	205 €	Back (long)	(31 00 003)	205 €		<table><tr><td></td><th>Height</th></tr><tr><td>Back (short)</td><td>28 - 33 cm</td></tr><tr><td>Back (long)</td><td>30 - 35 cm</td></tr><tr><td colspan="2">Price / piece</td></tr><tr><td>Back (short)</td><td>(31 00 005)</td><td>231 €</td></tr><tr><td>Back (long)</td><td>(31 00 006)</td><td>231 €</td></tr></table>		Height	Back (short)	28 - 33 cm	Back (long)	30 - 35 cm	Price / piece		Back (short)	(31 00 005)	231 €	Back (long)	(31 00 006)	231 €
	Height	Width																																	
Back (short)	28 - 33 cm	24 cm																																	
Back (long)	30 - 35 cm	24 cm																																	
Price / piece																																			
Back (short)	(31 00 002)	205 €																																	
Back (long)	(31 00 003)	205 €																																	
	Height																																		
Back (short)	28 - 33 cm																																		
Back (long)	30 - 35 cm																																		
Price / piece																																			
Back (short)	(31 00 005)	231 €																																	
Back (long)	(31 00 006)	231 €																																	
Lateral pads – padded and covered, incl. bracket	Lateral pads – long (trunk guidance) padded and covered, incl. bracket	Back without pad and cover – for the installation of custom cushions	Name embroidery on the back pad please fill out the form on page 137! An embroidery of Flexi fabrics is not possible.																																
																																			
<table><tr><td>Height</td><td>7 cm</td></tr><tr><td>Depth</td><td>9 cm</td></tr><tr><td>Chest width</td><td>14 - 20 cm</td></tr><tr><td colspan="2">Price / pair</td></tr><tr><td>31 00 009</td><td>119 €</td></tr></table>	Height	7 cm	Depth	9 cm	Chest width	14 - 20 cm	Price / pair		31 00 009	119 €	<table><tr><td>Height</td><td>14 cm</td></tr><tr><td>Depth</td><td>6 cm</td></tr><tr><td>Trunk width</td><td>14 - 20 cm</td></tr><tr><td colspan="2">Price / pair</td></tr><tr><td>31 00 007</td><td>188 €</td></tr></table>	Height	14 cm	Depth	6 cm	Trunk width	14 - 20 cm	Price / pair		31 00 007	188 €	31 00 008	Netto -50 €												
Height	7 cm																																		
Depth	9 cm																																		
Chest width	14 - 20 cm																																		
Price / pair																																			
31 00 009	119 €																																		
Height	14 cm																																		
Depth	6 cm																																		
Trunk width	14 - 20 cm																																		
Price / pair																																			
31 00 007	188 €																																		

Configuration options - seat

Seat pad standard, incl. cover 	Seat pad – with anatomic leg guidance, incl. cover 		Dysplasia seat – with a seat plate (up to 10 cm abductable) in the front sitting area, internal and external leg guidance (approx. 30° abduction per side), incl. pad and cover and a 10 cm wider adjustable footplate
Price / piece 3100011 incl.	Upcharge / piece 3100012 242 €		Upcharge / piece 3100013 359 €
Thigh guides – padded and covered, angle, width, depth adjustable (included in dysplasia seat) 	Internal guide – abductable (included in dysplasia seat) 	Incontinence coating in the seat area 	Seat without pad and cover for the installation of custom cushions
Price / pair 13 cm (3100016) 216 € 16 cm (3100017) 216 €	Price / set 3100015 231 €	Price / piece 3100014 217 €	Price reduction 3100018 Netto -50 €

Accessories

Therapy tray – made of wood with rail, pluggable 	Table upholstery covered with black Flexi cover, height: 1cm 	4 Castors – 2 with locking devices 	4 Castors – all with locking devices 
Price / piece 3100022 280 €	Price / piece 3100023 138 €	Upcharge / set 75 mm (7300058) incl. 100 mm (7301058) 64 €	Upcharge / set 75 mm (7300060) 42 € 100 mm (7301060) 101 €

Accessories

Harness

Harness assembly kit for the seat

consisting of 2 brackets with clip buckle and 2 clip buckles (drilled)

Harness assembly kit for the back

consisting of 2 threaded plates with clip buckles and 2 clip buckles (drilled)

The madita-fun mini. harness is assembled with the harness assembly sets

Price / set

3100030

48 €


Price / set

3100031

48 €

Harness

Chest-shoulder strap – made of elastic, breathable and moisture repellent material, three straps with clips, sewn strap to the left shoulder; colour: black



	Size mini	Size 0
A	8,5cm	13cm
B	17,5cm	17,5cm
C	14cm	15cm
D	3,5cm	3,5cm
E	24,5cm	25,5cm
Strap width	2cm	2cm

Price / piece


Size mini (73 05 032)

197 €

Size 0 (3100032)

197 €

T-shape groin harness – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black



	Size mini	Size 0
A	18cm	23cm
B	10cm	15cm
C	25cm	33cm
D	19cm	24cm
Strap width	2cm	2cm

Price / piece

Size mini (73 05 033)

165 €

Size 0 (73 00 033)

165 €

Harness

Pelvic harness – with pad and adjustable clip on both sides; colour: black



	Size mini	Size 0
Pad width	3,5cm	4cm
Pad length per side	10cm	16cm
Strap length per side	20cm	30cm
Strap width	2cm	2cm

Price / piece


Size mini (73 05 035)

99 €

Size 0 (73 00 035)

99 €

4-point pelvic harness – with pad and adjustable clip on both sides; colour: black



	Size mini	Size 0
Pad width	3,5cm	4cm
Pad length per side	10cm	16cm
Strap length per side	20cm	30cm
Strap width	2cm	2cm

Price / piece

Size mini (73 05 036)


172 €

Size 0 (73 00 036)

172 €


Harness

Ankle cuffs – with velcro fastener and idler pulley; colour: black (not possible in combination with padded cover for footplate)



	Ankle circumference	Width Gurt
Size 1	12 - 19cm	25mm
Size 2	15 - 23cm	25mm
		Price / pair
Size 1 (73 01 061)		236 €
Size 2 (73 02 061)		236 €

Foot positioner – with straps, incl. velcro fastener; colour: black (not possible in combination with padded cover for footplate)



	A	Strap width
Size 0	5cm	2,5cm
Size 1	6,5cm	2,5cm
		Price / pair
Size 0 (78 00 034)		217 €
Size 1 (78 01 034)		217 €