



madita-fun.

15

Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The modern and integrative therapy chair.

- Incl. armrests, footplate with slots for harness, push handle, standard seat pad with cover
- Height adjustment via gas pressure spring (1200 N for size 0 + 1; 1600 N for size 1b + 2)
- 4 castors (75 mm), 2 with locking devices
- Adjustment options:
 - Tilt-in-space (via gas pressure spring 500 N)
 - Back angle
 - Seat depth and seat width
 - Angle, height and width adjustable armrests
 - Height adjustable backrest
 - Folding footplate, knee angle and height adjustable

HMV-No.: 26.11.05.0022




Basic model

madita-fun. size 0	madita-fun. size 1	madita-fun. size 1b	madita-fun. size 2
Price / piece	Price / piece	Price / piece	Price / piece
79 00 000 2.711 €	79 01 000 2.711 €	79 01 200 2.922 €	79 02 000 3.049 €

Frame colours

Anthracite grey	Telegrey 4 (white grey)
	
Price	Price
92 00 017 incl.	92 00 011 incl.

Presettings



Presettings on the therapy chair at Schuchmann.
Please provide us with the following body measurements:

I Axilla height	_____	IV Pelvic width	_____
II Chest width	_____	V Seat depth	_____
III Back height	_____	VI Lower leg length	_____


74 00 080

Price

incl.

Dimensions

		Size 0	Size 1	Size 1b	Size 2
A	Seat depth*	15 - 25 cm	21 - 30 cm	26 - 36 cm	30 - 42 cm
B	Seat width	16 - 26 cm	20 - 30 cm		28 - 38 cm
C	Back height with extra short back	32 - 35 cm		35 - 43 cm	
C	Back height with short back	35 - 43 cm		42 - 52 cm	
C	Back height with long back	42 - 52 cm		53 - 63 cm	
C	Back height with Active back	35 - 43 cm		42 - 52 cm	
D	Seat height		25 - 62 cm		
	Seat height with electric height adjustment		32 - 62 cm		
	Seat height with hydraulic height adjustment		27 - 63 cm		
E	Lower leg length	12 - 22 cm** / 15 - 25 cm	16 - 29 cm** / 19 - 32 cm	19 - 33 cm** / 22 - 36 cm	25 - 43 cm
F	Total width	54 cm		64 cm	
G	Length of base frame	65 cm		75 cm	
	Total width with combi subframe	66 cm		76 cm	
	Total length with combi subframe	76 cm		86 cm	
	Tilt-in-space		(-)8° - 37°		
	Back angle		(-)5° - 25°		
	Max. load	30 kg		45 kg	50 kg
	Weight	19 kg	20 kg	22 kg	24 kg



* reduced by 1 cm in combination with adjustable pelvic guide

** this adjustment range can be achieved by repositioning the bracket.

Covers

Pink 2
(No. 115)


Price

9100 115

incl.

Deep sea blue
(No. 111)


Price

9100 111

incl.

Leaf green
(No. 123)


Price

9100 123

incl.

Slate grey
(No. 124)


Price

9100 124

incl.

Name embroidery

Name embroidery on the back pad

please fill out the form on page 137!



An embroidery of Flexi fabrics is not possible.

Covers

Flexi blue
(No. 44)


Upcharge

9100 044

272 €

Flexi light grey
(No. 50)


Upcharge

9100 050

272 €

Cover according to customer request

Upcharge

9100 999

278 €

Minor position

Without seat and back upholstery and cover
for the installation of custom cushions

Price reduction

73 00 059

Net -100 €

Headrests

Headrest (standard)
incl. standard holder, padded, incl. cover, height and angle adjustable


Price / piece

for size 0 **(78 00 009)**

281 €

for size 1 **(78 01 009)**

281 €

for size 2 **(78 02 009)**

281 €

Headrest (shell-shaped)
incl. standard holder, padded, incl. cover, height and angle adjustable


Price / piece

for size 0+1 **(78 01 010)**

399 €

for size 2 **(78 02 010)**

399 €

Headrest (lateral guide)
incl. standard holder, padded, incl. cover, height and angle adjustable


Price / piece

for size 0+1 **(78 01 011)**


377 €

for size 2 **(78 02 011)**

377 €

Headrests

Headrest (Lateral guide, angle adjustable) – incl. standard holder, padded, incl. cover, height, depth and angle adjustable, side guides are each angle adjustable



A = Width contact surface
B = Height headrest
C = Width lateral guide

	A	B	C
Size 1	13,5 cm	10,5 cm	10 cm
Size 2	15 cm	10,5 cm	13,5 cm
Size 3	17 cm	12,5 cm	16,5 cm
Price / piece			
Size 1 (78 01 012)	651 €		
Size 2 (78 02 012)	651 €		
Size 3 (78 03 012)	651 €		

Adapter for headrests for lateral adjustment

Double joint for headrest mounts to extend the depth adjustment range




Universal adapter for tailored headrests in the version 15 x 15 mm



Price / piece		Price / set		Price / piece	
78 01 041	154 €	78 00 041	98 €	78 00 013	95 €


Standard back (optional with screw-on upper arm guide)

Back – incl. pad and cover (short back is standard unless otherwise specified)



	for size 0	for size 1	for size 1b	for size 2
extra short back:	32 – 35 cm	32 – 35 cm	35 – 43 cm	35 – 43 cm
short back:	35 – 43 cm	35 – 43 cm	42 – 52 cm	42 – 52 cm
long back:	42 – 52 cm	42 – 52 cm	53 – 63 cm	53 – 63 cm
Price / piece				
extra short back:	78 00 001 205 €	78 01 001 205 €	78 01 201 205 €	78 02 001 231 €
short back:	78 00 002 205 €	78 01 002 205 €	78 01 202 205 €	78 02 002 231 €
long back:	78 00 003 205 €	78 01 003 205 €	78 01 203 205 €	78 02 003 231 €

Attachable upper arm guide (for standard back)



	for size 0	for size 1	for size 1b	for size 2
Price / piece				
for extra short back:	78 00 004 133 €	78 01 004 133 €	78 01 204 133 €	78 02 004 133 €
for short back:	78 00 005 133 €	78 01 005 133 €	78 01 205 133 €	78 02 005 133 €
for long back:	78 00 006 133 €	78 01 006 133 €	78 01 206 133 €	78 02 006 133 €

Active back (upper arm guidance not possible)

Active back
incl. pad and cover



for size 0	for size 1	for size 1b	for size 2
35 – 43 cm	35 – 43 cm	35 – 43 cm	42 – 52 cm
Price / piece			
78 00 007	261 €	78 01 007	261 €
		78 02 007	261 €

Configuration options - back

Back unit – fixed (standard configuration = with push handles, if nothing else stated)



Back unit – removable



Adjustable limitation of the back angle – to determine the individual max. back angle

Price / piece		Price / piece		Price / set	
for size 0+1 with push handles (73 01 004)	incl.	for size 0+1 with push handles (73 01 006)	411 €	73 01 111	69 €
for size 2 with push handles (73 02 004)	incl.	for size 2 with push handles (73 02 006)	411 €		
for size 0+1 with push loop (73 01 053)	incl.				

Armrest



Armrest pad size 1

(Standard configuration for basic model size 0, if nothing else stated)

Length: 20 cm Width: 6 cm

Price / pair	
78 01 008	incl.

Armrest pad size 2

(Standard configuration for basic model size 1+2, if nothing else stated)

Length: 27 cm Width: 7 cm

Price / pair	
78 02 008	incl.

Armrest pad size 3

Length: 34 cm Width: 8 cm

Price / pair	
78 03 008	incl.

Lateral pads (Fold-down lateral pads cannot be used in combination with the upper arm support.)


Lateral pads size 0
padded, incl. cover
(Standard configuration for **madita**.
size 0, if nothing else stated)

Height: 7 cm Depth: 9 cm
Chest width: 12 – 27 cm


Lateral pads size 1
padded, incl. cover
(Standard configuration for **madita**.
size 1+1b, if nothing else stated)

Height: 8 cm Depth: 13 cm
Chest width: 16 – 30 cm

Lateral pads size 2
padded, incl. cover
(Standard configuration for **madita**.
size 2, if nothing else stated)

Height: 11 cm Depth: 16 cm
Chest width: 19 – 36 cm

Price / pair		Price / pair		Price / pair	
Fixed	(7800045) 120 €	Fixed	(7801045) 120 €	Fixed	(7802045) 120 €
Fixed active back	(7800047) 120 €	Fixed active back	(7801047) 120 €	Fixed active back	(7802047) 120 €
Flip away	(7800049) 405 €	Flip away	(7801049) 405 €	Flip away	(7802049) 405 €

Lateral pads (Fold-down lateral pads cannot be used in combination with the upper arm support.)


**Lateral pads size 0, long
(trunk guidance)**
padded, incl. cover

Height: 12,5 cm Depth: 9 cm
Min. adjustable forearm height in
combination with adjustable pelvic
guide: 23 cm

**Lateral pads size 1, long
(trunk guidance)**
padded, incl. cover

Height: 15,5 cm Depth: 13 cm
Min. adjustable forearm height in
combination with adjustable pelvic
guide: 26 cm

**Lateral pads size 2, long
(trunk guidance)**
padded, incl. cover


Height: 18,5 cm Depth: 16 cm
Min. adjustable forearm height in
combination with adjustable pelvic
guide: 30 cm

Price / pair		Price / pair		Price / pair	
Fixed	(7800046) 185 €	Fixed	(7801046) 185 €	Fixed	(7802046) 185 €
Fixed active back	(7800048) 185 €	Fixed active back	(7801048) 185 €	Fixed active back	(7802048) 185 €
Flip away	(7800050) 470 €	Flip away	(7801050) 470 €	Flip away	(7802050) 470 €

i Lateral pads (long) with flip away bracket can be folded 20° to the side.

Cover

**Cover for upper arm
system/pad holder**



Price / pair
7801040 41 €

Seat angle adjustment

Trigger lever – mounted
on the push bar/handle
(standard configuration, if
nothing else stated)
This article is not possible with
the removable backrest. For
the removable backrest, please
select the „**Trigger lever with
detachable Bowden cable**
– mounted on the push bar/
handle“.

Trigger lever – mounted
directly on the gas spring


**Trigger lever with detachable
Bowden cable** – mounted on
the push bar/handle (for
the removable back unit)

Safety bolt – for release
lever, prevents unintentional
release of height and angle
adjustments
(Not possible with release
lever „**directly on the gas
spring**“)


Price / piece	Price / piece	Price / piece	Price / piece
7401085 incl.	7400085 incl.	7402085 39 €	7301086 82 €

Seat option


Standard seat pad – incl. cover
(standard configuration, if nothing
else stated)



Seat pad – incl. cover with
anatomic leg guidance




Adjustable pelvic guidance – padded, incl. cover



Price / piece		Upcharge / piece		Price / piece	
for size 0	(7800017) incl.	for size 0	(7800018) 242 €	for size 0 – seat width A = 20 – 24 cm	(7800020) 425 €
for size 1	(7801017) incl.	for size 1	(7801018) 242 €	for size 1+1b – seat width A = 22 – 26 cm	(7801020) 425 €
for size 1b	(7801217) incl.	for size 1b	(7801218) 242 €	for size 2 – seat width A = 28 – 34 cm	(7802020) 452 €
for size 2	(7802017) incl.	for size 2	(7802018) 242 €		


Seat option				Accessories		
Abductable leg guidance with leg length compensation, padded, incl. cover (already includes thigh guides)				Knee abduction wedge pluggable using quick-release axle, padded, incl. cover (not usable with abductable leg guidance)		
	for size 0	for size 1+1b	for size 2		Width	Depth
A	0° – 20°	0° – 20°	0° – 20°	Size 0	7 cm	10 cm
B	7 – 10 cm	9 – 12 cm	10 – 14 cm	Size 1	8 cm	14 cm
C	3 cm	3 cm	3 cm	Size 2	10 cm	17 cm
Price / piece				Price / piece		
for size 0	(78 000 019)		576 €	Size 0	(73 000 023)	249 €
for size 1	(78 010 019)		576 €	Size 1	(73 010 023)	249 €
for size 1b	(78 012 019)		576 €	Size 2	(73 020 023)	249 €
for size 2	(78 020 019)		598 €			






A = Abduction angle
B = Thigh width
C = Leg length compensation



Thigh guides					
		Thigh guides for size 0 padded, incl. cover, adjustable angle, depth and width		Thigh guides for size 1 + 1b padded, incl. cover, adjustable angle, depth and width	
		Price / pair		Price / pair	
		12 cm	(73 00 114) 213 €	14 cm	(73 01 114) 213 €
		14 cm	(73 00 115) 213 €	17 cm	(73 01 115) 213 €
		17 cm		22 cm	(73 02 114) 218 €
					(73 02 115) 218 €

Footplate						Calf plate			
				Width	Depth	Lower leg length			
			for size 0	28,4 cm	17,5 cm	12 – 22 cm* / 15 – 25 cm			
			for size 1	33,4 cm	20,5 cm	16 – 29 cm* / 19 – 32 cm			
			for size 1b			19 – 33 cm* / 22 – 36 cm			
			for size 2	40,3 cm	24,5 cm	25 – 43 cm			
						* This adjustment range can be achieved by repositioning the bracket.			
Footplate standard – without locking mechanism, foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps (Standard configuration, unless otherwise specified)			Footplate with lock – foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps.			Footplate with angle adjustment foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps.			
						Foot angle = (-)10° – 10°			
Price / piece			Upcharge / piece			Upcharge / piece			
for size 0	(78 000 021)	incl.	for size 0	(78 000 022)	58 €	for size 0	(78 000 023)	106 €	
for size 1	(78 010 021)	incl.	for size 1	(78 010 022)	58 €	for size 1	(78 010 023)	106 €	
for size 1b	(78 012 021)	incl.	for size 1b	(78 012 022)	58 €	for size 1b	(78 012 023)	106 €	
for size 2	(78 020 021)	incl.	for size 2	(78 020 022)	58 €	for size 2	(78 020 023)	106 €	
						Calf plate for footplates covered with black Flexi fabric			
									
							for size 0	for size 1+1b	for size 2
						Width	27 cm	32 cm	38 cm
						Height	6,5 cm	6,5 cm	6,5 cm
						Price / piece			
						for size 0	(78 000 024)	228 €	
						for size 1 + 1b	(78 010 024)	228 €	
						for size 2	(78 020 024)	228 €	

Configuration possibilities footplates					
Step 1: Select the bracket Knee angle to accommodate the footplates individually adjustable in height, width and angle			Step 2: Select the footplates		
			Split footplates – foldable, fixed mounted on the footplate holder, with non-slip surface and each 4 slots for the for the fixation of straps		
					
			Foot angle = (-)15° – 15°		
					
			Foot cups – Foldable, fixed mounted on the footplate holder, with magnetic lock		
Price / piece			Price / piece		
for size 0	(78 000 025)	316 €	for size 0	(78 000 024)	228 €
for size 1	(78 010 025)	316 €	for size 1 + 1b	(78 010 024)	228 €
for size 1b	(78 012 025)	316 €	for size 2	(78 020 024)	228 €
for size 2	(78 020 025)	316 €			

Accessories

Therapy tray – made of wood with rail, flip away



	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

Therapy tray – transparent, with rail, flip away



	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

Therapy tray – made of wood with rail, prefixable, fixed on both sides



	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

Therapy tray – transparent, with rail, prefixable, fixed on both sides



	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

	Price / piece
for size 0 (78 00 028)	276 €
for size 1 (78 01 028)	276 €
for size 2 (78 02 028)	276 €

	Price / piece
for size 0 (78 00 029)	377 €
for size 1 (78 01 029)	377 €
for size 2 (78 02 029)	377 €

	Price / piece
for size 0 (78 00 038)	294 €
for size 1 (78 01 038)	294 €
for size 2 (78 02 038)	294 €

	Price / piece
for size 0 (78 00 039)	421 €
for size 1 (78 01 039)	421 €
for size 2 (78 02 039)	421 €

Add-on option:

Protection against unintentional folding up **(78 01 036)** 108 €

Accessories

4 Castors – 2 with locking devices



	Upcharge / set
75 mm (73 00 058)	incl.
100 mm (73 01 058)	64 €
125 mm (73 02 058)	101 €

4 Castors – all with locking devices



	Upcharge / set
75 mm (73 00 060)	42 €
100 mm (73 01 060)	101 €
125 mm (73 02 060)	138 €

Accessories

Central parking brake for rear castors – incl. 4 dual castors (100 mm), operated by using the brake lever which is mounted to the seat system



Conversion kit for combi subframe

complete with 12 1/2" wheels, screw-in axles, drum brake, steering wheels 6" and push bar



	Price / set
for size 0+1 Parking lever mounted left (79 00 030)	467 €
for size 0+1 Parking lever mounted right (79 01 030)	467 €
for size 1b+2 Parking lever mounted left (79 02 030)	467 €
for size 1b+2 Parking lever mounted right (79 03 030)	467 €
for size 0+1 Parking lever mounted left – castors front with brakes (79 00 031)	504 €
for size 0+1 Parking lever mounted right – castors front with brakes (79 01 031)	504 €
for size 1b+2 Parking lever mounted left – castors front with brakes (79 02 031)	504 €
for size 1b+2 Parking lever mounted right – castors front with brakes (79 03 031)	504 €

	Price / set
for size 0+1 (74 01 001)	933 €
for size 1b+2 (74 02 001)	933 €

Accessories

Hydraulic height adjustment

Retrofitting not possible! Not possible together with the conversion kit for combi subframe!



Electric height adjustment

by a 24 Volt electric motor; incl. battery and hand control



Additional gas dampers (2 pcs)

to support the lift function

(not possible in combination with the hydraulic or electric height adjustment!)

2x 300 N	empfohlen ab 30 kg Belastung
2x 500 N	empfohlen ab 40 kg Belastung

Anti tip – adjustable in height



	Upcharge / set	Upcharge / set	Price / set	Price / piece
for size 0+1 (74 01 050)	589 €	74 00 102	1.034 €	for standard & electric height adjustment (74 00 004) 281 €
for size 1b+2 (74 02 050)	589 €			for hydr. height adjustment (74 01 004) 281 €

Padding

Padded cover for footplate

covered with black Flexi cover, height: 1cm (not possible in combination with foot straps)

Padded cover for split footplate

covered with black Flexi cover, height: 1cm (not possible in combination with foot straps)



Table upholstery – covered with black Flexi cover, height: 1cm

	Price / piece	Price / pair	Price / piece
for size 0 (78 00 032)	133 €	for size 1 (78 01 033)	191 €
for size 1+1b (78 01 032)	133 €	for size 2 (78 02 033)	191 €
for size 2 (78 02 032)	133 €		

Accessories



Incontinence coating in the seat area	Wheel base extension inserts at rear, length = 10 cm (not possible together with central parking brake)	Harness assembly kit for the seat consisting of 2 brackets with clip buckle and 2 clip buckles (drilled)	Harness assembly kit for the back consisting of 2 threaded plates with clip buckles and 2 clip buckles (drilled)
Price / piece	Price / pair	Price / set	Price / set
73 00 020 217 €	74 00 005 191 €	25 mm (73 01 025) 48 €	25 mm (73 01 026) 48 €

Harness



Chest-shoulder strap – made of elastic, breathable and moisture repellent material, three straps with clips, sewn strap to the left shoulder; colour: black					Positioning jacket – with zipper, made of elastic, breathable material, shoulder straps sewn at back, pelvis straps with clips; colour: black						
	Size 0	Size 1	Size 2	Size 3			Size 0	Size 1	Size 2		
A	13cm	13cm	18cm	18cm		A	11cm	11cm	18cm		
B	17,5cm	19cm	20cm	21cm		B	21cm	21cm	22cm		
C	15cm	15cm	19cm	28cm		C	12cm	12cm	15cm		
D	3,5cm	5cm	5cm	5cm		D	2,5cm	3,5cm	3,5cm		
E	25,5cm	28cm	34cm	44cm		E	25,5cm	28cm	34cm		
Strap width	2cm	2cm	2,5cm	2,5cm		Strap width	2cm	2cm	2,5cm		
Price / piece						Price / piece					
Size 0 (73 00 032)					197 €	Size 0 (73 00 056)					228 €
Size 1 (73 01 032)					197 €	Size 1 (73 01 056)					228 €
Size 2 (73 02 032)					197 €	Size 2 (73 02 056)					228 €
Size 3 (73 03 032)					197 €						

Harness

Harness

Pelvic harness – with pad and adjustable clip on both sides; colour: black					4-point pelvic harness – with pad and adjustable clip on both sides; colour: black				
									
	Size 0	Size 1	Size 2			Size 0	Size 1	Size 2	
Pad width	4 cm	4 cm	5 cm		Pad width	4 cm	4 cm	5 cm	
Pad length per side	16 cm	18 cm	28 cm		Pad length per side	16 cm	18 cm	28 cm	
Strap length per side	30 cm	30 cm	40 cm		Strap length per side	30 cm	30 cm	40 cm	
Strap width	2 cm	2,5 cm	2,5 cm		Strap width	2 cm	2,5 cm	2,5 cm	
Price / piece					Price / piece				
Size 0 (73 00 035) 99 €					Size 0 (73 00 036) 172 €				
Size 1 (73 01 035) 99 €					Size 1 (73 01 036) 172 €				
Size 2 (73 02 035) 99 €					Size 2 (73 02 036) 172 €				

Harness

Ankle cuffs – with velcro fastener and idler pulley, colour: black (not possible in combination with padded cover for footplate/split footplate)					Foot positioner – with straps, incl. velcro fastener, colour: black (not possible in combination with padded cover for footplate/split footplate)				
									
		Ankle circumference	Strap width				A	Strap width	
Size 1	12 - 19 cm	2,5 cm			Size 1	6,5 cm	2,5 cm		
Size 2	15 - 23 cm	2,5 cm			Size 2	8 cm	2,5 cm		
Size 3	19 - 27 cm	4 cm			Size 3	9,5 cm	2,5 cm		
Price / pair					Price / pair				
Size 1 (73 01 061) 236 €					Size 1 (78 01 034) 217 €				
Size 2 (73 02 061) 236 €					Size 2 (78 02 034) 217 €				
Size 3 (73 03 061) 236 €					Size 3 (78 03 034) 217 €				