



Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The mobile therapy chair which grows with the child.

- Incl. armrests, footplate with slots for harness, push handle, standard seat pad with cover
- Adjustable seat height with gas pressure spring support
- 4 castors (75 mm), 2 with locking devices
- Adjustment options:
 - Tilt-in-space (via gas pressure spring 500N)
 - Back angle
 - Seat depth and seat width
 - Angle, height and width adjustable armrests
 - Height adjustable backrest
 - Folding footplate, knee angle and height adjustable

HMV-No.: 26.11.05.0021




Basic model

madita. size 0	madita. size 1	madita. size 1b	madita. size 2
Price / piece	Price / piece	Price / piece	Price / piece
78 00 000 2.001 €	78 01 000 2.001 €	78 01 200 2.144 €	78 02 000 2.288 €

Frame colours

Anthracite grey		Telegrey 4 (white grey)	
Price	Price	Price	Price
92 00 017 incl.		92 00 011 incl.	

Presettings



Presettings on the therapy chair at Schuchmann.
Please provide us with the following body measurements:

I Axilla height	_____	IV Pelvic width	_____
II Chest width	_____	V Seat depth	_____
III Back height	_____	VI Lower leg length	_____


73 00 080	Price
	incl.

Dimensions

		Size 0	Size 1	Size 1b	Size 2
A	Seat depth*	15 - 25 cm	21 - 30 cm	26 - 36 cm	30 - 42 cm
B	Seat width	16 - 26 cm	20 - 30 cm		28 - 38 cm
C	Back height with extra short back	32 - 35 cm		35 - 43 cm	
C	Back height with short back	35 - 43 cm		42 - 52 cm	
C	Back height with long back	42 - 52 cm		53 - 63 cm	
C	Back height with Active back	35 - 43 cm			42 - 52 cm
D	Seat height	47,5 - 58 cm			
E	Lower leg length	12 - 22 cm** / 15 - 25 cm	16 - 29 cm** / 19 - 32 cm	19 - 33 cm** / 22 - 36 cm	25 - 43 cm
F	Total width	45 cm		53 cm	
G	Length of base frame	66 cm		77 cm	
	Tilt-in-space	(-)5° - 37°			
	Back angle	(-)5° - 25°			
	Max. load	30 kg		45 kg	50 kg
	Weight	16 kg	17 kg	18 kg	19 kg

* reduced by 1 cm in combination with adjustable pelvic guide

** this adjustment range can be achieved by repositioning the bracket.



Covers

Pink 2
(No. 115)



Deep sea blue
(No. 111)



Leaf green
(No. 123)



Slate grey
(No. 124)



Name embroidery

Name embroidery on the back pad

please fill out the form on page 137!



An embroidery of Flexi fabrics is not possible.

9100 115	Price	9100 111	Price	9100 123	Price	9100 124	Price
	incl.		incl.		incl.		incl.

Covers

Flexi blue
(No. 44)



Flexi light grey
(No. 50)



new

Cover according to customer request

9100 044	Upcharge	9100 050	Upcharge	9100 999	Upcharge
	272 €		272 €		278 €

Minor position

Without seat and back upholstery and cover
for the installation of custom cushions

73 00 059	Price reduction
	Net -100 €

Headrests

Headrest (standard)

incl. standard holder, padded, incl. cover, height and angle adjustable



Headrest (shell-shaped)

incl. standard holder, padded, incl. cover, height and angle adjustable



Headrest (lateral guide)


incl. standard holder, padded, incl. cover, height and angle adjustable



Price / piece		Price / piece		Price / piece	
for size 0 (78 00 009)	281 €	for size 0+1 (78 01 010)	399 €	for size 0+1 (78 01 011)	377 €
for size 1 (78 01 009)	281 €	for size 2 (78 02 010)	399 €	for size 2 (78 02 011)	377 €
for size 2 (78 02 009)	281 €				

Headrests

Headrest (Lateral guide, angle adjustable) – incl. standard holder, padded, incl. cover, height, depth and angle adjustable, side guides are each angle adjustable



A = Width Anlagefläche
B = Height Kopfstütze
C = Width Seitenführung

	A	B	C
Size 1	13,5 cm	10,5 cm	10 cm
Size 2	15 cm	10,5 cm	13,5 cm
Size 3	17 cm	12,5 cm	16,5 cm
Price / piece			
Size 1 (78 01 012)	651 €		
Size 2 (78 02 012)	651 €		
Size 3 (78 03 012)	651 €		

Adapter for headrests for lateral adjustment

Double joint for headrest mounts to extend the depth adjustment range




Universal adapter for tailored headrests 15 x 15 mm



Price / piece		Price / set		Price / piece	
7801041	154 €	7800041	98 €	7800013	95 €


Standard back (optional with screw-on upper arm guide)

Back – incl. pad and cover (short back is standard unless otherwise specified)



	for size 0	for size 1	for size 1b	for size 2
extra short back:	32 – 35 cm	32 – 35 cm	35 – 43 cm	35 – 43 cm
short back:	35 – 43 cm	35 – 43 cm	42 – 52 cm	42 – 52 cm
long back:	42 – 52 cm	42 – 52 cm	53 – 63 cm	53 – 63 cm
Price / piece		Price / piece	Price / piece	Price / piece
extra short back:	78 00 001 205 €	78 01 001 205 €	78 01 201 205 €	78 02 001 231 €
short back:	78 00 002 205 €	78 01 002 205 €	78 01 202 205 €	78 02 002 231 €
long back:	78 00 003 205 €	78 01 003 205 €	78 01 203 205 €	78 02 003 231 €

Attachable upper arm guide (for standard back)



	for size 0	for size 1	for size 1b	for size 2
Price / piece	Price / piece	Price / piece	Price / piece	Price / piece
for extra short back:	78 00 004 133 €	78 01 004 133 €	78 01 204 133 €	78 02 004 133 €
for short back:	78 00 005 133 €	78 01 005 133 €	78 01 205 133 €	78 02 005 133 €
for long back:	78 00 006 133 €	78 01 006 133 €	78 01 206 133 €	78 02 006 133 €

Active back (upper arm guidance not possible)

Active back
incl. pad and cover



for size 0	for size 1	for size 1b	for size 2
35 – 43 cm	35 – 43 cm	35 – 43 cm	42 – 52 cm
Price / piece		Price / piece	Price / piece
78 00 007	261 €	78 01 007	261 €
			78 02 007 261 €

Configuration options - back

Back unit – fixed (standard configuration = with push handles, if nothing else stated)



Back unit – removable



Adjustable limitation of the back angle
to determine the individual max. back angle

Price / piece		Price / piece	Price / set
for size 0+1 with push handles (73 01 004)	incl.	for size 0+1 with push handles (73 01 006)	411 €
for size 2 with push handles (73 02 004)	incl.	for size 2 with push handles (73 02 006)	411 €
for size 0+1 with push loop (73 01 053)	incl.		73 01 111 69 €

Armrest



Armrest pad size 1

(Standard configuration for basic model size 0, if nothing else stated)

Length: 20 cm Width: 6 cm

Armrest pad size 2

(Standard configuration for basic model size 1+2, if nothing else stated)

Length: 27 cm Width: 7 cm

Armrest pad size 3


Length: 34 cm Width: 8 cm


Price / pair		Price / pair	Price / pair
78 01 008	incl.	78 02 008	incl.
		78 03 008	incl.

Lateral pads (Fold-away lateral pads cannot be used in combination with the upper arm support.)

	Lateral pads size 0 padded, incl. cover (Standard configuration for madita. size 0, if nothing else stated)	Lateral pads size 1 padded, incl. cover (Standard configuration for madita. size 1+1b, if nothing else stated)	Lateral pads size 2 padded, incl. cover (Standard configuration for madita. size 2, if nothing else stated)
	Height: 7 cm Depth: 9 cm Chest width: 12 – 27 cm	Height: 8 cm Depth: 13 cm Chest width: 16 – 30 cm	Height: 11 cm Depth: 16 cm Chest width: 19 – 36 cm
	Price / pair	Price / pair	Price / pair
	Fixed (7800045) 120 €	Fixed (7801045) 120 €	Fixed (7802045) 120 €
	Fixed active back (7800047) 120 €	Fixed active back (7801047) 120 €	Fixed active back (7802047) 120 €
	Flip away (7800049) 405 €	Flip away (7801049) 405 €	Flip away (7802049) 405 €

Lateral pads (Fold-away lateral pads cannot be used in combination with the upper arm support.)

	Lateral pads size 0, long (Trunk guidance) padded, incl. cover Height: 12,5 cm Depth: 9 cm Chest width: : 12 – 27 cm Min. adjustable forearm height in combination with adjustable pelvic guide: 23 cm	Lateral pads size 1, long (Trunk guidance) padded, incl. cover Height: 15,5 cm Depth: 13 cm Chest width: 16 – 30 cm Min. adjustable forearm height in combination with adjustable pelvic guide: 26 cm	Lateral pads size 2, long (Trunk guidance) padded, incl. cover Height: 18,5 cm Depth: 16 cm Chest width: 19 – 36 cm Min. adjustable forearm height in combination with adjustable pelvic guide: 30 cm
	i Lateral pads (long) with flip away bracket can be folded 20° to the side.		
	Price / pair	Price / pair	Price / pair
	Fixed (7800046) 185 €	Fixed (7801046) 185 €	Fixed (7802046) 185 €
	Fixed active back (7800048) 185 €	Fixed active back (7801048) 185 €	Fixed active back (7802048) 185 €
	Flip away (7800050) 470 €	Flip away (7801050) 470 €	Flip away (7802050) 470 €

Cover
Cover for upper arm system/pad holder


Price / pair
7801040 41 €

Seat angle adjustment

Trigger lever – mounted on the push bar/handle (standard configuration, if nothing else stated)

This is not possible with the removable backrest. For the removable backrest, please select the „**Trigger lever with detachable Bowden cable – mounted on the push bar/handle**“.

Trigger lever – mounted directly on the gas spring

Trigger lever with detachable Bowden cable – mounted on the push bar/handle (for the removable back unit)

Safety bolt – for release lever, prevents unintentional release of height and angle adjustments (Not possible with release lever „**directly on the gas spring**“)

Price / piece
7301085 incl.


Price / piece
7300085 incl.

Price / piece
7302085 39 €


Price / piece
7301086 82 €

Seat option


Standard seat pad – incl. cover (standard configuration, if nothing else stated)



Seat pad – incl. cover with anatomic leg guidance



Adjustable pelvic guidance – padded, incl. cover



Price / piece

Upcharge / piece

Price / piece

for size 0 **(7800017)** incl.
for size 1 **(7801017)** incl.
for size 1b **(7801217)** incl.
for size 2 **(7802017)** incl.

for size 0 **(7800018)** 242 €
for size 1 **(7801018)** 242 €
for size 1b **(7801218)** 242 €
for size 2 **(7802018)** 242 €

for size 0 – seat width A = 20 – 24 cm **(7800020)** 425 €
for size 1+1b – seat width A = 22 – 26 cm **(7801020)** 425 €
for size 2 – seat width A = 28 – 34 cm **(7802020)** 452 €

Seat option


Accessories

Abductable leg guidance

with leg length compensation, padded, incl. cover (already includes thigh guides)

Knee abduction wedge

pluggable using quick-release axle, padded, incl. cover (not usable with abductible leg guidance)



A = Abduction angle
B = Thigh width
C = Leg length compensation

	Width	Depth
Size 0	7 cm	10 cm
Size 1	8 cm	14 cm
Size 2	10 cm	17 cm
Price / piece		
Size 0	(73 00 023)	249 €
Size 1	(73 01 023)	249 €
Size 2	(73 02 023)	249 €

	for size 0	for size 1+1b	for size 2
A	0° – 20°	0° – 20°	0° – 20°
B	7 – 10 cm	9 – 12 cm	10 – 14 cm
C	3 cm	3 cm	3 cm
Price / piece			
for size 0	(78 00 019)		576 €
for size 1	(78 01 019)		576 €
for size 1b	(78 01 219)		576 €
for size 2	(78 02 019)		598 €

Thigh guides



Thigh guides for size 0

padded, incl. cover, adjustable angle, depth and width

Thigh guides for size 1 + 1b

padded, incl. cover, adjustable angle, depth and width

Thigh guides for size 2

padded, incl. cover, adjustable angle, depth and width

	Price / pair			Price / pair			Price / pair	
12 cm	(73 00 114)	213 €	14 cm	(73 01 114)	213 €	17 cm	(73 02 114)	218 €
14 cm	(73 00 115)	213 €	17 cm	(73 01 115)	213 €	22 cm	(73 02 115)	218 €

Footplate

Wadenanlage



	Width	Depth	Lower leg length
for size 0	28,4 cm	17,5 cm	12 – 22 cm* / 15 – 25 cm
for size 1	33,4 cm	20,5 cm	16 – 29 cm* / 19 – 32 cm
for size 1b			19 – 33 cm* / 22 – 36 cm
for size 2	40,3 cm	24,5 cm	25 – 43 cm

* This adjustment range can be achieved by repositioning the bracket.

Footplate standard – without locking mechanism, foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps (Standard configuration, unless otherwise specified)


Footplate with lock – foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps.

Footplate with angle adjustment foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps.

Foot angle = (-)10° – 10°

Calf plate for footplates

covered with black Flexi fabric




Price / piece			Upcharge / piece			Upcharge / piece			Price / piece		
for size 0	(78 00 021)	incl.	for size 0	(78 00 022)	58 €	for size 0	(78 00 023)	106 €	for size 0	(78 00 024)	228 €
for size 1	(78 01 021)	incl.	for size 1	(78 01 022)	58 €	for size 1	(78 01 023)	106 €	for size 1 + 1b	(78 01 024)	228 €
for size 1b	(78 01 221)	incl.	for size 1b	(78 01 222)	58 €	for size 1b	(78 01 223)	106 €	for size 2	(78 02 024)	228 €
for size 2	(78 02 021)	incl.	for size 2	(78 02 022)	58 €	for size 2	(78 02 023)	106 €			

	for size 0	for size 1+1b	for size 2
Width	27 cm	32 cm	38 cm
Height	6,5 cm	6,5 cm	6,5 cm
Price / piece			
for size 0	(78 00 024)		228 €
for size 1 + 1b	(78 01 024)		228 €
for size 2	(78 02 024)		228 €

Configuration possibilities footplates

Step 1: Select the bracket


Knee angle to accommodate the footplates individually adjustable in height, width and angle



Price / piece		
for size 0	(78 00 025)	316 €
for size 1	(78 01 025)	316 €
for size 1b	(78 01 225)	316 €
for size 2	(78 02 025)	316 €

Step 2: Select the footplates


Split footplates – foldable, fixed mounted on the footplate holder, with non-slip surface and each 4 slots for the for the fixation of straps



Foot angle = (-)15° – 15°

	Front width	Rear width	Length
Size 1	11,4 cm	10 cm	18,4 cm
Size 2	13,8 cm	13,2 cm	21,4 cm
Price / pair			
Size 1	(78 01 026)		incl.
Size 2	(78 02 026)		incl.

Foot cups – Foldable, fixed mounted on the footplate holder, with magnetic lock



	min. front width	max. front width	min. rear width	max. rear width	Length
Size 1	8,7 cm	10,4 cm	5,7 cm	7,4 cm	17,4 cm
Size 2	9,5 cm	12 cm	6,7 cm	9,2 cm	20,1 cm
Size 3	11,5 cm	14 cm	8,1 cm	10,6 cm	23,8 cm
Upcharge / pair					
Size 1	(78 01 027)				186 €
Size 2	(78 02 027)				186 €
Size 3	(78 03 027)				186 €

Accessories

Central parking brake for rear castors

incl. 4 dual castors (100 mm), operated by using the brake lever which is mounted to the seat system



4 Castors – 2 with locking devices



4 Castors – all with locking devices



Upcharge / set

Parking lever mounted left	(78 00 030)	467 €
Parking lever mounted right	(78 01 030)	467 €
Parking lever mounted left; castors front with brakes	(78 00 031)	504 €
Parking lever mounted right; castors front with brakes	(78 01 031)	504 €

Upcharge / set

75 mm	(73 00 058)	incl.
100 mm	(73 01 058)	64 €
125 mm	(73 02 058)	101 €

Upcharge / set

75 mm	(73 00 060)	42 €
100 mm	(73 01 060)	101 €
125 mm	(73 02 060)	138 €

Accessories

Therapy tray – made of wood, with rail, flip away



Therapy tray – transparent, with rail, flip away



Therapy tray – made of wood, with rail, prefixable, fixed on both sides



Therapy tray – transparent, with rail, prefixable, fixed on both sides



	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

	Price / piece
for size 0 (78 00 028)	276 €
for size 1 (78 01 028)	276 €
for size 2 (78 02 028)	276 €

	Price / piece
for size 0 (78 00 029)	377 €
for size 1 (78 01 029)	377 €
for size 2 (78 02 029)	377 €

	Price / piece
for size 0 (78 00 038)	294 €
for size 1 (78 01 038)	294 €
for size 2 (78 02 038)	294 €

	Price / piece
for size 0 (78 00 039)	421 €
for size 1 (78 01 039)	421 €
for size 2 (78 02 039)	421 €

Add-on option:

Protection against unintentional folding up	(78 01 036)	108 €
---	-------------	-------

Accessories

Incontinence coating in the seat area



Wheel base extension

inserts at rear, length = 10 cm (not possible together with central parking brake)

Price / piece

73 00 020	217 €
-----------	-------

Price / pair

73 01 028	191 €
-----------	-------

Padding

Padded cover for footplate

covered with black Flexi cover, height: 1cm (not possible in combination with foot straps)

Padded cover for split footplate

covered with black Flexi cover, height: 1cm (not possible in combination with foot straps)

Table upholstery – covered with black Flexi cover, height: 1cm



Price / piece

for size 0	(78 00 032)	133 €
for size 1+1b	(78 01 032)	133 €
for size 2	(78 02 032)	133 €

Price / pair

for size 1	(78 01 033)	191 €
for size 2	(78 02 033)	191 €

Price / piece

for size 0	(73 00 054)	138 €
for size 1	(73 01 054)	138 €
for size 2	(73 02 054)	138 €

Harness assembly kits

Harness assembly kit for the seat

consisting of 2 brackets with clip buckle and 2 clip buckles (drilled)

Harness assembly kit for the back


consisting of 2 threaded plates with clip buckles and 2 clip buckles (drilled)

	Price / set		Price / set
25mm (73 01025)	48 €	25mm (73 01026)	48 €

Harness


Chest-shoulder strap – made of elastic, breathable and moisture repellent material, three straps with clips, sewn strap to the left shoulder; colour: black

	Size 0	Size 1	Size 2	Size 3
A	13cm	13cm	18cm	18cm
B	17,5cm	19cm	20cm	21cm
C	15cm	15cm	19cm	28cm
D	3,5cm	5cm	5cm	5cm
E	25,5cm	28cm	34cm	44cm
Strap width	2cm	2cm	2,5cm	2,5cm



Positioning jacket – with zipper, made of elastic, breathable material, shoulder straps sewn at back, pelvis straps with clips; colour: black

	Size 0	Size 1	Size 2
A	11cm	11cm	18cm
B	21cm	21cm	22cm
C	12cm	12cm	15cm
D	2,5cm	3,5cm	3,5cm
E	25,5cm	28cm	34cm
Strap width	2cm	2cm	2,5cm




	Price / piece		Price / piece
Size 0 (73 00032)	197 €	Size 0 (73 00056)	228 €
Size 1 (73 01032)	197 €	Size 1 (73 01056)	228 €
Size 2 (73 02032)	197 €	Size 2 (73 02056)	228 €
Size 3 (73 03032)	197 €		

Harness


T-shape groin harness – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black

	Size 0	Size 1	Size 2
A	23cm	25cm	28cm
B	15cm	17cm	20cm
C	33cm	35cm	38cm
D	24cm	26cm	28cm
Strap width	2cm	2cm	2,5cm



Groin harness (standard) – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black

	Size 1	Size 2	Size 3
A	24cm	29cm	33cm
B	14cm	17cm	22cm
C	35cm	43cm	49cm
Strap width	2cm	2,5cm	2,5cm



	Price / piece		Price / piece
Size 0 (73 00033)	165 €	Size 1 (73 01034)	165 €
Size 1 (73 01033)	165 €	Size 2 (73 02034)	165 €
Size 2 (73 02033)	165 €	Size 3 (73 03034)	165 €


Harness

Pelvic harness – with pad and adjustable clip on both sides; colour: black



	Size 0	Size 1	Size 2
Pad width	4cm	4cm	5cm
Pad length per side	16cm	18cm	28cm
Strap length per side	30cm	30cm	40cm
Strap width	2cm	2,5cm	2,5cm

4-point pelvic harness – with pad and adjustable clip on both sides; colour: black




	Size 0	Size 1	Size 2
Pad width	4cm	4cm	5cm
Pad length per side	16cm	18cm	28cm
Strap length per side	30cm	30cm	40cm
Strap width	2cm	2,5cm	2,5cm

	Price / piece		Price / piece
Size 0 (73 00035)	99 €	Size 0 (73 00036)	172 €
Size 1 (73 01035)	99 €	Size 1 (73 01036)	172 €
Size 2 (73 02035)	99 €	Size 2 (73 02036)	172 €

Harness


Ankle cuffs – with velcro fastener and idler pulley; colour: black (not possible in combination with padded cover for footplate/split footplate)



	Ankle circumference	Strap width
Size 1	12 – 19cm	2,5cm
Size 2	15 – 23cm	2,5cm
Size 3	19 – 27cm	4cm

	Price / pair
Size 1 (73 01061)	236 €
Size 2 (73 02061)	236 €
Size 3 (73 03061)	236 €

Foot positioner – with straps, incl. velcro fastener, colour: black (not possible in combination with padded cover for footplate/split footplate)



Velcro strip to open

	A	Strap width
Size 1	6,5cm	2,5cm
Size 2	8cm	2,5cm
Size 3	9,5cm	2,5cm

	Price / pair
Size 1 (78 01034)	217 €
Size 2 (78 02034)	217 €
Size 3 (78 03034)	217 €