



malte outdoor.

57

Invoice address:

Dispatch address:

See invoice address

Order Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The compact and versatile posterior-walker for outdoor use.

- Lightweight and folding aluminium frame
- Individual height and width adjustable handles
- Drum brake incl. parking handbrake
- EVA back wheels
- EVA front wheels with direction locks

HMV-No.: 10.46.02.3052




new

Now also
available
in size 4!

Dimensions

		Size 1	Size 2	Size 3	Size 4
A1	Handle height	34 - 58 cm	51 - 75 cm	68 - 94 cm	87 - 113 cm
A2	Height of forearm supports	46 - 64 cm	63 - 87 cm	79 - 103 cm	98 - 124 cm
B	Handle spacing	22 - 36 cm	30 - 44 cm	38 - 52 cm	43 - 57 cm
C	Seat height	33 cm	45 cm	60 cm	70 cm
D	Total width	56 cm	66 cm	73 cm	80 cm
E	Total length	66 cm	75 cm	82 cm	92 cm
	Frame height (folded)	25 cm	32 cm	34 cm	37 cm
	Transport dimensions (L x W x H)	66 x 56 x 25 cm	75 x 66 x 32 cm	82 x 73 x 34 cm	92 x 80 x 37 cm
	Max. load	30 kg	50 kg	70 kg	100 kg
	Weight	7,5 kg	9,9 kg	10,6 kg	11,6 kg



Basic model

new

malte outdoor. size 1 6" tyres (front) and 8" (rear)	malte outdoor. size 2 8" tyres (front) and 10" (rear)	malte outdoor. size 3 8" tyres (front) and 10" (rear)	malte outdoor. size 4 8" tyres (front) and 10" (rear)
Price / piece	Price / piece	Price / piece	Price / piece
75 01 100 1.086 €	75 02 100 1.102 €	75 03 100 1.215 €	75 04 100 1.455 €

Frame colours

Colour selection for size 1	Colour selection for size 2	Colour selection for size 3 + 4	Coated with an individual colour (on request) Please provide the requested colour here: RAL _____
			
Piece	Piece	Piece	Upcharge
Telemagenta (Pink) (92 00 012) incl. Sky blue (92 00 013) incl.	Lasurit blue metallic (92 00 010) incl. Traffic red (92 00 007) incl.	Anthracite metallic (92 00 014) incl. Pure white (92 00 100) incl.	75 00 090 278 €

Accessories

Combined pelvic and buttock support pad – height, width and depth adjustable, PU padded 		4-point pelvic harness – for width adjustable pelvic support pad 				Folding seat – padded and covered 
	Width		Outer part width	Centre section width	Outer part length	Centre section length
Size 1	20 – 28 cm	Size 1	7 cm	7 cm	12 cm	12 cm
Size 2	25 – 33 cm	Size 2	9 cm	8,5 cm	13 cm	14 cm
Size 3	30 – 36 cm	Size 3	11,5 cm	11 cm	15 cm	16 cm
Piece / set		Price / piece				Price / piece
Size 1 (75 01 005)	433 €	Size 1 (75 01 020)	193 €			for size 1 (75 01 007) 344 €
Size 2 (75 02 005)	433 €	Size 2 (75 02 020)	193 €			for size 2 (75 02 007) 344 €
Size 3 (75 03 005)	433 €	Size 3 (75 03 020)	193 €			for size 3 + 4 (75 03 007) 344 €

Accessories

Stops for the hand 	Forearm armrests with handles – adjustable in depth, height and angle, handles adjustable in depth and angle 	Handbrake – with extra small brake handles at reduced grip width with separate locking function 	Storage basket 
	Length	Width	
	Size 1 18 cm	6,5 cm	
	Size 2 25 cm	8 cm	
Piece / set	Piece / set		Upcharge / Set
75 00 004 145 €	Size 1 (75 01 001) 373 €	Size 2 (75 02 001) 373 €	for size 1 (75 01 016) 176 € for size 2 (75 02 016) 176 € for size 3 (75 03 016) 176 € for size 4 (75 04 016) 176 €
			75 00 011 73 €

Accessories

Push bar for attendant 	Anti tip – incl. support roller 	Backwards stop 
Price / piece	Piece / set	Upcharge / pair
for right (75 00 018) 200 € for left (75 01 018) 200 €	75 00 019 145 €	for size 1 (75 01 013) 264 € for size 2 – 4 (75 02 013) 264 €