



**Order
Cost estimate**

Invoice address:

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Dispatch address:

See invoice address

Commission:

The compact and versatile posterior-walker.


- Lightweight and folding aluminium frame (size 0 cannot be folded)
- Individual height and width adjustable handles
- EVA rear wheels with integrated backwards stop function in the hub, can be switched off
- EVA front wheels with direction locks

HMV-No.: 10.46.02.3029



Dimensions

		Size 0 frame can't be folded	Size 1	Size 2	Size 3	Size 4
A1	Handle height	29 - 53 cm	34 - 58 cm	51 - 75 cm	68 - 94 cm	87 - 113 cm
A2	Height of forearm supports	-	46 - 64 cm	63 - 87 cm	79 - 103 cm	98 - 124 cm
B	Handle spacing	22 - 36 cm		30 - 44 cm	38 - 52 cm	43 - 57 cm
C	Seat height	-	33 cm	45 cm	60 cm	70 cm
D	Total width	52 cm	51 cm	61 cm	68 cm	75 cm
E	Total length	58 cm		70 cm	77 cm	87 cm
	Frame height (folded)	37 cm	25 cm	32 cm	34 cm	37 cm
	Transport dimensions (L x W x H)	58 x 52 x 37 cm	58 x 51 x 25 cm	70 x 61 x 32 cm	77 x 68 x 34 cm	87 x 75 x 37 cm
	Max. load	20 kg	30 kg	50 kg	70 kg	100 kg
	Weight	3,9 kg	6,4 kg	7,8 kg	8,4 kg	9,5 kg



Basic model				
malte mini. size 0 Frame colour: Lasurit blue metallic (frame not foldable)	malte. size 1	malte. size 2	malte. size 3	malte. size 4
Price / piece	Price / piece	Price / piece	Price / piece	Price / piece
35 00 100 891 €	75 01 000 923 €	75 02 000 944 €	75 03 000 1.056 €	75 04 000 1.284 €

Frame colours			
Colour selection for size 1 	Colour selection for size 2 	Colour selection for size 3+4 	Coated with an individual colour (on request) Please provide the requested colour here: RAL _____
Piece	Piece	Piece	Upcharge
Telemagenta (Pink) (92 00 012) incl. Sky blue (92 00 013) incl.	Lasurit blue metallic (92 00 010) incl. Traffic red (92 00 007) incl.	Anthracite metallic (92 00 014) incl. Pure white (92 00 100) incl.	75 00 090 278 €

Accessories					
Forearm armrests with handles – adjustable in depth, height and angle, handles adjustable in depth and angle (for size 0 on request)		Forearm armrests with grip balls – adjustable in depth, height and angle, handles adjustable in depth and angle (for size 0 on request)		Grip bar – pivoting, as a cross-bar between the handles (for sizes 0+4 and in combination with adjacent forearm armrests not possible)	
					
Size 1	Length 18 cm Width 6,5 cm	Size 1	Length 18 cm Width 6,5 cm		
Size 2	Length 25 cm Width 8 cm	Size 2	Length 25 cm Width 8 cm		
Piece / set		Piece / set		Price / piece	
Size 1 (75 01 001) 373 € Size 2 (75 02 001) 373 €		Size 1 (75 01 002) 373 € Size 2 (75 02 002) 373 €		for size 1 (75 01 003)	305 €
				for size 2 (75 02 003)	305 €
				for size 3 (75 03 003)	305 €

Accessories			
Stops for the hand	Parking handbrake (not possible for size 0)	Handbrake – with extra small brake handles at reduced grip width with separate locking function (not possible for size 0)	
			
Piece / set	Piece / set	Assembly in combination with a forearm rest (not possible in combination with the grip ball)	
75 00 004 145 €	for size 1 (75 01 008) 264 € for size 2 (75 02 008) 264 € for size 3 (75 03 008) 264 € for size 4 (75 04 008) 264 €	Price / set	
		for size 1 (75 01 015)	441 €
		for size 2 (75 02 015)	441 €
		for size 3 (75 03 015)	441 €
		for size 4 (75 04 015)	441 €

Accessories

Combined pelvic and buttock support pad for sizes 1-4 – height, width and depth adjustable, PU padded



	Width
Size 1	20 – 28 cm
Size 2	25 – 33 cm
Size 3	30 – 36 cm

Combined pelvic and buttock support pad for size 0 – height, width and depth adjustable, PU padded



	Width
for size 0	20 – 28 cm

4-point pelvic harness
for width adjustable pelvic support pad



	Outer part width	Centre section width	Outer part length	Centre section length
Size 1	7 cm	7 cm	12 cm	12 cm
Size 2	9 cm	8,5 cm	13 cm	14 cm
Size 3	11,5 cm	11 cm	15 cm	16 cm

	Piece / set		Piece / set		Price / piece
Size 1 (75 01 005)	433 €	for size 0 (35 01 103)	433 €	Size 1 (75 01 020)	193 €
Size 2 (75 02 005)	433 €			Size 2 (75 02 020)	193 €
Size 3 (75 03 005)	433 €			Size 3 (75 03 020)	193 €

Accessories

Folding seat – padded and covered (not possible for size 0)



Friction brake



Push bar for attendant
(not possible for size 0)



	Price / piece		Piece / set		Price / piece
for size 1 (75 01 007)	344 €	for size 0+1 (75 00 009)	86 €	for rechte (75 00 018)	200 €
for size 2 (75 02 007)	344 €	for size 2-4 (75 01 009)	86 €	for links (75 01 018)	200 €
for size 3+4 (75 03 007)	344 €				

Accessories

Anti tip – incl. support roller
(for size 0 without support roller)



Deflector set
(not possible for size 0)



Bumper wheels – for size 0



Storage basket



	Piece / set		Piece / set		Piece / set		Price / piece
for size 0 (35 00 005)	145 €	75 00 012	213 €	35 00 010	213 €	for size 0 (35 00 011)	73 €
for size 1-4 (75 00 019)	145 €					for size 1-4 (75 00 011)	73 €