



**Order
Cost estimate**

Invoice address:

Dispatch address:

See invoice address

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The anterior-walker for more independence and mobility.


- Lightweight aluminium frame (sizes 2+3 can be folded)
- Individual height and width adjustable handles
- EVA tyres (150 mm for sizes 0+1, 200 mm for sizes 2+3)

HMV-No.: 10.46.02.3049



Dimensions

		Size 0	Size 1	Size 2	Size 3
A1	Handle height	30 – 47 cm	38 – 59 cm	53 – 75 cm	63 – 88 cm
A2	Height of forearm supports	41 – 58 cm	49 – 70 cm	63 – 85 cm	73 – 98 cm
B	Handle spacing	18 – 32 cm	22 – 36 cm	30 – 44 cm	36 – 50 cm
C	Seat height	25 cm	30 cm	45 cm	45 cm
D	Total width	44 cm	48 cm	62 cm	71 cm
E	Total length	36,5 cm	42,5 cm	57,5 cm	61 cm
E	Total length (telescopic wheelbase)	39 – 47 cm	45 – 55 cm	60 – 74 cm	65 – 79 cm
	Max. load	20 kg	30 kg	50 kg	65 kg
	Weight	4,2 kg	5 kg	7,4 kg	9,1 kg



Basic model

marcy. size 0-3

	Price / piece
Size 0 (7600 000)	791 €
Size 1 (7601 000)	791 €
Size 2 (7602 000)	791 €
Size 3 (7603 000)	833 €



**marcy. size 0-3
(telescopic)**

	Price / piece
Size 0 (7600 100)	1.029 €
Size 1 (7601 100)	1.029 €
Size 2 (7602 100)	1.029 €
Size 3 (7603 100)	1.114 €

Basic features
with additional,
adjustable
wheelbase



