



momo tricycle.

Order Cost estimate

Invoice address:

Dispatch address:

See invoice address

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The tricycle as therapy aid for children, teenagers and adults.


- Aluminium frame with extra-low access point
- Rim brake with separate parking brake
- Angle adjustable handlebar stem
- Steering damper for stabilisation when moving straight on (not applicable for steering for attendant)
- Drum brake in front wheel
- incl. guard plates and basket
- Lighting system according to StVZO with hub dynamo (from 24/20")



Dimensions

		12"	16/12"	16"	20"	24/20"	24"	26"	26" (XL)
A*	Pedal to upper edge of saddle	30 - 43 cm	35 - 48 cm	42 - 62 cm	49 - 71 cm	53 - 75 cm	55 - 77 cm	65 - 87 cm	72 - 94 cm
A**	Pedal to upper edge of saddle	33 - 40 cm	38 - 45 cm	45 - 60 cm	52 - 67 cm	56 - 71 cm	58 - 81 cm	68 - 91 cm	75 - 98 cm
B	Saddle support tube to front tube	33 cm	40 cm	39 cm	45 cm	54 cm	52 cm	58 cm	
C	Wheelbase	68 cm	75 cm	80 cm	91 cm	99 cm	104 cm	112 cm	
D	Access point height	16 cm		18 cm	20 cm		22 cm	23 cm	
E	Total length	102 cm	110 cm	120 cm	140 cm	153 cm	165 cm	178 cm	
	Total width	58 cm	65 cm		69 cm		75 cm		
	Max. load	40 kg	50 kg	60 kg	80 kg	90 kg	100 kg	120 kg	
	Weight	12 kg	13,5 kg	14,5 kg	17 kg	18,5 kg	19,5 kg	21,5 kg	21,7 kg

* Dimensions with standard saddle post / ** Dimensions with T-saddle post



Basic model

momo tricycle. 12" (HMV-No.: 22.51.02.0053)	momo tricycle. 16/12" (HMV-No.: 22.51.02.0054)	momo tricycle. 16" (HMV-No.: 22.51.02.0054)	momo tricycle. 20" (HMV-No.: 22.51.02.0055)
Basic features, without drum brake and rear guard plate, additionally with pennant	Basic features, additionally with pennant	Basic features, additionally with pennant	Basic features, additionally with pennant
Price / piece	Price / piece	Price / piece	Price / piece
37 01 000 974 €	37 02 100 1.454 €	37 02 000 1.454 €	37 03 000 1.587 €

Basic model

momo tricycle. 24/20" (HMV-No.: 22.51.02.0051)	momo tricycle. 24" (HMV-No.: 22.51.02.0051)	momo tricycle. 26" (HMV-No.: 22.51.02.0052)	momo tricycle. 26" (XL) (HMV-No.: 22.51.02.0052)
Basic features	Basic features	Basic features	Basic features
Price / piece	Price / piece	Price / piece	Price / piece
37 04 100 1.842 €	37 04 000 1.842 €	37 05 000 1.987 €	37 06 000 2.167 €

Frame colours

Colour selection 12" + 16" + 16/12"	Colour selection 20" + 24/20"	Colour selection 24" + 26" + 26" XL
		
Piece	Piece	Piece
Window grey Pure white (92 00 200) incl. Telemagenta (Pink) Pure white (92 00 201) incl. Turquoise blue Pure white (92 00 202) incl.	Window grey Pure white (92 00 200) incl. Pastel blue Pure white (92 00 203) incl. Kiwisorbet green Pure white (92 00 204) incl.	Window grey Pure white (92 00 200) incl. Ivory Pure white (92 00 205) incl. Black matt Pure white (92 00 206) incl.

Configuration possibilities - drive

	Art.-No.	Drive	Upcharge / piece
	37 00 001	Rigid sprocket (standard configuration, if nothing else stated)	incl.
	37 00 004	Freewheel brake hub (with back pedal brake)	192 €
	37 00 503	Change from rigid sprocket to neutral hub	437 €
	37 00 002	Extra light transmission (only for above drive types)	incl.
	37 00 603	Switchover from rigid sprocket, neutral hub and 3-gears-freewheel brake hub with reverse option	683 €
	37 00 005	3-gear freewheel brake hub (with back pedal brake)	253 €
	37 00 007	3-gear freewheel brake hub with reverse option	424 €
	37 00 008	7-gear freewheel brake hub (with back pedal brake)	479 €
	37 00 202	7-gear-freewheel hub with roll brake (brake lever mounted on the handlebar, without back pedal brake)	604 €

Saddle

Standard saddle (standard configuration, if nothing else stated)				Gel saddle			
							
	Width		Length		Rear width	Front width	Length
Size 1	15 cm		21 cm	Size 1	19 cm	4 cm	24 cm
Size 2	15 cm		24 cm	Size 2	24 cm	7 cm	27 cm
Size 3	18 cm		26 cm				
Price / piece				Upcharge / piece			
Size 1 (37 01 001)			incl.	Size 1 (37 01 024)			64 €
Size 2 (37 02 001)			incl.	Size 2 (37 02 024)			64 €
Size 3 (37 03 001)			incl.				

Saddle

Unicycle saddle – banana shape (caution: min. inside leg length increases by 2 cm)



Rear width	Middle width	Front width	Length
11 cm	6 cm	8 cm	25 cm

Upcharge / piece

37 00 010 119 €

Moped saddle – (caution: min. inside leg length increases by 2 cm)



Rear width	Front width	Rear length	Length
26 cm	9,5 cm	12 cm	25 cm

Upcharge / piece

37 00 023 176 €

Saddle supports

Standard saddle support – (standard configuration, if nothing else stated)



	for
Size 1	12" + 16"/12"
Size 2	16"
Size 3	20" – 26"

Price / piece

Size 1	(37 01 004)	incl.
Size 2	(37 02 004)	incl.
Size 3	(37 03 004)	incl.

T-saddle support – for horizontal saddle adjustment



	for
Size 1	12" + 16"/12"
Size 2	16" + 20"
Size 3	24" + 24"/20" + 26"


Upcharge / piece

Size 1	(37 01 005)	93 €
Size 2	(37 02 005)	93 €
Size 3	(37 03 005)	93 €

Holding brackets

Holding brackets incl. holder – only for use together with the T-saddle support

	max. Height*	Depth
Size 1	20 cm	12 cm
Size 2	30 cm	12 cm
Size 3	37 cm	12 cm
Size 4	53 cm	12 cm
Tailored	____ cm (max. 85 cm)	____ cm (max. 30 cm)



* Max. height: measured from the standard saddle to the top edge of the respective support pad

Headrest

Headrest – only for use together with a holding bracket



	Width	Height
Size 1	20 cm	15 cm
Size 2	23 cm	18 cm

Price / set

Price / piece

Size 1	(37 01 007)	143 €	Size 3	(37 03 007)	143 €
Size 2	(37 02 007)	143 €	Size 4	(37 04 007)	143 €
			Tailored	(37 09 007)	304 €

Size 1	(37 01 029)	384 €
Size 2	(37 02 029)	384 €

Torso supports

Dynamic back support pad padded with strap and bracket (only for use together with a holding bracket)



	Width
Size 1	20 – 30 cm
Size 2	25 – 35 cm
Size 3	30 – 40 cm

Price / piece

Size 1	(37 01 008)	271 €
Size 2	(37 02 008)	271 €
Size 3	(37 03 008)	271 €

Width adjustable back support pad padded, incl. bracket (only for use together with a holding bracket)



	Width
Size 1	20 – 28 cm
Size 2	25 – 33 cm
Size 3	30 – 36 cm

Price / piece

Size 1	(37 01 025)	362 €
Size 2	(37 02 025)	362 €
Size 3	(37 03 025)	362 €

Pelvic supports

Dynamic pelvic guide support pad padded with strap and bracket (only for use together with a holding bracket)



	Width
Size 1	20 – 30 cm
Size 2	25 – 35 cm
Size 3	30 – 40 cm

Price / piece

Size 1	(37 01 009)	271 €
Size 2	(37 02 009)	271 €
Size 3	(37 03 009)	271 €

Width adjustable pelvic guide support pad – padded, incl. bracket (only for use together with a holding bracket)



	Width
Size 1	20 – 28 cm
Size 2	25 – 33 cm
Size 3	30 – 36 cm

Price / piece

Size 1	(37 01 026)	362 €
Size 2	(37 02 026)	362 €
Size 3	(37 03 026)	362 €

Harness system

Chest strap – for width adjustable back support pad



	Width	Length
Size 1	5 cm	18,5 cm
Size 2	6,5 cm	25,5 cm
Size 3	7 cm	30,5 cm

Price / piece

Size 1	(37 01 050)	109 €
Size 2	(37 02 050)	109 €
Size 3	(37 03 050)	109 €

Positioning jacket – for width adjustable back support pad, incl. universal bracket



	Inside width	Total width	Side height	Total height
Size 1	8 cm	25 cm	10 cm	28 cm
Size 2	9 cm	30 cm	12 cm	30 cm
Size 3	12 cm	34 cm	16 cm	35 cm

Price / piece

Size 1	(37 01 051)	271 €
Size 2	(37 02 051)	271 €
Size 3	(37 03 051)	271 €

Harness system

T-shape groin harness – for width adjustable pelvic support pad

	Seat width	Front width	Seat length	Total length
Size 1	24 cm	23 cm	15 cm	32 cm
Size 2	26 cm	26 cm	17 cm	34 cm
Size 3	28 cm	28 cm	19 cm	38 cm



4-point pelvic harness – for width adjustable pelvic support pad

	Outer part width	Centre section width	Outer part length	Centre section length
Size 1	7 cm	7 cm	12 cm	12 cm
Size 2	9 cm	8,5 cm	13 cm	14 cm
Size 3	11,5 cm	11 cm	15 cm	16 cm



Piece / Stück

Price / piece

Size 1 **(37 01 052)**
Size 2 **(37 02 052)**
Size 3 **(37 03 052)**

186 €
186 €
186 €

Size 1 **(37 01 053)**
Size 2 **(37 02 053)**
Size 3 **(37 03 053)**

191 €
191 €
191 €

Accessories

Universal bracket – for harness used without support pads (width = 27,5 cm)



Price / piece

37 02 055

124 €

Hand positioner – for all handlebars

	Front width
Size 0	8,5 cm
Size 1	10 cm
Size 2	12 cm
Size 3	14 cm



Piece / pair

Size 0 **(37 00 057)**
Size 1 **(37 01 057)**
Size 2 **(37 02 057)**
Size 3 **(37 03 057)**

156 €
156 €
156 €
156 €

Steering stop – adjustable up to direction lock (incl. for steering for attendant)



	for
Size 1	12"
Size 2	16/12" – 26"

Price / piece

Size 1 **(37 01 006)**
Size 2 **(37 03 006)**

143 €
143 €

Handlebar extension



	for	Length
Size 1	12"	10 cm
Size 2	16/12" – 26"	10 cm

Price / piece

Size 1 **(37 01 022)**
Size 2 **(37 02 022)**

32 €
32 €

Handlebar

Classic handlebar – black (standard configuration, if nothing else stated); When selecting the 7-speed gear system, available only from size 3 onwards.



	Standard for	Width
Size 1	12"	47 cm
Size 2	16/12" + 16"	49 cm
Size 3	20"	52 cm
Size 4	24/20" – 26"	58 cm

Price / piece

Size 1 **(37 01 011)** incl.
Size 2 **(37 02 011)** incl.
Size 3 **(37 03 011)** incl.
Size 4 **(37 04 011)** incl.

Touring handlebar



	Width	Depth
Size 0	50 cm	17 cm
Size 1	55 cm	17 cm
Size 2	58 cm	17 cm

Upcharge / piece

Size 0 **(37 00 112)** 125 €
Size 1 **(37 01 012)** 125 €
Size 2 **(37 02 012)** 125 €

Round handlebar



	Width	Depth
Size 1	40 cm	17 cm
Size 2	43 cm	21 cm
Size 3	48 cm	25 cm

Upcharge / piece

Size 1 **(37 01 013)** 125 €
Size 2 **(37 02 013)** 125 €
Size 3 **(37 03 013)** 125 €

Multifunctional handlebar suitable for 20" – 26"



	Width	Depth
Size 1	61,5 cm	16,5 cm

Upcharge / piece

37 02 014 181 €

Handlebar brake

Round handlebar with brake lever ring



	Width	Depth
Size 1	40 cm	17 cm
Size 2	43 cm	21 cm
Size 3	48 cm	25 cm

Upcharge / piece

Size 1 **(37 01 028)** 367 €
Size 2 **(37 02 028)** 367 €
Size 3 **(37 03 028)** 367 €

Round handlebar with brake function – the brake function is simply activated by pressing down the round handlebars



	Width	Depth
Size 1	40 cm	17 cm
Size 2	43 cm	21 cm
Size 3	48 cm	25 cm

Upcharge / piece

Size 1 **(37 01 031)** 367 €
Size 2 **(37 02 031)** 367 €
Size 3 **(37 03 031)** 367 €

Accessories

Push bar – height adjustable and removable with integrated twist lock



Price / piece

37 00 011

200 €

Parking brake for attendant
(can only be used in connection with a push bar)



Price / piece

37 00 012

133 €

Steering for attendant – via deflection bar, incl. steering stop, in this case the steering damper is omitted (can only be used in connection with a push bar)



Price / piece

for 12"	(37 01 010)	484 €
for 16/12"	(37 02 110)	484 €
for 16"	(37 02 010)	484 €
for 20"	(37 03 010)	484 €
for 24/20"	(37 04 110)	484 €
for 24"	(37 04 010)	484 €
for 26"	(37 05 010)	484 €

Crank

Progressively adjustable crank reducer

for 12"-20" there is no widening of the pedal spacing and for 24"+26" there is a widening of the pedal spacing of approx. 18 mm on each side



Piece / pair

for 12"-20"+16/12"+24/20"
for 24"+26"

(37 00 016)
(37 04 016)

197 €
197 €

Crank reducer – minimal pedal distance widening - reduction in size by 2.5 or 5 cm respectively



Piece / pair

for 16"+20"+24/20"
for 24"+26"

(37 02 015)
(37 05 015)

102 €
102 €

Special crank for knee contracture – for 20" - 26"



Price / piece

mounted on right
mounted on left

(37 03 015)
(37 04 015)

422 €
422 €

Pedals

Foot fixing pedals – with racing hooks and compensation weights



Depth

13 cm


Home exerciser pedals



Width

12 cm

Ergotec pedals – with foot positioning



Velcro strip to open

	A	Pedal width
Size 1	6,5 cm	9 cm
Size 2	8 cm	9 cm



Upcharge / pair

37 00 013

102 €

Upcharge / pair

37 00 014

69 €

Upcharge / pair

Size 1 **(37 01 056)**
Size 2 **(37 02 056)**

191 €
191 €

Foot cups

Foot cups – width adjustable, with fixing straps & with magnetic lock



Accessories

Walking stick holder



Price / piece

37 00 040

132 €

	Min. front width	Max. front width	Min. rear width	Max. rear width	Length
Size 1	8,7 cm	10,4 cm	5,7 cm	7,4 cm	17,4 cm
Size 2	9,5 cm	12 cm	6,7 cm	9,2 cm	20,1 cm
Size 3	11,5 cm	14 cm	8,1 cm	10,6 cm	23,8 cm
Piece / pair					
Size 1 (37 01 037)	238 €				
Size 2 (37 02 037)	238 €				
Size 3 (37 03 037)	238 €				

Foot cups with leg guidance


Foot cups with leg guidance – width adjustable and with fixing straps & with magnetic lock

	Min. front width	Max. front width	Min. rear width	Max. rear width	Length	Min. leg guide height	Max. leg guide height
Size 1	8,7 cm	10,4 cm	5,7 cm	7,4 cm	17,4 cm	15 cm	18,5 cm
Size 2	9,5 cm	12 cm	6,7 cm	9,2 cm	20,1 cm	18 cm	22 cm
Size 3	11,5 cm	14 cm	8,1 cm	10,6 cm	23,8 cm	21 cm	26 cm
Size 4	11,5 cm	14 cm	8,1 cm	10,6 cm	23,8 cm	28 cm	36 cm
							Piece / pair
Size 1 (37 01 038)							404 €
Size 2 (37 02 038)							404 €
Size 3 (37 03 038)							404 €
Size 4 (37 04 038)							404 €



Foot cups with dynamic leg guidance

Foot cups with dynamic leg guidance
width adjustable and with fixing straps
& with magnetic lock



	Min. front width	Max. front width	Min. rear width	Max. rear width	Length	Min. leg guide height	Max. leg guide height
Size 1	8,7 cm	10,4 cm	5,7 cm	7,4 cm	17,4 cm	15 cm	18,5 cm
Size 2	9,5 cm	12 cm	6,7 cm	9,2 cm	20,1 cm	18 cm	22 cm
Size 3	12 cm	15 cm	8,8 cm	11,8 cm	24 cm	20,5 cm	25,5 cm
							Piece / pair
Size 1 (37 01 039)							528 €
Size 2 (37 02 039)							528 €
Size 3 (37 03 039)							528 €

Configuration possibilities - lighting system

	Art.-No.	Lighting system	for	Upcharge / Set
	3700021	Upgrading of a battery light system attachable incl. reflectors according to StVZO	12" - 26"	189 €
	3701021	Upgrade with a lighting system with hub dynamo and reflectors according to StVZO	16/12"+ 16"+20"	189 €

momo motion. - The tricycle with E-drive

A motor support when using a tricycle is useful for many users and necessary due to the indication. Thus, for some users, the starting situation is often the only obstacle to independent use of the bike. With clinical pictures such as dystrophies, permanent support is required to compensate the relative loss of muscle strength.

In order to make the medical benefits of riding a tricycle accessible to these people, we have developed our **momo motion.**

The corresponding order form for our **momo motion.** as well as further information on the possibilities of motor support can be found on page 94 or on our homepage at www.schuchmann.de/en/momo-motion.

