



Invoice address:

Dispatch address: See invoice address

Order
Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____


Commission:

- The pivoting car seat.**
- Seat and back plate made of aluminium with back angle adjustment
 - Removable pads and fixed lateral pads
 - Incl. height adjustable and padded headrest and 5-point brace straps
 - Three-point strap guidance, frame colour: black – silky matt, RAL 9005
 - Tested in accordance with ECE-44/03; for use in accordance with the third ordinance on exemptions from road traffic regulations

HMV-NR.: 26.11.06.0009



Dimensions							
		Size 1		Size 2		Size 3	
		a	b	a	b	a	b
A	Seat depth	25cm	30cm		35cm		40cm
B	Seat width	28cm		32cm		36cm	
C	Back height	40cm		48cm		56cm	
D	Headrests height	59 - 69cm		65 - 77cm		76 - 86cm	
E	Lateral pad height	25 - 33cm		32 - 40cm		36 - 44cm	
	Tilt-in-space*	0° / 15°			0° / 10°		
F	Chest width	14 - 25cm		16 - 29cm		16 - 32cm	
G	Back angle			90° - 108°			
	Max. load	30 kg		45 kg		60 kg	
* when using the footplate and tilt adjustment							




Basic model					
sascha. size 1		sascha. size 2		sascha. size 3	
Price / piece		Price / piece		Price / piece	
Size 1a (14 01 001)	2.189 €	Size 2a (14 02 001)	2.273 €	Size 3a (14 03 001)	2.320 €
Size 1b (14 01 100)	2.189 €	Size 2b (14 02 100)	2.273 €	Size 3b (14 03 100)	2.320 €

Covers					
Fabric design Black / Red (No. 105 + 91)		Fabric design Black / Light blue (No. 105 + 106)		Fabric design Black / Spacer fabric (No. 105 + 47)	
					
Piece		Piece		Piece	
9100206		9100230		9100231	
incl.		incl.		incl.	

Lateral pads					
Standard setting range for fixed lateral pads The max. chest width reduces by 2 cm on each side when using the pluggable lateral pads For size 1 a / b Height from seat to armpit 25 – 33 cm Width of chest area 14 – 25 cm For size 2 a / b Height from seat to armpit 32 – 40 cm Width of chest area 16 – 29 cm For size 3 a / b Height from seat to armpit 36 – 44 cm Width of chest area 16 – 32 cm		Setting range for tailored lateral pads outside of the adjacent dimensions Width of chest area _____ cm Width of chest area outside of the adjacent dimensions on request.		Lateral pads – pluggable, easier to get in and out, incl. cover 	
				Upcharge / piece	
				for size 1, mounted on right (14 01 007)	248 €
				for size 1, mounted on left (14 01 008)	248 €
				for size 2, mounted on right (14 02 007)	248 €
				for size 2, mounted on left (14 02 008)	248 €
				for size 3, mounted on right (14 03 007)	248 €
				for size 3, mounted on left (14 03 008)	248 €
				Upcharge	
				14 00 009	253 €

Strap guide		Accessories			
Shoulder strap guide 14 00 002: Shoulder strap guide mounted on right (standard configuration, if nothing else stated) 14 00 003: Shoulder strap guide mounted on left 14 00 004: Shoulder strap guide mounted on left		Lateral guide and torso support pad as safety impact protection, incl. cover and fixing strap 		Abduction wedge – pluggable, with bracket and cover 	
				ISOFIX bracket 	
Price / piece		Price / piece		Price / piece	
right (14 00 002)	incl.	for size 1 (14 01 005)	807 €	for size 1 (14 01 006)	366 €
left (14 00 003)	incl.	for size 2 (14 02 005)	807 €	for size 2 (14 02 006)	366 €
both sides (14 00 004)	100 €	for size 3 (14 03 005)	807 €	for size 3 (14 03 006)	366 €
				Piece / set	
				for size 1 (14 01 010)	260 €
				for size 2 (14 02 010)	260 €
				for size 3 (14 03 010)	260 €

Pivoting device					
Pivoting device with tilt adjustment incl. footplate support and retainer bar  		Right pivoting		Left pivoting	
		Piece / set		Piece / set	
		X = 394 mm for size 1 (14 01 002)		X = 394 mm for size 1 (14 01 102)	
		502 €		502 €	
		X = 454 mm for size 2 (14 02 002)		X = 454 mm for size 2 (14 02 102)	
		502 €		502 €	
		X = 494 mm for size 3 (14 03 002)		X = 494 mm for size 3 (14 03 102)	
		502 €		502 €	

Accessories												
Tilt adjustment on base plate incl. footplate support		Footplate – removable (only for use in combination with the pivoting device or tilt adjustment)		Footplate – for installation in the footwell, adjustable in height and angle								
												
				<table><tr><th>Width</th><th>Height</th><th>Depth</th></tr><tr><td>36 cm</td><td>19–31 cm</td><td>26 cm</td></tr></table>	Width	Height	Depth	36 cm	19–31 cm	26 cm		
Width	Height	Depth										
36 cm	19–31 cm	26 cm										
Price / piece		Price / piece		Price / piece								
for size 1 (1401003)	188 €	for size 1 (1401004)	177 €	1400012 348 €								
for size 2 (1402003)	188 €	for size 2 (1402004)	177 €									
for size 3 (1403003)	188 €	for size 3 (1403004)	177 €									