



Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The convenient seating aid for everyday use.


- Sturdy wooden frame with integrated aluminium profiles
- Incl. two castors
- Incl. footrest
- Incl. push handle
- Incl. seat and back pad
- Height and angle adjustable backrest
- Height adjustable armrests
- Material: varnished beech
- **Base model with quick release: easy and fast adjustable without tools:**
 - height-, angle- and depth-adjustable seat
 - height- and depth adjustable footplate



HMV-No.: 26.11.05.0034

Dimensions

			Dimensions				
			Size 0 Nursery	Size 0	Size 1 Nursery	Size 1	Size 2
A	Seat width		17 – 28 cm		20 – 35 cm		30 – 40 cm
B	Seat depth		16 – 26 cm		25 – 35 cm		30 – 40 cm
C	Max. seat height		42 cm	62 cm	42 cm	62 cm	
D*	Back height		25 cm		29 cm		34 cm
E	Back angle		(-) 3° – 17,5°				
F	Seat angle		(-) 5° – 15°				
G	Footrest angle		(-) 5° – 15°				
H**	Footrest quick height adjustment min.	Standard	12 cm				
		with 50 mm dual castors	17 cm				
		with 75 mm dual castors	20 cm				
H**	Height of the footrest at floor level min.	Standard	2,5 cm				
		with 50 mm dual castors	7,5 cm				
		with 75 mm dual castors	10,5 cm				
I	Lower leg height with quick height-adjustable footrest min.		14 cm				
I	Lower leg height with footrest low to the ground min.		18,6 cm				
	Total width		46 cm	47 cm	53 cm	54 cm	59 cm
	Length base frame		63 cm	69 cm***	63 cm	69 cm***	69 cm***
	Max. load		30 kg		45 kg		60 kg
	Weight		9,2 kg	9,7 kg	9,6 kg	10,2 kg	12 kg



* Minimum back height measured from seat surface to upper edge of back = D + 6cm** adjustable lower leg length max. = selected/required seat height - footrest height min. (H)
***long stand + 10 cm

Basic model <u>WITHOUT</u> quick-release heights adjustment					
smilla. size 0		smilla. size 1		smilla. size 2	
Price / piece		Price / piece		Price / piece	
Size 0 Nursery (13 00 500)	1.310 €	Size 1 Nursery (13 01 500)	1.310 €	Size 2 (13 02 400)	1.530 €
Size 0 (13 00 400)	1.310 €	Size 1 (13 01 400)	1.310 €		



Basic model <u>WITH</u> quick-release heights adjustment					
smilla. size 0		smilla. size 1		smilla. size 2	
Price / piece		Price / piece		Price / piece	
Size 0 Nursery (13 00 300)	1.526 €	Size 1 Nursery (13 01 300)	1.526 €	Size 2 (13 02 200)	1.748 €
Size 0 (13 00 200)	1.526 €	Size 1 (13 01 200)	1.526 €		

Frame colours		
Beech natural	Beech white	Beech slate grey
		
Piece	Piece	Piece
92 00 300 incl.	92 00 301 incl.	92 00 302 incl.

Covers			
Pink 2 (No. 115)	Deep sea blue (No. 111)	Leaf green (No. 123)	Slate grey (No. 124)
			
Price	Price	Price	Price
91 00 115 incl.	91 00 111 incl.	91 00 123 incl.	91 00 124 incl.

Covers			
Flexi blue (No. 44)	Flexi light grey (No. 50)	Cover according to customer request	Name embroidery on the back pad please fill out the form on page 137
			
Upcharge	Upcharge	Upcharge	An embroidery of Flexi fabric is not possible.
13 00 044 272 €	13 00 050 272 €	13 00 999 278 €	

Accessories						
Headrest – padded and height adjustable		Lateral pads – padded and adjustable		Thigh guides – padded and adjustable		
						
Adjustable height (in 3 cm steps)		Width 1	Width 2	Length	Height	
for size 0+1	12 – 21cm	for size 0	13 – 19cm 17 – 23cm	Size 0	17,5cm	7cm
for size 1		for size 1	20 – 26cm 24 – 30cm	Size 1	20 cm	8,5cm
for size 2	20,5 – 32,5cm	for size 2	24 – 30cm 28 – 34cm	Size 2	22,5cm	10cm
Price / piece		Price / pair		Price / pair		
for size 0 (13 00 101)	301 €	Width 1 (13 01 002)	192 €	Size 0 (13 00 004)	259 €	
for size 1 (13 01 101)	301 €	Width 2 (13 01 003)	192 €	Size 1 (13 01 004)	259 €	
for size 2 (13 02 101)	301 €			Size 2 (13 02 004)	259 €	

Configuration options - seat			
Standard seat pad – incl. cover (standard configuration, if nothing else stated)		Anatomical seat pad – incl. cover with anatomic leg guidance	
			
Price / piece		Upcharge / piece	
for size 0 (13 00 021)	incl.	for size 0 (13 00 022)	243 €
for size 1 (13 01 021)	incl.	for size 1 (13 01 022)	243 €
for size 2 (13 02 021)	incl.	for size 2 (13 02 022)	243 €

Accessories

Wooden therapy tray – plug-on



	Width	Depth
for size 0	48cm	37cm
for size 1	52cm	39cm
for size 2	57cm	39cm
Price / piece		
for size 0 (1300007)		288 €
for size 1 (1301007)		288 €
for size 2 (1302007)		288 €

Large Wooden therapy tray – plug-on, with rail



	Width	Depth
for size 0	48cm	45cm
for size 1	52cm	49cm
for size 2	57cm	49cm
Price / piece		
for size 0 (1300008)		320 €
for size 1 (1301008)		320 €
for size 2 (1302008)		320 €

Abduction wedge – padded and attachable at the front



Price / piece	
for size 0 (1900008)	252 €
for size 1+2 (1901008)	252 €

Footrest for Basic model WITHOUT quick-release heights adjustment

Footrest height quickly adjustable

height and depth adjustable, can be lowered to 12 cm above the floor



Upcharge / piece	
for size 0 (1300041)	89 €
for size 1 (1301041)	89 €
for size 2 (1302041)	89 €

Footrest angle and height quickly adjustable

angle, height and depth adjustable, can be lowered up to 12 cm above the floor



Upcharge / piece	
for size 0 (1300042)	205 €
for size 1 (1301042)	205 €
for size 2 (1302042)	205 €

Footrest floor level – height and depth adjustable, can be lowered to the floor



Price / piece	
for size 0 (1300040)	incl.
for size 1 (1301040)	incl.
for size 2 (1302040)	incl.

No footrest

Price reduction	
1300037	-220 €

Footrest for Basic model WITH quick-release heights adjustment

Footrest height quickly adjustable

height and depth adjustable, can be lowered to 12 cm above the floor



Price / piece	
for size 0 (1300032)	incl.
for size 1 (1301032)	incl.
for size 2 (1302032)	incl.

Footrest angle and height quickly adjustable

angle, height and depth adjustable, can be lowered up to 12 cm above the floor



Upcharge / piece	
for size 0 (1300038)	124 €
for size 1 (1301038)	124 €
for size 2 (1302038)	124 €

Footrest floor level – height and depth adjustable, can be lowered to the floor



Price / piece	
for size 0 (1300031)	w/o upcharge
for size 1 (1301031)	w/o upcharge
for size 2 (1302031)	w/o upcharge

No footrest

Price reduction	
1300033	-220 €

Seat

Seat without quick height adjustment

(standard for base model w/o quick-release heights adjustment)



Price / piece	
for size 0 (1300035)	incl.
for size 1 (1301035)	incl.
for size 2 (1302035)	incl.


Quick height adjustment seat

(upcharge for base model w/o quick-release heights adjustment)



Upcharge / set	
for size 0 (1300036)	136 €
for size 1 (1301036)	136 €
for size 2 (1302036)	136 €

Pedestal – when the child is able to transfer independently over the footrest into the smilla., we recommend the longer pedestal for size 1+2 (not possible for the KITA version!).



Upcharge / set	
Standard for size 0-2	(1300034) incl.
Long (10 cm longer) for size 0-2	(1301034) 108 €

Accessories



Footrest pad – covered with black Flexi fabric, height: 1 cm



Footrest pad – covered with black Flexi fabric, height: 3 cm



Preis / Stück	
for size 0 (1300046)	133 €
for size 1 (1301046)	133 €
for size 2 (1302046)	133 €

Preis / Stück	
for size 0 (1300045)	165 €

Castors		Harness	
Dual castors – 4 pieces 		Harness assembly kit for the seat – consisting of 2 brackets with 2 clip buckles 	
Harness assembly kit for the back – consisting of belt holder with T-groove incl. 2 clip buckles and 2 clip buckles (drilled) 			
Price / set		Price / set	Price / set
50 mm – 2 with locking devices	(13 01 006) 320 €	25 mm (13 01 023)	48 €
50 mm – all with locking devices	(13 01 016) 360 €		
75 mm – 2 with locking devices	(13 00 006) 320 €		
75 mm – all with locking devices	(13 00 016) 360 €		
		25 mm (13 01 024)	164 €

Harness				
Chest-shoulder strap – made of elastic, breathable and moisture repellent material, three straps with clips, sewn strap to the left shoulder, colour: black				
	Size 0	Size 1	Size 2	Size 3
A	13 cm	13 cm	18 cm	18 cm
B	17.5 cm	19 cm	20 cm	21 cm
C	15 cm	15 cm	19 cm	28 cm
D	3.5 cm	5 cm	5 cm	5 cm
E	25.5 cm	28 cm	34 cm	44 cm
Strap width	2 cm	2 cm	2.5 cm	2.5 cm
				
Price / piece				
Size 0 (13 00 025) 197 €				
Size 1 (13 01 025) 197 €				
Size 2 (13 02 025) 197 €				
Size 3 (13 03 025) 197 €				
Positioning jacket – with zipper, made of elastic, breathable material, shoulder straps sewn at back, pelvis straps with clips, colour: black				
	Size 0	Size 1	Size 2	
A	11 cm	11 cm	18 cm	
B	21 cm	21 cm	22 cm	
C	12 cm	12 cm	15 cm	
D	2.5 cm	3.5 cm	3.5 cm	
E	25.5 cm	28 cm	34 cm	
Strap width	2 cm	2 cm	2.5 cm	
				
Price / piece				
Size 0 (13 00 026) 228 €				
Size 1 (13 01 026) 228 €				
Size 2 (13 02 026) 228 €				

Harness				
T-shape groin harness – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back, colour: black				
	Size 0	Size 1	Size 2	
A	23 cm	25 cm	28 cm	
B	15 cm	17 cm	20 cm	
C	33 cm	35 cm	38 cm	
D	24 cm	26 cm	28 cm	
Strap width	2 cm	2 cm	2.5 cm	
				
Price / piece				
Size 0 (13 00 027) 165 €				
Size 1 (13 01 027) 165 €				
Size 2 (13 02 027) 165 €				
Groin harness (standard) – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back, colour: black				
	Size 1	Size 2	Size 3	
A	24 cm	29 cm	33 cm	
B	14 cm	17 cm	22 cm	
C	35 cm	43 cm	49 cm	
Strap width	2 cm	2.5 cm	2.5 cm	
				
Price / piece				
Size 1 (13 01 028) 165 €				
Size 2 (13 02 028) 165 €				
Size 3 (13 03 028) 165 €				

Harness				
Pelvic harness – with pad and adjustable clip on both sides; colour: black				
				
	Size 0	Size 1	Size 2	
Pad width	4 cm	4 cm	5 cm	
Pad length per side	16 cm	18 cm	28 cm	
Strap length per side	30 cm	30 cm	40 cm	
Strap width	2 cm	2.5 cm	2.5 cm	
Price / piece				
Size 0 (13 00 029) 99 €				
Size 1 (13 01 029) 99 €				
Size 2 (13 02 029) 99 €				
4-point pelvic harness – with pad and adjustable clip on both sides; colour: black				
				
	Size 0	Size 1	Size 2	
Pad width	4 cm	4 cm	5 cm	
Pad length per side	16 cm	18 cm	28 cm	
Strap length per side	30 cm	30 cm	40 cm	
Strap width	2 cm	2.5 cm	2.5 cm	
Price / piece				
Size 0 (13 00 030) 172 €				
Size 1 (13 01 030) 172 €				
Size 2 (13 02 030) 172 €				