



Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:


The height adjustable supine standing frame.

- Height and angle adjustment (0-90°) via gas pressure spring for sizes 1+2
- Height and angle adjustment (0-90°) via electric power for sizes 3+4 (incl. safety shut-off and manual operation)
- Frame shaped to work with mobile hoists
- Upper body support adjustable in length, incl. padding elements to adjust the length of the upper body support
- Incl. chest strap and pelvic harness
- Adapter for leg rests (standard)
- Leg rests (standard), adjustable in length, with up to 5 cm pelvic rotation adjustment
- Split footplates (standard), which can be lowered to the floor via height adjustment of the center column
- Base frame with 5 (for size 1-3) or 6 (for size 4) castors (75 mm) with total locking devices, one of them with directional lock
- Frame colour: Telegrey 4 (white grey), attachment parts: black
- Incl. tool for all settings (adjustments can be made without tools)

HMV-No.: 28.29.02.1057




Dimensions



			Size 1	Size 2	Size 3	Size 4
A	Body size		75 - 115cm	100 - 135 cm	120 - 160 cm	145 - 185 cm
B	Upper body support	Hip/thorax width	14 - 27 cm	17 - 31cm	20 - 37 cm	27 - 46 cm
C		Length	30,5 - 43 cm	38 - 48 cm	49 - 61 cm	58 - 69 cm
D	Leg rests	Foot to knee	20,5 - 30,5 cm	27,5 - 37 cm	32,5 - 45 cm	42 - 54 cm
E		Inside leg to knee	12 - 18,5 cm	16,5 - 24 cm	20 - 29,5 cm	24,5 - 34 cm
		Knee angle	0° - 65°			
		Hip angle	0° - 45°			
		Abduction (each leg)	0° - 30°			
F		Pelvic rotation adjustment	0 - 5 cm			
G	Feet	Dorsiflexion	0° - 20°			
G		Plantarflexion	0° - 15°			
H		Inversion/Eversion	0° - 15° each direction			
I	Height of lying surface (for transfer)		53,5 - 73,5 cm	60 - 85 cm	71,5 - 93,5 cm	82,5 - 104,5 cm
J	Tilt		0° - 90°			
	Transport dimensions (L x W x H)		95 x 51 x 75 cm	105 x 61 x 85 cm	125 x 76 x 97 cm	152 x 77 x 110 cm
	Max. load		30 kg	60 kg	80 kg	100 kg
	Weight		30 kg	35 kg	80 kg	90 kg

Settings



The following measurements should be pre-set on the standing frame:

I Body size	_____	VI Pelvic width	_____
II Axilla height (Centre lateral pads)	_____	VII Hip height (Centre pelvic pads)	_____
III Shoulder height	_____	VIII Knee height	_____
IV Inside leg height	_____	IX Neck width	_____
V Chest width	_____	X Tray height	_____

42 00 020

Piece

incl.

Basic model			
till. size 1 for users with a body size of 75 – 115 cm		till. size 2 for users with a body size of 100 – 135 cm	
with gas spring supported height and angle adjustment	with gas spring supported angle adjustment and power supported height adjustment (mains supply)	with gas spring supported height and angle adjustment	with gas spring supported angle adjustment and power supported height adjustment (mains supply)
Price / piece	Price / piece	Price / piece	Price / piece
42 01 000	2.829 €	42 01 100	4.757 €
		42 02 000	3.049 €
		42 02 100	4.989 €

Basic model	
till. size 3 for users with a body size of 120 – 160 cm	
with power supported height and angle adjustment (wired)	with power supported height and angle adjustment and contourable upper body support (divided into shoulder, trunk and pelvis areas – areas can be individually fixed in the desired position with locking levers)
Price / piece	Price / piece
42 03 000	5.009 €
	42 03 100
	5.991 €

Basic model	
till. size 4 for users with a body size of 145 – 185 cm	
with power supported height and angle adjustment (wired)	with power supported height and angle adjustment and contourable upper body support (divided into shoulder, trunk and pelvis areas – areas can be individually fixed in the desired position with locking levers)
Price / piece	Price / piece
42 04 000	5.257 €
	42 04 100
	6.356 €



Contourable upper body support




Akku	
	Rechargeable battery (Li-Ion 2.25 Ah) – for basic models with power supported height adjustment (Note: can only be retrofitted by trained personnel)
	Price / piece 42 01 026 514 €

Covers for size 1 + 2	
Black / Deep sea blue (No. 105 + 111)	Black / Pink 2 (No. 105 + 115)
	
Piece	Piece
for size 1+2 (9100239)	for size 1+2 (9100234)
incl.	incl.

Covers for size 3 + 4	
Black / Petrol (No. 105 + 112)	Black / Dark grey (No. 105 + 92)
	
Piece	Piece
for size 3+4 (9100237)	for size 3+4 (9100236)
incl.	incl.


Covers	
Cover according to customer request	
	Upcharge / piece
42 00 999	278 €

Headrests		
Headrest (Standard) – padded, incl. cover, height, depth and angle adjustment	Headrest (Shell-shaped) – padded, incl. cover, height, depth and angle adjustment	Headrest (Lateral guide) – padded, incl. cover, height, depth and angle adjustment
		
Price / piece	Price / piece	Price / piece
for size 1 (42 01 001)	for size 1 (42 01 003)	for size 1 (42 01 002)
for size 2 (42 02 001)	for size 2 (42 02 003)	for size 2 (42 02 002)
for size 3+4 (42 03 001)	for size 3+4 (42 03 003)	for size 3+4 (42 03 002)
287 €	392 €	392 €
287 €	392 €	392 €
287 €	392 €	392 €

Headrests						
Headrest (Lateral guide, angle adjustable) – padded, incl. cover, height, depth and angle adjustable, side guides are each angle adjustable		Adapter for lateral adjustment of headrests		Universal adapter for tailored headrests – 15x15 mm, height and depth adjustable		
						
		A	B	C		
Size 1	13,5cm	10,5cm	10cm			
Size 2	15cm	10,5cm	13,5cm			
Size 3	17cm	12,5cm	16,5cm			
		Price / piece		Price / piece		
Size 1 (42 01 033)	651 €	42 02 004		154 €	42 01 004	221 €
Size 2 (42 02 033)	651 €					
Size 3 (42 03 033)	651 €					


Support pads

Lateral pads for lying surface (standard) – incl. upper arm support, flip away, stepless adjustment for height and width, adjustable for rotation



	Height (A)	Depth (C)	suitable for	Chest width (B)
Size 1	10,5cm	9 cm	till. Size 1	14 – 27 cm
			till. Size 2	17 – 31 cm
Size 2	11,5cm	11 cm	till. Size 2	17 – 31 cm
			till. Size 3	20 – 37 cm
Size 3	12,5cm	16 cm	till. Size 3	20 – 37 cm
			till. Size 4	29 – 46 cm
Size 4	13,5cm	18 cm	till. Size 3	20 – 37 cm
			till. Size 4	27 – 46 cm

Pelvic pads for lying surface (standard) – incl. forearm rest, flip away, stepless adjustment for height and width, adjustable for rotation



	Height (A)	Depth (C)	suitable for	Chest width (B)
Size 1	12 cm	9 cm	till. Size 1	14 – 27 cm
			till. Size 2	17 – 31 cm
Size 2	14 cm	10,5 cm	till. Size 2	17 – 31 cm
			till. Size 3	20 – 37 cm
Size 3	17 cm	15,5 cm	till. Size 3	20 – 37 cm
			till. Size 4	29 – 46 cm
Size 4	21 cm	17,5 cm	till. Size 3	20 – 37 cm
			till. Size 4	27 – 46 cm



The **standard version** of the lateral and pelvis pads is only suitable for the **till. without** contourable upper body support!

The **depth adjustable version** of the lateral and pelvic pads is only suitable for the **till. (size 3 & 4) with** contourable upper body support!

Piece / pair	Piece / pair	Piece / pair	Piece / pair
Size 1 (42 01 005) 325 €		Size 1 (42 01 007) 325 €	
Size 2 (42 02 005) 325 €	Size 2, depth adjustable (42 02 006) 325 €	Size 2 (42 02 007) 325 €	Size 2, depth adjustable (42 02 008) 325 €
Size 3 (42 03 005) 325 €	Size 3, depth adjustable (42 03 006) 325 €	Size 3 (42 03 007) 325 €	Size 3, depth adjustable (42 03 008) 325 €
Size 4 (42 04 005) 325 €	Size 4, depth adjustable (42 04 006) 325 €	Size 4 (42 04 007) 325 €	Size 4, depth adjustable (42 04 008) 325 €


Pads

Knee support pads with patella recess

individual three-dimensional adjustment, with positioning aid




	Inside width (A)
Size 1	8 cm
Size 2	10 cm
Size 3	12 cm



Flexible knee positioning (FKP) – Individually three- dimensionally adjustable, with length-ad-justable contact straps that adapt to the shape of the knee



	Inside width (A)
Size 1	8 cm
Size 2	10 cm
Size 3	12 cm



Padding – for additional support of the knee, to prevent hyperextension (knee overextended)



	Length	Ø
Size 1	8 cm	4 cm
Size 2	8 cm	6 cm
Size 3	10 cm	8 cm
Size 4	12 cm	10 cm

Piece / pair	Piece / pair	Piece / pair
Size 1 (42 01 035) 576 €	Size 1 (42 01 036) 649 €	Size 1 (42 01 011) 179 €
Size 2 (42 02 035) 576 €	Size 2 (42 02 036) 649 €	Size 2 (42 02 011) 179 €
Size 3 (42 03 035) 576 €	Size 3 (42 03 036) 649 €	Size 3 (42 03 011) 179 €
		Size 4 (42 04 011) 179 €

Leg configuration

Adapter for leg rests (standard)
(standard configuration, if nothing else stated)



Piece / pair

for size 1+2 (**42 01 014**) incl.
for size 3+4 (**42 03 014**) incl.

Adapter for leg rests (abduction)
enables individual and continuous abduction of each leg from 0°–30° (Note: not retrofittable!)



Upcharge / pair

for size 1+2 (**42 01 015**) 881 €
for size 3+4 (**42 03 015**) 881 €

Leg rests (standard) – length adjustable, with 5 cm pelvic rotation adjustment (standard configuration if nothing else stated)



Piece / pair

for size 1 (**42 01 016**) incl.
for size 2 (**42 02 016**) incl.
for size 3 (**42 03 016**) incl.
for size 4 (**42 04 016**) incl.

Leg rests (flex) – adjustable in length and angle, with 5 cm pelvic rotation adjustment (Note: not retrofittable!)



Upcharge / pair

for size 1 (**42 01 017**) 533 €
for size 2 (**42 02 017**) 533 €
for size 3 (**42 03 017**) 533 €
for size 4 (**42 04 017**) 533 €

Foot configuration

Split footplates (standard)
standard configuration if nothing else stated (Note: not retrofittable!)



Piece / pair

for size 1 (**42 01 018**) incl.
for size 2 (**42 02 018**) incl.
for size 3 (**42 03 018**) incl.
for size 4 (**42 04 018**) incl.

Split footplates (flex) – individually adjustable, dorsiflex (20°), plantarflex (15°) and inversion/eversion (+/-15°) (Note: not retrofittable!)



Upcharge / pair

for size 1 (**42 01 019**) 269 €
for size 2 (**42 02 019**) 269 €
for size 3 (**42 03 019**) 269 €
for size 4 (**42 04 019**) 269 €

Accessories

Heel edges – incl. belt + magnetic lock (Note: cannot be combined with simple foot strap!)



Piece / pair

for size 1 (**42 01 021**) 195 €
for size 2 (**42 02 021**) 195 €
for size 3 (**42 03 021**) 195 €
for size 4 (**42 04 021**) 195 €

Simple foot strap – with magnetic lock (Note: cannot be combined with heel edges!)




Piece / pair

for size 1 (**42 01 022**) 119 €
for size 2 (**42 02 022**) 119 €
for size 3 (**42 03 022**) 119 €
for size 4 (**42 04 022**) 119 €

Therapietisch

Wood therapy tray (standard) – with rail, removable, depth, height and angle adjustable



Price / piece

for size 1 (**42 01 024**) 832 €
for size 2 (**42 02 024**) 832 €
for size 3 (**42 03 024**) 832 €
for size 4 (**42 04 024**) 832 €

Wood therapy tray (with quick adjustment) – with rail, removable, depth, height and angle adjustable via table linkage and additionally via quick adjustment



Price / piece

for size 1 (**42 01 034**) 1.003 €
for size 2 (**42 02 034**) 1.003 €
for size 3 (**42 03 034**) 1.003 €
for size 4 (**42 04 034**) 1.003 €

Padded tray cover – covered with black wipe clean material, height: 1 cm



Price / piece

for size 1 (**42 01 025**) 174 €
for size 2 (**42 02 025**) 174 €
for size 3 (**42 03 025**) 174 €
for size 4 (**42 04 025**) 174 €

Accessories

for size 1+2



for size 3+4



Push bar

Price / piece

for size 1+2 (**42 01 012**) 342 €
for size 3+4 (**42 03 012**) 342 €



Safety bolt – for release lever, prevents unintentional release of height and angle adjustments (Note: only for size 1+2)

Price / piece

42 00 013 84 €

Harness

Chest strap – with pad and adjustable clip on both sides, colour: black (standard configuration if nothing else stated: for Basic model size 1 = belt size 1, for Basic model size 2 = belt size 2, and for Basic model sizes 3 & 4 = belt size 3.)



	Size 1	Size 2	Size 3
Pad width	4 cm	6 cm	6 cm
Pad length per side	18 cm	25 cm	30 cm
Strap length per side	25 cm	45 cm	65 cm
Strap width	2,5 cm	2,5 cm	2,5 cm
Price / piece			
Size 1 (42 01 029)			incl.
Size 2 (42 02 029)			incl.
Size 3 (42 03 029)			incl.

Wide Chest strap – padded, colour: black




	Size 1	Size 2	Size 3
Pad width	8 cm	10 cm	12 cm
Pad length per side	10 cm	25 cm	35 cm
Strap length per side	25 cm	45 cm	65 cm
Strap width	2,5 cm	2,5 cm	2,5 cm
Upcharge / piece			
Size 1 (42 01 030)			63 €
Size 2 (42 02 030)			63 €
Size 3 (42 03 030)			63 €

Harness

Positioning jacket – made of elastic, breathable and moisture repellent material, shoulder straps and pelvis straps (25 mm) with clips, colour: black

	A	B	C	D
Size 1	14 cm	26,5 cm	14 cm	33 cm
Size 2	15 cm	29,5 cm	16 cm	36 cm
Size 3	17 cm	35,5 cm	19 cm	41 cm
Size 4	18 cm	38,5 cm	21 cm	43,5 cm



Extension for positioning jacket – for the belt attachment of the positioning jacket

	recommended from a neck width of
for size 1	13 cm
for size 2	17 cm
for size 3	21 cm
for size 4	30 cm



	Upcharge / piece	Piece / set
Size 1 (42 01 027)	251 €	
Size 2 (42 02 027)	251 €	
Size 3 (42 03 027)	251 €	
Size 4 (42 04 027)	251 €	
for size 1+2 (42 01 028)		61 €
for size 3+4 (42 03 028)		61 €

Harness

Pelvic harness –with pad and adjustable clip on both sides, colour: black (standard configuration if nothing else stated: for Basic model size 1 = belt size 1, for Basic model size 2 = belt size 2, and for Basic model sizes 3 & 4 = belt size 3.)



Pelvic harness – by „Pelvi.Loc“, colour: black



	Size 1	Size 2	Size 3		Size 1 (XS)	Size 2 (S)	Size 3 (M)	Size 4 (L)
for pelvic width	16 – 25 cm	25 – 40 cm	38 – 46 cm	for pelvic width	16 – 25 cm	25 – 34 cm	32 – 40 cm	38 – 46 cm
Pad width	4 cm	6 cm	6 cm	Pad width	6,5 cm	7 cm	7,5 cm	8,5 cm
Pad length per side	18 cm	25 cm	30 cm	Pad length per side	16 cm	19 cm	22 cm	24 cm
Strap length per side	25 cm	45 cm	65 cm	Strap length per side	25 cm	25 cm	45 cm	45 cm
Strap width	2,5 cm	2,5 cm	2,5 cm	Strap width	2,5 cm	2,5 cm	2,5 cm	2,5 cm
Price / piece				Upcharge / piece				
Size 1 (42 01 031)			incl.	Size 1 (XS) (42 01 032)				325 €
Size 2 (42 02 031)			incl.	Size 2 (S) (42 02 032)				325 €
Size 3 (42 03 031)			incl.	Size 3 (M) (42 03 032)				325 €
				Size 4 (L) (42 04 032)				325 €

Accessories



75 mm Castors –5 (for sizes 1 – 3) or 6 (for size 4) bra-kable castors with total locking, one of which only has directional locking. Size 1+2 = double castors; Size 3+4 = single castors (standard configuration if nothing else stated)

Piece / set
for size 1–4 (42 00 038)



100 mm castors – only possible with size 3 + 4, 5 (for size 3) or 6 (for size 4) braked castors with total locking, one with directional locking

Upcharge / Set
for size 3+4 (42 01 038)