



Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The standing frame – from prone to supine.


- Basic frame with height and depth adjustable footplate
- Angular adjustment of the central column via gas pressure spring from (-)15° – 30°
- Incl. chest support pad with lateral guide and buttocks support pad
- Height and depth adjustable pelvis frame, pelvis position adjustable via a crank
- Incl. three-dimensionally adjustable knee support pads with patella cutout
- Incl. therapy tray
- 4 castors with locking devices (75 mm) and total lock

HMV-No.: 28.29.01.1059



Please send us the following measurements for optimum presetting of the standing frame

Settings



The following measurements should be pre-set on the standing frame:

- | | |
|--------------------------------|-------------------------------|
| I Body size _____ | V Chest width _____ |
| II Sternum height _____ | VI Pelvic width _____ |
| III Elbow height _____ | VII Knee width _____ |
| IV Crotch height _____ | VIII Knee height _____ |


48 00 020

Piece

incl.

Dimensions

		Size 1	Size 2	Size 3
	Body size	75 – 110 cm	100 – 140 cm	120 – 160 cm
A	Height from footplate to tray	48,5 – 80,5 cm	71 – 95 cm	80 – 112 cm
B	Pelvis depth	12 – 22 cm	16 – 28 cm	20 – 35 cm
	Depth adjustment of pelvis frame (via crank)	6 cm		
C	Total length	84 cm	90 cm	100 cm
D	Total width	57 cm	65 cm	71 cm
E	Abduction divided footplate	0° – 27°		
F	Pelvis height	45 – 60 cm	57 – 80 cm	65 – 90 cm
G	Sternum height	55 – 85 cm	75 – 110 cm	95 – 125 cm
	Max. load	40 kg	60 kg	80 kg
	Weight	27 kg	29 kg	32 kg



Basic model

		
tim. size 1 for users with a body size 75 – 110 cm Melon yellow (92 00 006)	tim. size 2 for users with a body size 100 – 140 cm Dormant Red (92 00 003)	tim. size 3 for users with a body size 120 – 160 cm Lasurit blue metallic (92 00 010)
Price / piece	Price / piece	Price / piece
48 01 100 3.690 €	48 02 100 3.690 €	48 03 100 3.994 €

Frame colours

Coated with an individual colour (on request) Please provide the requested colour here:

RAL _____

Upcharge

48 00 090 278 €

Coated covers

Flexi Black



Piece

91 00 043 incl.

Flexi Light grey



Piece

91 00 050 incl.

Lateral pads

Chest support pad with lateral guide
can be adjusted in height, width
and depth



Chest width

Size 1	16 – 24 cm
Size 2	18 – 26 cm
Size 3	24 – 34 cm
Size 4	30 – 42 cm

Chest support pad with lateral guide
can be adjusted in height, width
and depth, additionally angle
adjustable



Chest width

Size 1	16 – 24 cm
Size 2	18 – 26 cm
Size 3	24 – 34 cm
Size 4	30 – 42 cm

Price / piece

Upcharge / piece

Size 1 (48 01 101)	incl.	Size 1 (48 01 102)	122 €
Size 2 (48 02 201)	incl.	Size 2 (48 02 202)	122 €
Size 3 (48 02 101)	incl.	Size 3 (48 02 102)	122 €
Size 4 (48 03 101)	incl.	Size 4 (48 03 102)	122 €

Accessories

Chin pad – height adjustable



Inside width

Size 1	6 cm
Size 2	8 cm
Size 3	10 cm

Price / piece

Size 1 (48 01 012)	271 €
Size 2 (48 02 012)	271 €
Size 3 (48 03 012)	271 €

Back strap



Pad width

Pad height

Size 1	17 cm	8,5 cm
Size 2	28 cm	10 cm
Size 3	39 cm	11,5 cm


Price / piece

Size 1 (48 01 013)	162 €
Size 2 (48 02 013)	162 €
Size 3 (48 03 013)	162 €

Combinations - pelvic frame

Option 1:

Pelvis support pad with lateral guide (mounted at the front of the pelvis frame) + buttocks support pad (mounted at the back of the pelvis frame)



Pelvis support pad with lateral guide width adjustable



	Pelvic width
Size 1	16 – 24 cm
Size 2	18 – 26 cm
Size 3	24 – 34 cm
Size 4	30 – 42 cm
Price / piece	
Size 1 (48 01 005)	426 €
Size 2 (48 02 105)	426 €
Size 3 (48 02 005)	426 €
Size 4 (48 03 005)	426 €


Buttocks support pad




	Width x Height
Size 1	16 x 10 cm
Size 2	20 x 10 cm
Size 3	25 x 14 cm
Price / piece	
Size 1 (48 01 003)	incl.
Size 2 (48 02 003)	incl.
Size 3 (48 03 003)	incl.

Option 2:

Pelvic positioner (mounted at the front of the pelvis frame) + pelvis support pad with lateral guide (mounted at the back of the pelvis frame)



Pelvic positioner



	Width x Height
Size 1	16 x 10 cm
Size 2	20 x 10 cm
Size 3	25 x 14 cm
Price / piece	
Size 1 (48 01 004)	162 €
Size 2 (48 02 004)	162 €
Size 3 (48 03 004)	162 €


Pelvis support pad with lateral guide width adjustable




	Pelvic width
Size 1	16 – 24 cm
Size 2	18 – 26 cm
Size 3	24 – 34 cm
Size 4	30 – 42 cm
Upcharge / piece	
Size 1 (48 01 006)	263 €
Size 2 (48 02 106)	263 €
Size 3 (48 02 006)	263 €
Size 4 (48 03 006)	263 €

Option 3:

Pelvic positioner and combined spine and pelvic support pads (mounted at the front of the pelvis frame) + buttocks support pad (mounted at the back of the pelvis frame)




Pelvic positioner



	Width x Height
Size 1	16 x 10 cm
Size 2	20 x 10 cm
Size 3	25 x 14 cm
Price / piece	
Size 1 (48 01 004)	162 €
Size 2 (48 02 004)	162 €
Size 3 (48 03 004)	162 €

Combined spine and pelvic support pads – width, depth and height adjustable



	Pelvic width
for size 1	14 – 24 cm
for size 2	22 – 32 cm
for size 3	32 – 42 cm
Piece / pair	
for size 1 (48 01 007)	426 €
for size 2 (48 02 007)	426 €
for size 3 (48 03 007)	426 €

Buttocks support pad



	Width x Height
Size 1	16 x 10 cm
Size 2	20 x 10 cm
Size 3	25 x 14 cm
Price / piece	
Size 1 (48 01 003)	incl.
Size 2 (48 02 003)	incl.
Size 3 (48 03 003)	incl.

Accessories

Back support pad – mounted at the back of the pelvis frame, height and angle adjustable

	Width	Height
for size 1	20 cm	23 cm
for size 2	20 cm	23 cm
for size 3	25 cm	28 cm
Price / piece		
for size 1 (48 01 010)		452 €
for size 2 (48 02 010)		452 €
for size 3 (48 03 010)		452 €



Headrest – height, depth and angle adjustable, only for use together with the back support pad

	Price / piece
for size 1 (48 01 111)	348 €
for size 2 (48 02 111)	348 €
for size 3 (48 03 111)	348 €




Configuration options - knee pads

Knee support pads with patella recess individual three-dimensional adjustment



	Inside width		Inside width
Size 1	8 cm	Size 1	8 cm
Size 2	10 cm	Size 2	10 cm
Size 3	12 cm	Size 3	12 cm

Flexible knee positioning (FKP) – Individually three-dimensionally adjustable, with length-adjustable contact straps that adapt to the shape of the knee



Knee pad support – 2 different options are available to accommodate the adjacent knee pads

Rigid
angle, depth and height adjustable

Ball head
angle-, depth-, height- and additionally tilt angle adjustable

(Rigid = Standard configuration unless otherwise specified)

Please select one of the adjacent knee pads.

Piece / pair		Upcharge / pair		Upcharge / pair	
Size 1 (48 01 027)	incl.	Size 1 (48 01 028)	122 €	Rigid (48 00 026)	incl.
Size 2 (48 02 027)	incl.	Size 2 (48 02 028)	122 €	Ball head (48 01 026)	160 €
Size 3 (48 03 027)	incl.	Size 3 (48 03 028)	122 €		

Accessories

Therapy tray



Therapy tray – with adjustable armrests



Therapy tray – with adjustable armrests and bowl (incl. cover)



Price / piece		Upcharge / piece		Upcharge / piece	
for size 1 (48 01 114)	incl.	for size 1 (48 01 116)	232 €	for size 1 (48 01 117)	502 €
for size 2 (48 02 114)	incl.	for size 2 (48 02 116)	232 €	for size 2 (48 02 117)	502 €
for size 3 (48 03 114)	incl.	for size 3 (48 03 116)	232 €	for size 3 (48 03 117)	502 €

Accessories

Pad for therapy tray without armrests
covered with black Flexi cover, height: 1 cm



Pad for therapy tray with armrests
covered with black Flexi cover, height: 1 cm



Angular adjustment of the center column by gas spring from (-) 15° - 30°
(standard configuration = with foot pedal, if nothing else stated)



Price / piece		Price / piece		Price / piece	
for size 1 (4801018)	181 €	for size 1 (4801019)	232 €	by foot pedal on the base frame (4800025)	incl.
for size 2 (4802018)	181 €	for size 2 (4802019)	232 €	by manual release under the tray (4801025)	incl.
for size 3 (4803018)	181 €	for size 3 (4803019)	232 €		

Accessories

Footplate – height and depth adjustable



Individual footplates – height, width and angle adjustable (abduction of legs)



Heel edges with ankle straps – angle, width and depth adjustable, for mounting on the footplate(s) with magnetic locks

