



Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The sporty and agile standing frame.


- **Basic frame (standard)** with height and depth adjustable footplate and four height adjustable support castors (75 mm), rear with springs and lock
- **Basic frame (with abduction)** with individual adjustable split footplate, two height adjustable support castors (75 mm) in the front and one support castor (50 mm) in the back
- Incl. chest support pad, rear support pad and three-dimensionally adjustable knee support pads with patella recess
- Mechanical angle adjustment of the central column from 0 - 15°
- Incl. drum brakes, easy operation via locking lever
- Incl. tray support and storage net (for size 3+4)
- Adjustable camber

HMV-No.: 28.29.01.3011



Please send us the following measurements for optimum presetting of the standing frame

Settings



The following measurements should be pre-set on the standing frame


- | | |
|--------------------------------|-------------------------------|
| I Body size _____ | V Knee height _____ |
| II Sternum height _____ | VI Chest width _____ |
| III Elbow height _____ | VII Pelvic width _____ |
| IV Crotch height _____ | VIII Knee width _____ |

49 00 020

Piece

incl.

Dimensions



		Size 1		Size 2			Size 3		Size 4			
	Body size	75 – 110 cm		100 – 130 cm			120 – 155 cm		150 – 175 cm			
	Wheel size	24"	28"	30"	32"	36"	32"	36"	36"	40"	42"	
A	Grip height: footplate – hand rim at top	40 – 48 cm	50 – 58 cm	57 – 65 cm	64 – 69 cm	72 – 80 cm	60 – 68 cm	71 – 79 cm	82 – 90 cm	87 – 95 cm		
	Recommended from elbow height	50 cm	60 cm	70 cm	80 cm	92 cm	82 cm	92 cm	92 cm	97 cm	99 cm	
B	Width of hand rim at top	0°	–	–	–	–	65 cm		67 cm			
		3°	–	–	54 cm	53 cm	52 cm	61 cm	60 cm	63 cm	62 cm	61 cm
		6°	–	–	51 cm	49 cm	47 cm	57 cm	55 cm	59 cm	57 cm	56 cm
		9°	45 cm	44 cm	48 cm	45 cm	42 cm	54 cm	51 cm	55 cm	53 cm	51 cm
		12°	42 cm	41 cm	45 cm	41 cm	38 cm	–	–	–	–	–
		15°	40 cm	–	–	–	–	–	–	–	–	–
C	Width of hand rim at the bottom (min. – max.)	68 – 74 cm	70 – 73 cm	71 – 83 cm		71 – 84 cm	73 – 83 cm		76 – 86 cm		76 – 90 cm	
D	Height from footplate to centre of knee support pads	19 – 30 cm		28 – 36 cm			34 – 44 cm		42 – 51 cm			
E	Height from footplate to support of buttocks support pad	26 – 42 cm		43 – 58 cm			53 – 69 cm		65 – 80 cm			
	Recommended for inside leg length	30 – 47 cm		48 – 63 cm			58 – 74 cm		70 – 84 cm			
F	Abduction split footplate	0° – 6°					0° – 30°					
F	Abduction split footplate (with abducted basic model)	0° – 40°					–		–			
G	Pelvis dept*	11 – 24 cm					15 – 28,5 cm					
	Pelvis dept** (adjustment range of pelvis position via crank max. 6 cm)	11 – 20 cm (+6 cm)					15 – 27 cm (+6 cm)					
H	Sternum height (height from footplate to centre of chest support pad)	54 – 80 cm		68 – 101 cm			92 – 123 cm		111 – 140 cm			
I	Tray height from footplate	43 – 66 cm		61 – 79 cm			72 – 96 cm		92 – 110 cm			
	Max. load	30 kg		45 kg			65 kg		80 kg			
	Weight	28 kg		31 kg			38 kg		42 kg			
* with bracket for buttocks support pad (standard) ** with bracket for buttocks support pad (flip away)												

* with bracket for buttocks support pad (standard) ** with bracket for buttocks support pad (flip away)

Basic model

todd. size 1 for users with a body size 75 – 110 cm		todd. size 2 for users with a body size 100 – 130 cm		todd. size 3 for users with a body size 120 – 155 cm		todd. size 4 for users with a body size 150 – 175 cm	
Price / piece		Price / piece		Price / piece		Price / piece	
Standard	(49 01 000) 3.126 €	Standard	(49 02 000) 3.275 €	Standard	(49 03 000) 3.935 €	Standard	(49 04 000) 4.447 €
Abductable	(49 01 100) 3.504 €	Abductable	(49 02 100) 3.658 €				

Basic model

Melon yellow	Anthracite grey	Pure white	Coating with a colour of your choice – (by arrangement) Please specify the desired colour here: RAL _____
			
Piece	Piece	Piece	Upcharge
92 00 006 incl.	92 00 017 incl.	92 00 100 incl.	49 00 090 278 €


Wipe-clean covers

Flexi Black	Flexi Light grey
	
Piece	Piece
49 00 043 incl.	49 00 050 incl.

Lateral pads
Accessories
Chest system – adjustable in height

Chest support pad with lateral guide
can be adjusted in height, width and depth

Chest support pad with lateral guide
can be adjusted in height, width and depth, additionally angle adjustable

Back strap – with magnetic lock for chest support pad with lateral guide


			Chest width		Chest width		Pad width		Pad height
	Inside width	Height	Size 1		Size 1		Size 1		
Size 1	27 cm	13 cm	Size 2	16 – 24 cm	Size 2	16 – 24 cm	Size 2	17 cm	8,5 cm
Size 2	36 cm	16 cm	Size 3	18 – 26 cm	Size 3	18 – 26 cm	Size 3	28 cm	10 cm
Size 3	43 cm	19 cm	Size 4	24 – 34 cm	Size 4	24 – 34 cm	Size 4	39 cm	11,5 cm
				30 – 42 cm		30 – 42 cm			
Price / piece			Upcharge / piece		Upcharge / piece		Price / piece		
Size 1 (49 01 001)		incl.	Size 1 (49 01 002)	382 €	Size 1 (49 01 003)	382 €	Size 1 (49 01 024)		162 €
Size 2 (49 02 001)		incl.	Size 2 (49 02 002)	382 €	Size 2 (49 02 003)	382 €	Size 2 (49 02 024)		162 €
Size 3 (49 03 001)		incl.	Size 3 (49 03 002)	382 €	Size 3 (49 03 003)	382 €	Size 3 (49 03 024)		162 €
			Size 4 (49 04 002)	382 €	Size 4 (49 04 003)	382 €			

Accessories
Back support pad

can be adjusted in height and angle

	Width	Height
for size 1+2	20 cm	23 cm
for size 3+4	25 cm	28 cm
Price / piece		
for size 1 (49 01 025)		452 €
for size 2 (49 02 025)		452 €
for size 3+4 (49 03 025)		452 €


Headrest – can be adjusted in height, depth and angle, only for use together with the back support pad

Price / piece	
for size 1 (49 01 026)	348 €
for size 2 (49 02 026)	348 €
for size 3 (49 03 026)	348 €
for size 4 (49 04 026)	348 €


Configuration options - Knee pads
Knee support pads with patella recess
individual three-dimensional adjustment

Flexible knee positioning (FKP)

Individually three-dimensionally adjustable, with length-adjustable contact straps that adapt to the shape of the knee


Knee pad support – 2 different options are available to accommodate the adjacent knee pads

Rigid
angle, depth and height adjustable

Ball head
angle-, depth-, height- and additionally tilt angle adjustable

(Rigid = standard configuration unless otherwise specified)

	Inside width		Inside width
Size 1	8 cm	Size 1	8 cm
Size 2	10 cm	Size 2	10 cm
Size 3	12 cm	Size 3	12 cm


Please select one of the adjacent knee pads.

Piece / pair		Upcharge / pair		Upcharge / pair	
Size 1 (49 01 031)	incl.	Size 1 (49 01 032)	122 €	Rigid (49 00 030)	incl.
Size 2 (49 02 031)	incl.	Size 2 (49 02 032)	122 €	Ball head (49 01 030)	160 €
Size 3 (49 03 031)	incl.	Size 3 (49 03 032)	122 €		

Configuration possibilities - Bracket for buttock and pelvic positioning

Configuration 1:



Bracket for buttocks support pad (standard)
depth adjustable via quick adjustment

	Depth adjustment Pelvic pad up to buttock pad	Standard configuration
Size 1	11 – 24 cm	for size 1+2
Size 2	15 – 28,5 cm	for size 3+4
Price / piece		
Size 1 (49 01 005)		incl.
Size 2 (49 03 005)		incl.

Configuration 2:



Bracket for buttocks support pad (flip away)
for step-by-step adjustment in the depth via hand crank and quick adjustment, fold-down mechanism for easier transfer (approx. 20 ° (10 cm) downwards)

	Depth adjustment Pelvic pad up to buttock pad	Depth adjustment via Crank	Standard configuration
Size 1	11 – 20 cm	6 cm	for size 1+2
Size 2	15 – 27 cm	6 cm	for size 3+4
Upcharge / Set			
Size 1 (49 01 004)			414 €
Size 2 (49 03 004)			414 €

Configuration 3:



Pelvic frame (with swivel arm) – variable depth adjustment, additional depth adjustment by hand crank for gradual and graded hip extension of the user, incl. holder for rear support pad and back support pad

 Please consider the combination options below and select a pad configuration on page 112.

	Width outside	Width inside	Depth adjustment Pelvic pad up to buttock pad	Depth adjustment via Crank
Size 1	37,5 cm	31 cm	11 – 22 cm	6 cm
Size 2	47,5 cm	41 cm	14 – 28 cm	6 cm
Size 3	57,5 cm	51 cm	18 – 35 cm	6 cm
Upcharge / Set				
Size 1 (49 01 029)				548 €
Size 2 (49 02 029)				548 €
Size 3 (49 03 029)				548 €

Combination options

The pelvic frame (with swivel arm) is available in three different sizes. Please take into account that not every camber setting is possible with every size of pelvic frame, as a collision with the drive wheels may occur above a certain camber setting. The following overview shows which size of pelvic frame can be combined with which basic frame.

	Camber	todd. (size 1)		todd. (size 2)			todd. (size 3)		todd. (size 4)		
Wheel size		24"	28"	30"	32"	36"	32"	36"	36"	40"	42"
	0°						● ■ ▲	● ■ ▲	● ■ ▲	● ■ ▲	● ■ ▲
	3°			● ■	● ■	● ■	● ■	● ■	● ■	● ■	● ■
	6°			●	●	●	● ■	● ■	● ■	● ■	● ■
	9°	●	●	●	●	●	● ■	●	● ■	● ■	●
	12°			●							
	15°										
● = Pelvic frame (size 1) ■ = Pelvic frame (size 2), ▲ = Pelvic frame (size 3))											

Configuration possibilities - Pelvic positioning

Option 1:

Pelvis support pad with lateral guide (mounted at the front of the pelvis frame) + buttocks support pad (mounted at the back of the pelvis frame)



i This combination is possible with pelvic frame (with swivel arm) as follows:

Pelvis support pad with lateral guide	Pelvic frame (with swivel arm)		
	Size 1	Size 2	Size 3
Size 1	16 - 24 cm	-	-
Size 2	-	18 - 25 cm	-
Size 3	-	24 - 32 cm	24 - 34 cm
Size 4	-	-	29 - 43 cm

Pelvis support pad with lateral guide
width and depth adjustable



	Pelvic width	Depth adjustment
Size 1	16 - 24 cm	5 cm
Size 2	18 - 26 cm	5 cm
Size 3	24 - 34 cm	5 cm
Size 4	30 - 42 cm	5 cm
		Price / piece
Size 1	(49 01 006)	426 €
Size 2	(49 02 106)	426 €
Size 3	(49 02 006)	426 €
Size 4	(49 03 006)	426 €


Buttocks support pad



	Width x Height
Size 1	16x10cm
Size 2	20x10cm
Size 3	25x14cm
Price / piece	
Size 1 (49 01 007)	incl.
Size 2 (49 02 007)	incl.
Size 3 (49 03 007)	incl.

Option 2:

Pelvic positioner (mounted at the front of the pelvis frame) + pelvis support pad with lateral guide (mounted at the back of the pelvis frame)



i This combination is possible with pelvic frame (with swivel arm) as follows:

Pelvis support pad with lateral guide	Pelvic frame (with swivel arm)		
	Size 1	Size 2	Size 3
Size 1	16 - 24 cm	-	-
Size 2	-	18 - 25 cm	-
Size 3	-	24 - 32 cm	24 - 34 cm
Size 4	-	-	29 - 43 cm

Pelvic positioner – depth adjustable



	Width x Height	Depth adjustment
Size 1	16 x 10 cm	5 cm
Size 2	20 x 10 cm	5 cm
Size 3	25 x 14 cm	5 cm
		Price / piece
Size 1	(49 01 008)	162 €
Size 2	(49 02 008)	162 €
Size 3	(49 03 008)	162 €


Pelvis support pad with lateral guide
width and depth adjustable



	Pelvic width
Size 1	16 – 24 cm
Size 2	18 – 26 cm
Size 3	24 – 34 cm
Size 4	29 – 43 cm
Upcharge / piece	
Size 1 (49 01 010)	263 €
Size 2 (49 02 110)	263 €
Size 3 (49 02 010)	263 €
Size 4 (49 03 010)	263 €

Option 3:

Pelvic positioner (mounted at the front of the pelvis frame) + pelvis support pad with lateral guide (mounted at the back of the pelvis frame)



Pelvic positioner depth adjustable




	Width x Height	Depth adjustment
Size 1	16 x 10 cm	5 cm
Size 2	20 x 10 cm	5 cm
Size 3	25 x 14 cm	5 cm
		Price / piece
Size 1	(49 01 008)	162 €
Size 2	(49 02 008)	162 €
Size 3	(49 03 008)	162 €

Combined spine and pelvic support pads – width, depth and height adjustable



	Pelvic width	Depth adjustment
Size 1	14 - 24 cm	5 cm
Size 2	22 - 32 cm	5 cm
Size 3	32 - 42 cm	5 cm
		Piece / pair
Size 1*	(49 01 009)	426 €
Size 2	(49 02 009)	426 €
Size 3	(49 03 009)	426 €

Buttocks support pad



	Width x Height
Size 1	16 x 10 cm
Size 2	20 x 10 cm
Size 3	25 x 14 cm
Price / piece	
Size 1 (49 01 007)	incl.
Size 2 (49 02 007)	incl.
Size 3 (49 03 007)	incl.

* For the combined spine and pelvic support pads, size 1, no pelvic positioner can be selected when used in combination.

Accessories



Therapy tray – made of wood, removable, with rail, incl. height adjustable attachment for chest support pad (The minimum setting „Table height from base plate“ (I) increases by 5 cm)

	Width	Depth
for size 3+4	70 cm	58 cm
		Price / piece
for size 3+4	(49 03 028)	414 €

for size 1+2



for size 3+4



Therapy tray with rail – made of Perspex, removable, not usable in combination with basket

	Width	Depth
for size 1+2	52 cm	40 cm
for size 3+4	62 cm	50 cm
		Price / piece
for size 1+2	(49 01 013)	414 €
for size 3+4	(49 03 027)	414 €

Accessories



Basket – removable, not usable in combination with therapy tray

	Price / piece
for size 1+2 (49 00 014)	181 €



Mounting position of the brake – the brake lever can be mounted under the tray on the left or right, as desired (standard configuration = under the tray on the right)

	Piece
under tray on the left (49 01 015)	incl.
under tray on the right (49 02 015)	incl.

Rear wheels

Drive wheels – with drum brake, push rims and black PU tyres



Drive wheels with closed hand rim

with drum brake and black PU tyres



		Piece / pair
24" for size 1	49 00 033	700 €
28" for size 1	49 01 033	852 €
30" for size 2	49 02 033	858 €
32" for size 2+3	49 03 033	884 €
36" for size 2-4	49 04 033	944 €
40" for size 4	49 05 033	1.097 €
42" for size 4	49 06 016	1.856 €

		Piece / pair
24" for size 1	49 00 034	918 €
28" for size 1	49 01 034	1.163 €
30" for size 2	49 02 034	1.170 €
32" for size 2+3	49 03 034	1.195 €
36" for size 2-4	49 04 034	1.264 €
40" for size 4	49 05 034	1.415 €

Rear wheels



Spoke guard – transparent

		Piece / pair
24" for size 1	(49 00 018)	239 €
28" for size 1	(49 01 018)	239 €
30" for size 2	(49 02 018)	239 €
32" for size 2+3	(49 03 018)	348 €
36" for size 3	(49 04 018)	348 €
40" for size 4	(49 05 018)	426 €



Wheel camber adjustment
to set the hand rim width

	Piece
0° for size 3+4 (49 00 019)	incl.
3° for size 2-4 (49 01 019)	incl.
6° for size 2-4 (49 02 019)	incl.
9° for size 1-4 (49 03 019)	incl.
12° for size 1+2 (49 04 019)	incl.
15° for size 1 (49 05 019)	incl.

Accessories

Footplate – height and depth adjustable



Split footplate – individual height, depth, width and angle adjustable

already included
in basic model
„with abduction“



Heel edges with ankle straps – angle, width and depth adjustable, for mounting on the footplate(s) with magnetic locks



	Front width
Size 1	8.5 – 11 cm
Size 2	10 – 12.5 cm

	Price / piece		Upcharge / pair		Piece / pair
for size 1 (49 01 020)	incl.	for size 1 (49 01 021)	439 €	Size 1 (49 01 027)	203 €
for size 2 (49 02 020)	incl.	for size 2 (49 02 021)	439 €	Size 2 (49 02 027)	203 €
for size 3 (49 03 020)	incl.	for size 3+4 (49 03 021)	439 €		
for size 4 (49 04 020)	incl.				