



Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The bathing and shower cradle.

- UV-resistant plastic frame with removable cover
- Seat and back section with adjustable angle
- Back angle adjustable with practical single hand operation
- Two heights possible via adjustable feet



HMV-No.: 04.40.04.0007



Basic model

wave. size 1-3
standard



Price / piece

Size 1 (58 01 000)	944 €
Size 2 (58 02 000)	1.000 €
Size 3 (58 03 000)	1.066 €

wave. size 1-3

with angle adjustable calf rest



Price / piece

Size 1 (58 01 100)	1046 €
Size 2 (58 02 100)	1.102 €
Size 3 (58 03 100)	1.178 €

Basic model incl. Accessoriespaket

wave size 1-3
standard, incl. headrest,
pelvis strap and leg straps



Price / piece

Size 1 (58 01 050)	1.056 €
Size 2 (58 02 050)	1.234 €
Size 3 (58 03 050)	1.306 €

wave size 1-3

with angle adjustable calf rest,
incl. headrest, pelvis strap and
leg straps



Price / piece

Size 1 (58 01 150)	1.158 €
Size 2 (58 02 150)	1.341 €
Size 3 (58 03 150)	1.413 €

Cover colours

Blue



Pink



Price / piece



Price / piece

9100 300

incl.

9100 301

incl.

Dimensions					
		Size 1	Size 2	Size 3	
A	Seat depth	27 cm	33 cm	43 cm	
B	Seat width	38 cm		44 cm	
C	Back height	58 cm	67 cm	81 cm	
D	Seat height	with plastic frame	5 cm / 18 cm		
		with bath cradle subframe	32 - 47 cm		
		with shower cradle subframe	70 cm		
		with transfer subframe	53 - 75 cm		
E	Width	46 cm		53 cm	
	Length of bathing cradle, completely flat, without calf pad	89 cm	104 cm	128 cm	
	Length of bathing cradle, completely flat, with calf pad	108 cm	124 cm	155 cm	
	Length, completely flat	75 cm			
	Height of bath subframe	20 cm			
	Length of bath subframe	99 cm			
	Width of bath subframe	44 cm			
	Height of shower subframe	56 cm			
	Width of shower subframe	69 cm			
	Length of shower subframe	86 cm			
	Max. load	34 kg	68 kg	114 kg	
	Back and seat angle	0°, 15°, 30°, 45°, 60°, 75°			

Subframe		
Bath subframe wave. to increase the bathing height 	Shower cradle subframe wave. mobile with 4 x 100 mm dual castors with brakes 	Transfer subframe wave. – mobile with 4 castors (2 with brakes), pivots to the right (by means of conversion also pivotable to the left), foldable, incl. wall hook to hang it on the wall. (we recommend the use in combination with a bathing and shower cradle with calf rest) 
Price / piece	Price / piece	Price / piece
58 00 500 664 €	58 00 510 1.412 €	58 00 520 1.636 €

Accessories			
Headrest / torso support 	Chest strap / pelvis strap 	Chest strap with lateral positioning 	Leg straps 
Price / piece	Price / piece	Price / piece	Price / piece
for size 1+2 (58 01 001) 149 € for size 3 (58 03 001) 149 €	for size 1 (58 01 002) 95 € for size 2 (58 02 002) 95 € for size 3 (58 03 002) 95 €	for size 1 (58 01 003) 118 € for size 2 (58 02 003) 118 € for size 3 (58 03 003) 118 €	for size 1 (58 01 004) 95 € for size 2 (58 02 004) 95 € for size 3 (58 03 004) 95 €