



50 **corzo xcountry go.**

Invoice address:

Dispatch address:

See invoice address

**Order
Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

I am your foldable rehab buggy.

- Tested according to the requirements of EN ISO 7176/19
- Steel / aluminum frame with footbrake
- Steerable and puncture-resistant front wheels
- Handles
- Adjustable seat depth and back height
- Height adjustable, foldable and removable footplate
- 5-point harness with padded crotch strap
- Lateral guide
- Basket
- Soft top
- Cover: black
- Frame: anthracite


HMV-No.: 18.99.01.1032

patron®



Dimensions

		Size 1	Size 2	Size 3	Size 4
A	Backrest height	55 cm	67 cm	65 / 85 cm	69 / 88 cm
B	Seat depth	28 - 34 cm	29 - 35 cm	32 - 39 cm	33 - 41 cm
	Seat width	30 cm	34 cm	38 cm	42 cm
C	Lower leg length	29 - 42 cm	32 - 46 cm	32 - 46 cm	33 - 50 cm
	Backrest angle	90° / 89°	94° / 91°	90° / 88°	94° / 92°
	Seat angle	20°	25°		21°
D	Total length	98 cm	102 cm		105 cm
	Total width	52 cm	57 cm	62 cm	65 cm
E	Total height	101 cm	105 cm	117 cm	121 cm
	Weight	12,7 kg	13 kg	14 kg	14,7 kg
	Max. load	45 kg	55 kg	75 kg	75 kg
	Length (folded)	117 cm	122 cm		126 cm
	Width (folded)	32 cm	31 cm	32 cm	37 cm
	Height (folded)	34 cm	36 cm		



Basic model

corzo xcountry go. size 1	corzo xcountry go. size 2	corzo xcountry go. size 3	corzo xcountry go. size 4
Price / piece	Price / piece	Price / piece	Price / piece
29 01 100 / 1 2.578 €	29 02 100 / 1 2.578 €	29 03 100 / 1 2.721 €	29 04 100 / 1 2.939 €

Accessories

Footplate – incl. calf strap, height adjustable, foldable and removable (standard configuration, if nothing else stated) 	Split footplates – incl. calf strap, height adjustable, foldable and removable 	Headrest 	Backrest extension 
Please choose either footplate or split footplates.			
Price / piece	Piece / Paar	Price / piece	Price / piece
for size 1 (29 01 001) incl. for size 2 (29 02 001) incl. for size 3 (29 03 001) incl. for size 4 (29 04 001) incl.	for size 1 (29 01 002) incl. for size 2-4 (29 02 002) incl.	29 00 003 203 €	for size 3 (29 03 012) 137 € for size 4 (29 04 012) 137 €

Accessories

Long lateral guide 	Chest-shoulder strap 	Positioning jacket 
Piece / Paar	Price / piece	Price / piece
for size 1 (29 01 105) 98 € for size 2-4 (29 02 105) 98 €	for size 1 (29 01 006) 162 € for size 2-4 (29 02 006) 162 €	29 00 007 197 €

Accessories

Pelvic harness 	Arm supports 	Rain cape – can only be used in conjunction with the soft top 
Price / piece	Piece / Paar	Price / piece
Size 1 (29 01 008) 104 € Size 2 (29 02 008) 104 €	29 00 013 207 €	for size 1+2 (29 01 010) 87 € for size 3+4 (29 02 010) 87 €