



marie integral.

Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

I am your safe companion on the go.

- 5-point belt incl. support pad and shoulder pad
- Pelvic support pad
- Headrest with additional wedges
- Height and width-adjustable lateral pads
- With ISOFIX bracket
- Innovative Side Impact Buddy for additional protection
- With four different sitting and resting positions
- Easily adjust the height of the headrest/backrest



HMV-No.: 26.11.06.0035

Approval no. 030245 S

**ECE R129/03
CERTIFIED**



**FOR THE SMALL AND
NOT SO SMALL**
61-105 cm, up to 16,5 kg



360° SECURITY
FORWARDS AND BACKWARDS




ECE R129/03
Certified for the highest safety
requirements



Dimensions

Recommended body size	61-105 cm
Recommended body weight	Max. 16.5 kg
Length	76 cm
Width	43 cm
Height	61 cm
Seat height	12.5 cm (lowest point to the horizontal plane in the seat)
Seat depth	28 cm
Seat width	28 cm
Armrest height	34 cm
Back height	Max. 41 cm
Backrest angle	Backwards: 34.5°-49.5°; Forward: 59.5° - 74.5° measured to horizontal plane
Headrest height	19 cm
Footrest height	-2 cm to + 16 cm to the horizontal surface
Weight	15.6 kg
Max. Load	16.5 kg
Fastening the safety device	5-point belt



Basic model

marie integral.

- 5-point belt incl. support pad und shoulder pad
- Pelvic support pad
- Headrest with additional wedges
- Height and width-adjustable lateral pads

Accessories

Abduction wedge



Arm support pad with additional abduction wedge



	Price / piece
12 01 200	2.480 €

	Price / piece	Piece / set
12 01 201	120 €	12 01 202 340 €

Accessories

Footrest – height adjustable



Sun protection



Car seat protector – Compatible with all car seats

- High-quality and non-slip neoprene
- Two practical, child-friendly front pockets
- Easy to attach to the headrest



Footmuff

- Water-repellent and windproof outer material
- Optimum fit up to 105 cm body height
- Easy to wash by hand wash



	Price / piece		Price / piece		Price / piece		Price / piece
12 01 203	280 €	12 01 204	80 €	12 01 205	89 €	12 01 206	48 €

Accessories

Seat minimiser – reduces the seat depth by 3 cm



	Price / piece
12 01 207	72 €