



Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The superior therapy chair.


- Incl. padding and cover
- Incl. back
- Incl. armrests
- Incl. 4 castors (100 mm) with directional locks at the front and brakes at the rear.
- Adjustment options:
 - Height, seat angle and backrest angle adjustment via electric motor with manual control
 - Seat depth and seat width
 - Height-, width- and depth-adjustable armrests
- Cover: Slate grey (No.124)
- Frame colour: Black grey

HMV-No.: 26.11.05.0051




Dimensions

		Size1	Size1b	Size 2	Size 3
A	Seat depth	30 - 42 cm	38 - 48 cm		45 - 55 cm
B	Seat width	28 - 38 cm		36 - 48 cm	
C	Short back height	42 - 60 cm		48 - 60 cm	
C	Long back height	52 - 70 cm		58 - 70 cm	
D	Seat height	40 - 55 cm			
E	Lower leg length	27 - 45 cm		35 - 52 cm	
F	Total width	60,5 cm			
G	Base frame length	84,5 cm			
	Seat angle	(-)20° - 20°			
	Back angle	(-)5° - 25°			
	max. load	125 kg			
	Weight	45 kg			



Presettings



Presettings on the therapy chair at Schuchmann.
Please provide us with the following body measurements:

I Axilla height _____

|| Chest width _____

III Back height _____

IV Pelvic width _____

V Seat depth _____

VI Lower leg length _____

Price

62 00 080

incl.

Basic model

maxi. size 1		maxi. size 1b 		maxi. size 2		maxi. size 3	
Price				Price		Price	
62 01 000	3,502 €	62 01 200	3,502 €	62 02 000	3,502 €	62 03 000	3,708 €

Covers


Minor position

Slate grey (No. 124) 	Soft grey (No. 127) 	Flexi light grey (No. 50) 	Without seat and back upholstery and cover – for the installation of custom cushions
Price	Price	Upcharge	Price reduction
9100 124 incl.	9100 127 incl.	6200 050 272 €	6200 001 Netto -100 €

Headrests

Headrest (standard) incl. standard holder, padded, incl. cover, height and angle adjustable		Headrest (shell-shaped) incl. standard holder, padded, incl. cover, height and angle adjustable		Headrest (lateral guide) incl. standard holder, padded, incl. cover, height and angle adjustable	
Price / piece		Price / piece		Price / piece	
for size 1+1b (62 01 005)	335 €	for size 1+1b (62 01 006)	453 €	for size 1+1b (62 01 007)	391 €
for size 2+3 (62 02 005)	335 €	for size 2+3 (62 02 006)	453 €	for size 2+3 (62 02 007)	391 €

Headrests

Headrest (Lateral guide, angle adjustable) – incl. standard holder, padded, incl. cover, height, depth and angle adjustable, side guides are each angle adjustable  A = Width contact surface B = Height headrest C = Width lateral guide <table border="1"> <thead> <tr> <th></th><th>A</th><th>B</th><th>C</th></tr> </thead> <tbody> <tr> <td>Size 1</td><td>13,5 cm</td><td>10,5 cm</td><td>10 cm</td></tr> <tr> <td>Size 2</td><td>15 cm</td><td>10,5 cm</td><td>13,5 cm</td></tr> <tr> <td>Size 3</td><td>17 cm</td><td>12,5 cm</td><td>16,5 cm</td></tr> <tr> <td colspan="4">Price / piece</td></tr> <tr> <td>Size 1 (62 01 008)</td><td></td><td></td><td>651 €</td></tr> <tr> <td>Size 2 (62 02 008)</td><td></td><td></td><td>651 €</td></tr> <tr> <td>Size 3 (62 03 008)</td><td></td><td></td><td>651 €</td></tr> </tbody> </table>		A	B	C	Size 1	13,5 cm	10,5 cm	10 cm	Size 2	15 cm	10,5 cm	13,5 cm	Size 3	17 cm	12,5 cm	16,5 cm	Price / piece				Size 1 (62 01 008)			651 €	Size 2 (62 02 008)			651 €	Size 3 (62 03 008)			651 €	Adapter for headrests for lateral adjustment Price / piece 62 02 009 154 €	Double joint for headrest mounts to extend the depth adjustment range  Price / set 62 00 009 98 €	Universal adapter for tailored headrests in the version 15 x 15 mm  Price / piece 62 01 009 95 €
	A	B	C																																
Size 1	13,5 cm	10,5 cm	10 cm																																
Size 2	15 cm	10,5 cm	13,5 cm																																
Size 3	17 cm	12,5 cm	16,5 cm																																
Price / piece																																			
Size 1 (62 01 008)			651 €																																
Size 2 (62 02 008)			651 €																																
Size 3 (62 03 008)			651 €																																

Accessories



Back – incl. pad and cover (short back is standard unless otherwise specified)

Back extension unit – incl. pad and cover (selectable for short back and for sizes 2 and 3)

	Price / piece		Price / piece
short for size 1+1b (62 01 002)	incl.	for size 2+3 62 02 023	132 €
short for size 2+3 (62 02 002)	incl.		
long for size 1+1b (62 01 003)	incl.		
long for size 2+3 (62 02 003)	incl.		

Accessories

Lateral pads – padded, incl. cover, angle-, depth and width adjustable



Lateral pads	Height	Depth
Size 1-3	12 cm	16 cm

Lateral pads, long (trunk guidance)
padded, incl. cover, angle-, depth and width adjustable



Lateral pads	Height	Depth
Size 1-3	20 cm	16 cm

	Width 1	Width 2		Width 1	Width 2
Chest width for maxi. size 1+1b	17 – 25 cm	25 – 33 cm	Chest width for maxi. size 1+1b	17 – 25 cm	25 – 33 cm
Chest width for maxi. size 2+3	25 – 33 cm	33 – 41 cm	Chest width for maxi. size 2+3	25 – 33 cm	33 – 41 cm
Price / pair			Price / pair		
Width 1 (62 01 022)	185 €		Width 1 (62 01 010)	185 €	
Width 2 (62 02 022)	185 €		Width 2 (62 02 010)	185 €	

Accessories


Thigh guides – padded, incl. cover, adjustable in angle, depth and width



Removable thigh guide – padded, incl. cover, adjustable in angle, depth and width



Knee abduction wedge – pre-pluggable, padded, incl. cover




	Price / pair		Price / pair		Price / piece
for size 1+1b (Length 22 cm)	(62 02 011) 218 €	for size 1+1b (Length 22 cm)	(62 01 043) 320 €	62 00 014	249 €
for size 2+3 (Length 25 cm)	(62 03 011) 218 €	for size 2+3 (Length 25 cm)	(62 02 043) 320 €		

Accessories

Adjustable pelvic guidance – padded, incl. cover, angle adjustable

	Seat width A
for size 1	28 – 38 cm
for size 2+3	36 – 48 cm



Therapy tray – made of wood with rail, prefixable, fixed on both sides (can only be used with armrests incl. therapy table insertion)



	Height	Depth		Width	Depth
Size 1	12 cm	12 cm		60 cm	50 cm
Size 2	15 cm	16 cm			
Price / piece			Price / piece		
for size 1+1b	62 01 012		452 €		62 00 013
for size 2+3	62 02 012		452 €		294 €

Accessories

Armrests standard – width and height adjustable



	Price / pair
62 00 004	incl.

Accessories

Armrests incl. therapy table insertion – width, height and angle adjustable



	Upcharge / pair
fixed (62 01 004)	247 €
fold-down (62 02 004)	371 €

Armrest padding, size 2 – for armrests incl. therapy table insertion.

(Standard configuration for the Basic model, size 1–2, unless otherwise specified.)

Length: 27 cm Width: 7 cm

62 02 046

Armrest padding, size 3 – for armrests incl. therapy table insertion.

(Standard configuration for the Basic model, size 3, unless otherwise specified.)

Length: 34 cm Width: 8 cm

62 03 046

Accessories


Plug-on footplate – height approx. 13 cm above the floor



Wooden footplate – height and knee angle adjustable, can be folded up, permanently attached to the seat unit and slots for attaching harnesses



Calf plate for footplate – covered with black Flexi fabric (not possible with plug-on footplate)



	Width	Depth		Width	Height
for size 1+1b	41cm	26cm	for size 1+1b	32 cm	6,5 cm
for size 2+3	45cm		for size 2+3	38 cm	
Price / piece			Price / piece		
for size 1+1b (6201016)		525 €	for size 1+1b (6201017)		228 €
for size 2+3 (6202016)		525 €	for size 2+3 (6202017)		228 €

Padding

Footrest padding for the attachable footrest

upholstered in black Flexi fabric, height: 1 cm (not compatible with the foot harness)

Footrest padding for height-adjustable footrest.

upholstered in black Flexi fabric, height: 1 cm (not compatible with the foot harness)

Price / piece		Price / piece	
62 00 038	165 €	for size 1+1b (62 01 038)	133 €
		for size 2+3 (62 02 038)	133 €

Accessories

Castors with locking devices – (standard with 2 locking devices at the rear, unless otherwise specified – with 4 locking devices, a directional lock is not possible)



Central parking brake for rear castors

parking lever optionally mounted on the right or left



Upcharge / set		Upcharge / set	
100 mm castors 2 with locking devices (62 00 018)	incl.	mounted left (100 mm castors)	(62 00 019) 474 €
100 mm castors 4 with locking devices (62 01 018)	101 €	mounted right (100 mm castors)	(62 01 019) 474 €
125 mm castors 2 with locking devices (62 00 041)	42 €	mounted left (100 mm front castors with locking device)	(62 02 019) 567 €
125 mm castors 4 with locking devices (62 01 041)	121 €	mounted right (100 mm front castors with locking device)	(62 03 019) 567 €
		mounted left (125 mm castors)	(62 00 042) 516 €
		mounted right (125 mm castors)	(62 01 042) 516 €
		mounted left (125 mm front castors with locking device)	(62 02 042) 609 €
		mounted right (125 mm front castors with locking device)	(62 03 042) 609 €

Accessories

Push bar – with rubber cover



Set for tool-free adjustment

Headrest, thoracic supports, back height, knee angle, seat depth, thigh guide, and lower leg length are adjustable without tools.

Price / piece		Price / set	
for size 1+1b (62 01 021)	216 €	62 00 045	78 €
for size 2+3 (62 02 021)	216 €		

Accessories

Harness assembly kit for the seat

consisting of 2 clip buckles and 2 clip buckles (drilled)

Harness assembly kit for the back

consisting of 2 threaded plates with clip buckles and 2 clip buckles (drilled)

Price / set

Price / set

25mm **(62 01 030)**

48 €


25mm **(62 01 031)**

48 €

Harness

Chest-shoulder strap – made of elastic, breathable and moisture repellent material, three straps with clips, sewn strap to the left shoulder; colour: black

	Size 1	Size 2	Size 3
A	13 cm	18 cm	18 cm
B	19 cm	20 cm	21 cm
C	15 cm	19 cm	28 cm
D	5 cm	5 cm	5 cm
E	28 cm	34 cm	44 cm
Strap width	2 cm	2,5 cm	2,5 cm



Price / piece

Size 1 **(62 01 032)**

197 €

Size 2 **(62 02 032)**


197 €

Size 3 **(62 03 032)**

197 €

Positioning jacket – with zipper, made of elastic, breathable material, shoulder straps sewn at back, pelvis straps with clips; colour: black

	Size 1	Size 2
A	11 cm	18 cm
B	21 cm	22 cm
C	12 cm	15 cm
D	3,5 cm	3,5 cm
E	28 cm	34 cm
Strap width	2 cm	2,5 cm



Price / piece

Size 1 **(62 01 033)**


228 €

Size 2 **(62 02 033)**

228 €

Harness

Pelvic harness – with pad and adjustable clip on both sides; colour: black



	Size 1	Size 2
Pad width	4 cm	5 cm
Pad length per side	18 cm	28 cm
Strap length per side	30 cm	40 cm
Strap width	2,5 cm	2,5 cm

Price / piece


Size 1 **(62 01 034)**

99 €

Size 2 **(62 02 034)**

99 €

4-point pelvic harness – with pad and adjustable clip on both sides; colour: black



	Size 1	Size 2
Pad width	4 cm	5 cm
Pad length per side	18 cm	28 cm
Strap length per side	30 cm	40 cm
Strap width	2,5 cm	2,5 cm

Price / piece

Size 1 **(62 01 035)**

172 €

Size 2 **(62 02 035)**

172 €

Harness

Pelvic harness – by „Pelvi.Loc“, colour: black



	Size 2 (S)	Size 3 (M)	Size 4 (L)
for pelvic width	25 – 34 cm	32 – 40 cm	38 – 46 cm
Padding width	7 cm	7,5 cm	8,5 cm
Padding Length per side	19 cm	22 cm	24 cm
Strap Length per side	25 cm	45 cm	45 cm
Strap width	2,5 cm	2,5 cm	2,5 cm

Price / piece

Size 2 (S) **(62 02 039)**

325 €


Size 3 (M) **(62 03 039)**

325 €

Size 4 (L) **(62 04 039)**

325 €

Groin harness (standard) – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black



	Size 2	Size 3	Size 4
A	29 cm	33 cm	36 cm
B	17 cm	22 cm	26 cm
C	43 cm	49 cm	55 cm
Strap width	2,5 cm	2,5 cm	2,5 cm

Size 2 **(62 02 040)**

165 €

Size 3 **(62 03 040)**


165 €

Size 4 **(62 04 040)**

165 €


Harness

Ankle cuffs – with velcro fastener and idler pulley, colour: black (not possible in combination with padded cover for footplate/split footplate)



	Ankle circumference	Strap width
Size 1	12 – 19 cm	2,5 cm
Size 2	15 – 23 cm	2,5 cm
Size 3	19 – 27 cm	4 cm
		Price / pair
Size 1 (62 01 036)		236 €
Size 2 (62 02 036)		236 €
Size 3 (62 03 036)		236 €

Foot positioner – with straps, incl. velcro fastener, colour: black (not possible in combination with padded cover for footplate/split footplate)



	A	Strap width
Size 1	6,5 cm	2,5 cm
Size 2	8 cm	2,5 cm
Size 3	9,5 cm	2,5 cm
		Price / pair
Size 1 (62 01 037)		217 €
Size 2 (62 02 037)		217 €
Size 3 (62 03 037)		217 €