



rifton stander. prone

117

Invoice address:

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Dispatch address:

See invoice address

Commission:

The versatile standing trainer.

- Open frame with height adjustment and push handle
- Angle adjustment of the main column from 0° - 90° using a gas spring
- Swivelling trunk laterals for easy transfer
- Leg support with hip abduction and individual leg length adjustment
- Incl. trunk strap and pelvic harness for secure positioning
- Use in the supine position thanks to the mountable upper body support



Dimensions

	Size 1	Size 2	Size 3
Total length	81,5cm	99 cm	109 cm
Total width	53,5 cm	68,5cm	73,5cm
Horizontal body support height above floor	56 - 71cm	56 - 91,5 cm	60 - 96,5cm
Sandal to knee height	15 - 28cm	25,5 - 38 cm	31 - 48 cm
Knee to hip height	15 - 25,5cm	20,5 - 33 cm	30,5 - 44 cm
Sandal to top of headrest	78,5 - 109cm	104 - 139,5 cm	124,5 - 169 cm
Distance between trunk lateral	12,5 - 25,5cm	18 - 28 cm	23 - 35,5 cm
Distance between thigh medial and lateral support	5 - 10 cm	7,5 - 12,5 cm	10 - 14 cm
Basic item weight	21kg	32,5kg	52 kg
Body size	71 - 106cm	96 - 134 cm	122 - 165 cm
Max load	34 kg	45kg	84 kg



Basic model

rifton stander, prone size 1 For users with a height of 71 – 106 cm Frame colour: red		rifton stander, prone size 2 For users with a height of 96 – 134 cm Frame colour: blue		rifton stander, prone size 3 For users with a height of 122 – 165 cm Frame colour: blue	
				 	
Price / piece		Price / piece		Price / piece	
45 01 200 4.190 €		45 02 200 4.890 €		45 03 200 5.690 €	

Accessories

Trunk strap – Standard, if not stated otherwise				Trunk strap with rotation control			
							
Price / piece initial equipment		Price / piece retrofitted		Upcharge / piece initial equipment		Price / piece retrofitted	
for size 1	45 01 117	incl.	45 01 117ET 383 €	for size 1	45 01 118	68 €	45 01 118ET 452 €
for size 2	45 02 117	incl.	45 02 117ET 395 €	for size 2	45 02 118	78 €	45 02 118ET 474 €
for size 3	45 03 117	incl.	45 03 117ET 405 €	for size 3	45 03 118	88 €	45 03 118ET 495 €

Accessories

Knee cuffs		Additional set of trunk laterals (can be used as laterals at chest or hip level)		Medial thigh supports		Lateral thigh supports	
							
Price / pair		Price / pair		Price / pair		Price / pair	
for size 1+2	(45 01 107) 258 €	for size 1+2	(45 01 104) 258 €	for size 1+2	(45 01 108) 258 €	for size 1+2	(45 01 109) 258 €
for size 3	(45 02 107) 279 €	for size 2+3	(45 02 104) 267 €	for size 2+3	(45 02 108) 267 €	for size 2+3	(45 02 109) 267 €

Accessories

Legrest incl. abduction adjustment Standard, if not stated otherwise				Legrest incl. abduction adjustment, knee contracture compensation and hip & knee flexion			
							
Price / pair initial equipment		Price / pair retrofitted		Upcharge / pair initial equipment		Price / pair retrofitted	
for size 1	45 01 101	incl.	45 01 101ET 826 €	for size 2	45 02 102	387 €	45 02 102ET 1.233 €
for size 2	45 02 101	incl.	45 02 101ET 851 €	for size 3	45 03 102	397 €	45 03 102ET 1.257 €
for size 2	45 03 101	incl.	45 03 101ET 879 €				

Accessories

Upper body support		Sandal lift kit (only for size 1 available)		Tray – angle and height adjustable, removable	
					
Price / piece		Price / pair		Price / piece	
for size 3	(45 03 222) 312 €	for size 1	(45 01 111) 126 €	for size 1+2	(45 01 110) 711 €
				for size 3	(45 03 110) 733 €

Extension for use in the supine position



Upper trunk support with a headrest – absolutely essential for extended use in the supine position!

Head laterals



	Price / piece		Price / pair
for size 1 (45 01 206)	435 €	for size 1+2 45 01 203	258 €
for size 2 (45 02 206)	447 €	for size 2+3 45 02 203	267 €
for size 3 (45 03 206)	459 €		

Extension for use in the supine position

Hip strap – Standard, if not stated otherwise



Hip strap with rotation control



One of these options is mandatory for extended use in the supine position!

	Price / piece		Price / piece
for size 1 (45 01 220)	383 €	for size 1 (45 01 221)	452 €
for size 2 (45 02 220)	395 €	for size 2 (45 02 221)	474 €
for size 3 (45 03 220)	405 €	for size 3 (45 03 221)	495 €

Extension for use in the supine position

Butterfly harness



	Price / piece
for size 1 (45 01 219)	452 €
for size 2 (45 02 219)	474 €
for size 3 (45 03 219)	495 €

Extension for use in the supine position

Knee supports



Swivelling knee supports
(only for size 2 available)



One of these options is mandatory for extended use in the supine position!

	Price / pair		Price / pair
for size 1 (45 01 205)	254 €	for size 2 (45 02 212)	532 €
for size 2 (45 02 205)	266 €	for size 3 (45 03 212)	562 €
for size 3 (45 03 205)	279 €		