



Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:


The dynamic standing trainer.

- Open central column with height adjustment via a hand-tightened screw
- Supplied with a flexible pelvic pad and a protective pad to cover the manual adjustment screw
- Dynamic centre column:
 - Tilt angle up to 10° in all directions
 - Individual resistance adjustment (from grade 1-4) on each side
- Universal mount for accessories
- Can be used in anterior and posterior positions
- 4 castors (75 mm) with brakes
- Base plate made of birch with profiles to accommodate the heel edges
- Optional mounting point for assisted communication devices.



HMV-No.: 28.29.01.1061

Dimensions



		Size 0 (mini)	Size 1	Size 2
	Body size	65 cm - 95 cm	85 cm - 115 cm	105 cm - 140 cm
A	Sternum height	42 cm - 62 cm	58 cm - 86 cm	76 cm - 100 cm
B	Table height from the base plate	48 cm - 72 cm	63 cm - 100 cm	80 cm - 115 cm
C	Total length	80 cm	88 cm	98 cm
D	Total width	50 cm	60 cm	70 cm
	Max. load	20 kg	40 kg	60 kg
	Weight with dynamic centre column	16 kg	19 kg	22 kg
	Weight with fixed centre column	13 kg	15 kg	18 kg

Basic model „rigid“		
trixi. size 0 for users with a height of 65 – 95 cm	trixi. size 1 for users with a height of 85 – 115 cm	trixi. size 2 for users with a height of 105 – 140 cm
Price / piece	Price / piece	Price / piece
Size 0 (63 00 000) 1.350 €	Size 1 (63 01 000) 1.450 €	Size 2 (63 02 000) 1.550 €

Basic model „dynamic“		
trixi. size 0 for users with a height of 65 – 95 cm	trixi. size 1 for users with a height of 85 – 115 cm	trixi. size 2 for users with a height of 105 – 140 cm
Price / piece	Price / piece	Price / piece
Size 0 (63 00 100) 1.780 €	Size 1 (63 01 100) 1.880 €	Size 2 (63 02 100) 1.980 €

Pelvic pads		
Pelvic pad (flex) – adjustable in height and width for rotation control (standard, unless otherwise specified: for basic model size 0 = strap size 0 for basic model size 1 = strap size 1 for basic model size 2 = strap size 2)	Pelvic pad – height and width adjustable with side guides	Retaining strap for pelvic pad with side support
		
Price / piece	Price / piece	
Size 0 (63 00 004) 298 €	Size 0 (63 00 003) 426	for size 0 (63 00 015) 172 €
Size 1 (63 01 004) 298 €	Size 1 (63 01 003) 426	for size 1 (63 01 015) 172 €
Size 2 (63 02 004) 298 €	Size 2 (63 02 003) 426	for size 2 (63 02 015) 172 €

Lateral pads		
Chest pad (flexible) – adjustable in height and width for rotation control	Chest pad – height and width adjustable with side guides	Retaining strap for chest pad with side support
		
Price / piece	Price / piece	
Size 0 (63 00 002) 298 €	Size 0 (63 00 001) 426 €	for sizes 0 (63 00 014) 172 €
Size 1 (63 01 002) 298 €	Size 1 (63 01 001) 426 €	for sizes 1 (63 01 014) 172 €
Size 2 (63 02 002) 298 €	Size 2 (63 02 001) 426 €	for sizes 2 (63 02 014) 172 €

Covers for chest and pelvic pads		
Slate grey (No. 124) 	Flexi Light Grey (No. 50) 	
Price		Upcharge
63 00 124 incl.	63 00 050	272 €

Accessories		
Vest – adjustable in size and height, including safety belt. Supports sensory perception. (only available in slate grey)	Depth adjustment for chest pads	Depth adjustment for pelvic pads
		
Price / piece	Price / piece	Price / piece
Size 0 (63 00 013) 480 €	up to 6 cm (63 00 005) 98 €	up to 6 cm (63 00 017) 98 €
Size 1 (63 01 013) 480 €	up to 10 cm (63 01 005) 98 €	up to 10 cm (63 01 017) 98 €
Size 2 (63 02 013) 480 €		

Accessories

Flexible Knee Positioner (FKP) Standard



Heel edges incl. strap



Price / pair

Price / pair

Size 0 **(63 00 106)**

490 €

Size 0 **(63 00 012)**

225 €

Size 1 **(63 01 106)**

490 €

Size 1 **(63 01 012)**

225 €

Size 2 **(63 02 106)**

490 €

Size 2 **(63 02 012)**

225 €

Accessories

Tray – height and depth adjustable, with edge


Price / piece

for sizes 0 **(63 00 008)**

425 €

for sizes 1+2 **(63 01 008)**

425 €

Armrests for tray



Price / piece

for size 0 **(63 00 010)**

215 €

for size 1+2 **(63 01 010)**

215 €

Accessories

Bracket for assisted communication for mounting on the tray


Price / piece

63 00 011

110 €

Bracket for assisted communication for mounting on the centre column


Price / piece

63 01 011

110 €

Dynamic foot cushion for sensory stimulation – incl. cover and fastening material


Price / piece

for size 0 **(63 00 018)**

725 €

for size 1 **(63 01 018)**

725 €

for size 2 **(63 02 018)**

725 €

Accessories

Angle adjustment of the base plate tilt adjustable in increments to 2.5°, 5° or 7.5°


Price / piece

63 00 019

380 €

stixi. – The innovative app with a standing diary for documentation and child-friendly game.

The standing diary enables structured recording of standing times.

The integrated game offers playful motivation for children.

EN-Version coming soon in the Apple App Store and Google Play Store for smartphones and tablets.

