

GMFCS  
Level III-V



ECE R129/03  
ZERTIFIZIERT

i-size  
universal  
ISOFIX

### I am your safe companion on the go.

And I am your best friend for relaxed journeys.

With my ergonomic seat, you can sit comfortably even on long trips. My unique EPP shell structure provides enhanced impact protection, and together with the side impact buddy, which absorbs up to 40% of impact energy, you are optimally protected. To make sure you can enjoy car journeys together for as long as possible, you can use me in both forward-facing and rear-facing positions. Thanks to my Isofix function, installing me on the back seat is child's play.

And your parents will love my 360-degree rotation – they can simply turn me to lift you in and out with ease.





**The EPP shell structure**  
of the seat offers improved  
impact protection.

**Lateral pads**  
are adjustable in width.

**Isofix mount**  
for use in  
in both directions.

**360° Rotation**  
for comfortable  
installation and removal.

**Headrest**  
can be adjusted in height.

**With the 5-point harness**  
your child will remain safely in place.

**Abduction wedge**  
enables optimum positioning.

**Adjustable tilt angle**  
for 4 different sitting and reclining positions.

**Foot rest**  
to relieve the legs,  
the back and neck.

**Height adjustment of the footrest**  
the car seat grows with the child.

**Rear facing  
travel position.**

**Forward facing  
travel position.**

**Comfortable and  
easy to get into  
and out of.**

**0 to  
16.5 kg**

**9 to  
16.5 kg**







chest pads  
are customisable.

| data                        | marie integral.                                                                          |
|-----------------------------|------------------------------------------------------------------------------------------|
| Recommended body size       | 61–105 cm                                                                                |
| Length                      | 76 cm                                                                                    |
| Width                       | 43 cm                                                                                    |
| Height                      | 61 cm                                                                                    |
| Seat height                 | 12.5 cm (lowest point to the horizontal plane in the seat)                               |
| Seat depth                  | 28 cm                                                                                    |
| Seat width                  | 28 cm                                                                                    |
| Arm rest height             | 34 cm                                                                                    |
| Back rest height            | Max. 41 cm                                                                               |
| Tilt-in-space               | Rear-facing: 34,5°–49,5°;<br>Forward-facing: 59,5° – 74,5° measured to horizontal plane. |
| Headrest height             | 19 cm                                                                                    |
| Footrest height             | –2 cm to + 16 cm from the horizontal surface                                             |
| Weight                      | 15,6 kg                                                                                  |
| Max. load                   | 16,5 kg                                                                                  |
| Fastening the safety device | 5-point harness                                                                          |
| HMV-No.                     | 26.11.06.0035                                                                            |



Optional abduction  
wedge for better leg  
positioning.



You can find more information on accessories  
in the price list.

