

GMFCS
Level III-V



Comfort and safety that last through your teen years.

And I am your companion for exciting trips, helping you become more independent along the way. I'm simple and intuitive to use, so you can quickly learn how to fasten your seat belt all by yourself. This not only boosts your confidence, but also means your parents can relax instead of trying to fasten your belt from an awkward angle! My unique EPP shell structure provides enhanced impact protection. And together with the Side Impact Buddy, which absorbs up to 40% of impact energy, you stay optimally protected wherever your adventures take you.


NEW:
now with
HMFV-No.



Version with
lap-belt and pelvic
support pad.



Version with
5-point harness.



Version with
5-point harness
with armrest and
abduction wedge.



charlie integral.

Head rest

can be adjusted in height with just one hand movement.

Lateral pads

stabilise the upper body.

5-point-harness

holds the child by its most stable points.

Side Impact Buddy

offers more side impact protection.

Arm support

including small armrests.

Isofix-mount

for particularly safe and easy use.

Foot rest

for a stable hold and a relaxed sitting position.



Guided belts

increase independence by allowing people to fasten their own seat belts.

Side Impact Buddy

absorbs 40– % of the force in a side impact.





data	charlie integral.
Recommended height	100–150 cm
Length	48 cm
Width	43–51 cm
Height	59–79 cm
Seat height	11 cm (lowest point to horizontal plane in seat)
Seat depth	27–32 cm (short / long seat)
Seat width	30–37 cm (back / front)
Back height	51–71 cm
Headrest height	20 cm
Lower leg length when using the footrest	21–42 cm
Weight	8 kg
Max. load.	45 kg
HMV-No.	26.11.06.0036



The optional footrest
is ergonomically shaped and
height-adjustable.



You can find more information on accessories
in the price list.

