

activity chair. Simply moving.

GMFCS
Level II-IV


max. 113




(-)15°-20°


(-)15°-25°


20-51 cm


24-61 cm


18-43 cm



I am your dynamic therapy chair.

I really hate sitting still! I'd much rather be changed into lots of different positions. I can tilt forwards or backwards, up or down - as dynamically as you require. Active sitting has a major advantage: it helps you to relieve physical and psychological tension. My wide range of accessories such as headrests, lateral pads or side guides provide the necessary stability in any situation.





activity chair.

Push handle
for convenient transfer.

Lateral pads
facilitate a stable and upright sitting position.

Groin harness
positions the pelvis and encourages proper alignment of the spine.

Sandals with foot straps
can be adjusted to the user's foot position for optimal joint alignment.

Dual castors with locking devices
facilitate transfer.

Headrest
can be adjusted in depth, height and angle.

Back angle
can be adjusted using a gas pressure spring.

Armrests
can be adjusted in height and angle.

Seat pads
are comfortable and avoid pressure points.

Sub-frame
is particularly sturdy and adjustable in height.

Foot support
can be folded away to the rear (incl. calf plate) and is therefore ideal for independent transfer.

Optional dynamic function

gently absorbs your movements and can be easily switched on and off using the rotary handle without requiring tools.



Frame colour  Effect black



		Size 0	Size 1	Size 2
Seat depth		20 - 30 cm	28 - 41 cm	38 - 51 cm
Seat width		30 cm	36 cm	43 cm
Seat width (with hip support pads)		18 - 23 cm	22 - 29 cm	28 - 36 cm
Chest width		17 - 29 cm	17 - 30 cm	24 - 36 cm
Back height		32 - 39 cm	39 - 48 cm	48 - 61 cm
Headrest height - seat		37 - 53 cm	44 - 62 cm	50 - 75 cm
Armrest height - seat		13 - 19 cm	18 - 27 cm	18 - 27 cm
Lower leg length		18 - 30 cm	20 - 38 cm	29 - 54 cm
Seat height over the floor	with short / with long feet	24 - 32 cm / 47 - 55 cm	32 - 42 cm / 48 - 58 cm	42 - 52 cm / 51 - 61 cm
	with short / with long feet and dual castors	34 - 42 cm / 47 - 55 cm	41 - 51 cm / 48 - 58 cm	42 - 52 cm / 51 - 61 cm
	with Hi/Lo base	25.5 - 63.5 cm	30.5 - 58.5 cm	40.5 - 61 cm
Frame width	with short / with long feet	53 / 58 cm	59 / 62 cm	68 / 68 cm
	with short / with long feet and dual castors	57 / 60 cm	62 / 65 cm	71 / 71 cm
	with Hi/Lo base	66 cm	68.5 cm	73.5 cm
Seat angle	with standard base / with Hi/Lo base	(-) 15° - 15° / (-) 15° - 25°		
Back angle		(-) 5° - 20°		
Total length/frame depth		54-73 cm	57-73 cm	77-90 cm
Total width		53-60 cm	59-65 cm	68-70,5 cm
Total height incl. headrest		61-108 cm	76-120 cm	92-136 cm
Total height without headrest		56-94 cm	71-106 cm	90-122 cm
Recommended body height		81-122 cm	102-142 cm	127-188 cm
max. load		34 kg	68 kg	113 kg
Weight of standard base frame		12,2 kg	15,9 kg	24,5 kg
Weight of Hi/Lo frame		25,4 Kg	39 Kg	44,5 kg
HMFV-No. standard base frame		26.11.05.0046		
HMFV-No. Hi/Lo frame		26.11.05.0047		

The dimensions given apply to both seat units with all bases. If certain dimensions differ between the base frames, these are listed specifically.



You can find more information on accessories & harnesses in the price list.

