

GMFCS
Level III-V




I am the therapy chair that grows with you.

I really love being by your side – whether you're a small child or a teenager. That's why I'm available in four different sizes. The great thing about me is that my seat can tilt forwards and backwards, and it's super easy to adjust. I can even tilt far enough for you to lie back and relax.

I can be raised or lowered so you're always right in the middle of whatever's happening.

I can even help you train your trunk muscles – just remove my backrest. It's quick and easy to do.







In just one single step, my back unit can simply be removed. This gives you more freedom to move.





My headrest with angle-adjustable side support is flexible, adaptable, practical and offers great therapeutic benefits.

data	Size 0	Size 1	Size 1b	Size 2
Seat depth*	15 - 25 cm	21 - 30 cm	26 - 36 cm	30 - 42 cm
Seat width	16 - 26 cm	20 - 30 cm		28 - 38 cm
Back height with extra short back	32 - 35 cm		35 - 43 cm	
Back height with short back	35 - 43 cm		42 - 52 cm	
Back height with long back	42 - 52 cm		53 - 63 cm	
Back height with Active back	35 - 43 cm			42 - 52 cm
Seat height	47,5 - 58 cm			
Lower leg length	12 - 22 cm** / 15 - 25 cm	16 - 29 cm** / 19 - 32 cm	19 - 33 cm** / 22 - 36 cm	25 - 43 cm
Total width	45 cm		53 cm	
Length of base frame	66 cm		77 cm	
Seat angle	(-)5° - 37°			
Back angle	(-)5° - 25°			
Max. load	30 kg		45 kg	50 kg
Weight	16 kg	17 kg	18 kg	19 kg
HMV-No.	26.11.05.0021			
* Reduced by 1- cm with adjustable pelvic guidance				
** This adjustment range can be achieved by moving the bracket.				



You can find more information on fabrics & embroidery
on pages 156 - 159.

You can find more information on accessories & harnesses
in the price list.

