

GMFCS
Levels I-III

kg
max. 60



°
(-) 3° - 17,5°

°
(-) 15° - 15°

cm
16 - 40

cm
17 - 40



I am your seating aid for everyday use.

I promote active sitting and support your development. The great thing is whatever size you are, I will grow with you! You can simply adjust my seat and footplate up or down so that you can sit comfortably at the table with your brothers or sisters. I'm made from beautifully crafted wood, with lots of attention to detail, so I fit perfectly with the furniture in your home. And with my wide range of accessories, I can be adapted to suit your individual needs.



You decide how I should look. Combine my frame in natural beech, white lacquer or slate grey with a selection of great coloured covers!



Height adjustment made easy!

Simply press the button and the seat and footplate can be moved up or down.



New armrests!

Made from glass fibre reinforced plastic, making them as strong as an elephant.



Headrest

supports the head control.

Backrest

can be adjusted in height and angle.

Armrests

can be adjusted in height.

Seat plate

can be adjusted in height, depth and angle.

Base foot

with plastic sliders.

Removable pad cover

made of breathable material.

Thorax pelotte pads

facilitate a stable and upright sitting position.

Thigh guide

provides the necessary support.

Wooden frame

with integrated aluminium profiles.

Foot plate

can be adjusted in height, depth and angle.

Fast height adjustment

via the aluminium profiles integrated in the frame.

**Frame colours**

- Natural beech
- White beech
- Beech slate grey



		Size 0 KITA	Size 0	Size 1 KITA	Size 1	Size 2
Seat width		17 - 28cm		20 - 35cm		30 - 40cm
Seat depth		16 - 26cm		25 - 35cm		30 - 40cm
Seat height max.		42 cm	62 cm	42 cm	62 cm	
Back height*		25 cm		29cm		34 cm
Back angle		(-) 3° - 17,5°				
Seat angle		(-) 5° - 15°				
Footrest quick height adjustment min.**	Standard	12cm				
	with 50 mm dual castors	17cm				
	with 75 mm dual castors	20cm				
Height of the footrest at floor level min.***	Standard	2,5cm				
	with 50 mm dual castors	7,5cm				
	with 75 mm dual castors	10,5cm				
Lower leg length with quick-adjust footrest min.		14 cm				
Lower leg height with footrest low to the ground min.		18,6cm				
Total width		46 cm	47 cm	53 cm	54 cm	59 cm
Length of base frame		63 cm	69 cm ***	63 cm	69 cm ***	69 cm ***
Max. load		30 kg		45 kg		60 kg
Weight		9,2kg	9,7 kg	9,6kg	10,2kg	12 kg
HMV-No.		26.11.05.0034				
* Minimum back height measured from seat surface to upper edge of back ** adjustable lower leg length max. = selected/required seat height - footrest height min. ***long stand + 10 cm						



You can find more information on fabrics & embroidery on pages 156 - 159.

You can find more information on accessories & harnesses in the price list.

