

GMFCS
Level III-V



I am the seat system which grows with you.

My dearest wish is to accompany you for a long time to come. For this reason, I am available in four different sizes which grow along with you. You can go anywhere with me, whatever the occasion - I am so flexible that I fit onto lots of different sub-frames. Due to my flexibility you can adjust my pelvic guidance as you wish. Or simply remove my back unit.

Headrest

can be adjusted in height, depth and angle.

Back unit

optionally removable, can be adjusted in angle and height.

Lateral pads

facilitate a stable and upright sitting position.

Armrests

can be adjusted in position.

Thigh guide

supports the movement of the upper thighs outwards.

Abducting leg guides

optionally available.



Frame colours

- Lazurite blue metallic
- Pure white
- White aluminium
- Anthracite metallic
- Dormant red



With diverse
accessories!



Data	Size 0	Size 1	Size 1b	Size 2
Seat depth*	15 - 25 cm	21 - 30 cm	26 - 36 cm	30 - 42 cm
Seat width	16 - 26 cm	20 - 30 cm		28 - 38 cm
Back height (short back)	35 - 43 cm		42 - 52 cm	
Back height (long back)	42 - 52 cm		53 - 63 cm	
Back angle	(-)5° - 25°			
max. load	30 kg		45 kg	50 kg
Weight	9 kg	10 kg		12 kg
Recommended for sub-frame(s)	freddy. size 0	freddy. size 0	freddy. Size 1 + Anti-tip mechanism, franka. Size 2, fridolin maxi.	freddy. Size 1 + Anti-tip mechanism, franka. Size 2, fridolin maxi.
HMV-Nr.	26.11.04.0xxx			

* reduced by 1 cm with adjustable pelvic guides



You can find more information on fabrics
on pages 176 - 179.

You can find more information on accessories & harnesses
in the pricelist



Scan this QR code
and watch videos about
the product directly on our
homepage.

