

GMFCS
Level IV-V


kg
max. 100


Home


(-)5° - 30°


44 - 56 cm

Support head

for all conventional seat pans and systems.


Foot pedal

facilitates inclination of the seat unit from (-)5° to 30°.

Dual castors

for easy transfer - with locking devices at the rear to prevent inadvertent rolling away.

Central column
can be adjusted up to 12 cm in height.



Height adjustment via adjustment wheel



12 cm

I am your robust sub-frame. As my name implies, I am the biggest of us three brothers. That's why it is really easy for me to bring adolescents and adults weighing up to 100 kg from A to B. To get you up to the right height, my seat height can be adjusted using the adjustment wheel.

Frame colour



Anthracite metallic



With a great range of accessories!

Two gas pressure springs
of 200 N each can be operated
by foot pedal



Castors
with 100 mm
width.



Footrest
can be plugged on and is
approx. 8 cm above the ground.



	fridolin maxi.
Height of seat pan support*	44-56 cm (46 - 58 cm**)
Seat angle	(-)5° - 30°
Length of base frame	80 cm
Total width	66 cm
Thickness of the gas pressure spring for tilt adjustment	2 x 200 N
max. load	100 kg
Weight	11,5 kg
HMV No.	26.99.01.0032

* with 75 mm castors ** incl. trapezoidal adapter



You can find more information on accessories
in the price list.

