

GMFCS
Level IV-V



I am your compact seat pan sub-frame. You can go anywhere with me, whether to school or to friends. I am really manoeuvrable and very lightweight. You can adjust my height precisely according to your wishes. If you like, you can rest your feet on my adjustable supports. The consciously designed short wheelbase creates maximum space for footrests and for your legs. And do you know what else is good about me? I have been tested in accordance with ISO 7176: Part 19 and therefore can also accompany you on car rides.



Individual footrests
with calf plates

can be adjusted in height, width and knee angle, and can also be swivelled away.



Divided footrests

are adjustable in height and width and can be swivelled away.





**The scale applied to
the inner tubes**
allows you to determine
the set seating height.

Back
can be folded down onto
the seat frame.

Spokes and clothing guard
prevent hands or clothing
getting into the spokes.

Rear wheels
available in 16", 20", 22"
and 24".

Back angle adjustment
takes place via a snap joint
or through central quick
adjustment.

**Scale on the
inner tubes**
for height adjustment
of the seat adaptor.

Trapezoidal adaptor
for all conventional seat pans.

Divided footrests
can be adjusted in height
and in angle.

Tilting device
makes it easy to
overcome kerbs.

**Swing-away
anti-tip mechanism**
prevents you from
tipping backwards.

Frame colours

 Lazurite blue metallic

 White aluminium

 Dormant red

 Pure white

 Anthracite metallic

 Black - satin finished

 Turquoise blue



Data		Size 0	Size 1	Size 2	Size 3
Seat pan support height	16"	45 – 55 cm (47 – 57 cm*)			
	22"	47 – 57 cm (49 – 59 cm*)			
	24"	50 – 60 cm (52 – 62 cm*)			
Short seat depth		26 – 28,5 cm	29 – 36,5 cm	-	-
Standard Seat depth		30,5 – 38 cm	33 – 40,5 cm	42,5 – 50 cm	
Tilt-in-space		(-) 5° – 30°			
Width between rear wheels (with clothing guard) = Max. outer width of seat pan		39 cm	44 cm	49 cm	54 cm
Frame width inside		28 cm	33 cm	38 cm	43 cm
Frame width outside		35 cm	40 cm	45 cm	50 cm
Seat frame width		25 cm	30 cm	35 cm	35 cm
Total length**		63 – 79 cm		66 – 82 cm	71 – 87 cm
Total width		58 cm	63 cm	68 cm	73 cm
Max. load		90 kg		120 kg	
Weight		ab 14 kg			
HMV No.		26.99.01.1033 (with 16" rear wheels) / 26.99.01.1079 (with 20" to 24" rear wheels and hand rim)			

* incl. trapezoidal adapter ** Depending on wheel size



You can find more information on accessories
in the price list.

