

GMFCS
Level III-V



I am your car seat that pivots.

I love taking a ride in cars. Do you too? Then please take me with you! With me, you are certain to arrive safely. I have been tested and found to be absolutely safe - on the passenger seat and on the rear seat. To make it very easy for your parents to position you on my seat, it is child's play to pivot me by 90 degrees out of the car. And my best feature of all is that I can be used as the basis for special constructions, which means I can be equipped according to your requirements. There aren't many other systems around which can do that.



With a great
range of
accessories!



Seat cover

can be removed.

Lateral pads

facilitate a stable and upright sitting position.

Back

can be tilted in three steps.

Swivel device

can be pivoted by 90°.

Padded headrest

can be adjusted in height.

Straps

help to support and position.

Anatomic side pads

keep the body upright.

Foot support

is adjustable in height and can be removed for transfer.



90°

swivel device
facilitates transfer.

**Removable
lateral pads**





	Size 1		Size 2		Size 3	
	a	b	a	b	a	b
Seat depth	25 cm	30 cm	35 cm	40 cm		
Seat width	28 cm	32 cm	36 cm			
Back height	40 cm	48 cm	56 cm			
Headrest height	59 - 69 cm	65 - 77 cm	76 - 86 cm			
Lateral pad height	25 - 33 cm	32 - 40 cm	36 - 44 cm			
Seat angle*	0° / 15°	0° / 10°				
Chest width	14 - 25 cm	16 - 29 cm	16 - 32 cm			
Back angle		90° - 108°				
Max. load	30 kg	45 kg	60 kg			
HMV No.	26.11.06.0009					

* When using the foot support and tilt adjustment



You can find more information on accessories
in the price list.

