

momo accessories.

Versatile.



Round handlebars with
brake lever ring



Steering movement
limited



Handlebar extension



Standard saddle



Gel-Saddle



Unicycle saddle



Moped saddle

	data	Size 1	Size 2.	Size 3
Standard saddle	Width	15 cm	15 cm	18 cm
	Length	21 cm	24 cm	26 cm
Gel saddle	Rear width	19 cm	24 cm	-
	Front width	4 cm	7 cm	-
	Length	24 cm	27 cm	-

	data	Size 1
unicycle saddle	Rear width	11 cm
	Width at centre	6 cm
	Front width	8 cm
	Length	25 cm
moped saddle	Rear width	26 cm
	Front width	9,5 cm
	Rear length	12 cm
	Total length	25 cm



Continuously adjustable
Crank shortener



Crank shortener



Special crank for
knee contracture

	Size 1	Size 2.	Size 3
Standard seat post	for 12" + 16 / 12"	for 16"	for 20" - 26"
T-saddle post	for 12" + 16 / 12"	for 16" + 20"	for 24" + 24 / 20" + 26"



Holding bracket
in four sizes



With appropriate mounts



data	Size 1	Size 2.	Size 3	Size 4
Max. height*	20 cm	30 cm	37 cm	53 cm
Depth	12 cm	12 cm	12 cm	12 cm

*Measured from a standard saddle up to the top edge of the respective support pad



Classic handlebars



Touring handlebars



Round handlebars



Multifunction handlebars

	data	Size 1	Size 2.	Size 3	Size 4
Classic handlebars	Width	47 cm	49 cm	52 cm	58 cm
Touring handlebars	Width	55 cm	58 cm	-	-
	Depth	17 cm	17 cm	-	-
Round handlebars	Width	40 cm	43 cm	48 cm	-
	Depth	17 cm	21 cm	25 cm	-
Multifunction handlebars	Width	61,5 cm	-	-	-
	Depth	16,5 cm	-	-	-



Foot fixation pedals



Exercise bike pedals



Ergotec pedals



Foot pan



Foot pan with leg guidance



Foot pan with dynamic leg guidance

Dynamic and width adjustable back and pelvic support pads



	data	Size 1	Size 2.	Size 3
Dynamic back support pad	Width	20 - 30 cm	25 - 35 cm	30 - 40 cm
Dynamic pelvic guidance support pad				
Width adjustable back support pad	Width	20 - 28 cm	25 - 33 cm	30 - 36 cm
Width adjustable pelvic guidance support pad				

hand positioning aid



chest strap



positioning vest



	data	Size 0	Size 1	Size 2	Size 3
Manual fixation	Front width	8.5 cm	10 cm	12 cm	14 cm
Chest strap	Width	-	5 cm	6.5 cm	7 cm
	Length	-	18.5 cm	25.5 cm	30.5 cm
Positioning jacket	Inside width	-	8 cm	9 cm	12 cm
	Total width	-	25 cm	30 cm	34 cm
	Side height	-	10 cm	12 cm	16 cm
	Total height	-	28 cm	30 cm	35 cm



Push bar



**Steering for companions
using steering linkage**



**Therapy bicycle frame as assembly set
for standard commercial bicycles**
(hub installation depth min. 120 mm)

