

GMFCS
Level III-IV

kg
max. 100



46–124 cm

29–113 cm



I am your versatile posterior walking aid.

I'm always there for you – even from your very first year of life. Walking is my favorite thing to do! No matter where you want to go, I'll follow you. When you hold onto my handles, you'll quickly feel how I help take some of your body weight. My special abducted frame shape also supports you so that your legs don't swing outwards when you walk. And my pelvic support pad helps you stay upright and walk straight. Pretty cool, aren't I?





Combined pelvic and buttock pad
can become a backrest in one
single easy step.

Forearm rests
with handles.

Padded seat
can easily be folded back.

Aluminium frame
can easily be folded up.

Direction locks
facilitate
walking straight.

Backwards stop
prevents inadvertent
rolling back.

EVA guide wheels
make the walking aid manoeuvrable.



Abducted frame shape
supports you when walking
and prevents you from
swinging outwards.



Aluminium frame
can easily be folded up via
folding joints.

Frame colour
size 0 mini

 Lazurite blue
metallic

Frame colours
size 1

 Telemagenta
 Sky blue

Frame colours
size 2

 Lazurite blue
metallic
 Traffic red

Frame colours
size 3 and size 4

 Anthracite
metallic
 Pure white



data	Size 0 mini*	Size 1	Size 2	Size 3	Size 4
Handle height	29 - 53 cm	34 - 58 cm	51 - 75 cm	68 - 94 cm	87 - 113 cm
Height of forearm supports	-	46 - 70 cm	63 - 87 cm	79 - 103 cm	98 - 124 cm
Handle spacing	22 - 36 cm		30 - 44 cm	38 - 52 cm	43 - 57 cm
Seat height	-	33 cm	45 cm	60 cm	70 cm
Overall width	52 cm	51 cm	61 cm	68 cm	75 cm
Total length	58 cm		70 cm	77 cm	87 cm
Frame height (folded)	37 cm	25 cm	32 cm	34 cm	37 cm
Transport dimensions (L x W x H)	58 x 52 x 37 cm	58 x 51 x 25 cm	70 x 61 x 32 cm	77 x 68 x 34 cm	87 x 75 x 37 cm
Max. load	20 kg	30 kg	50 kg	70 kg	100 kg
Weight	3,9 kg	7,1 kg	8,3 kg	8,9 kg	10,2 kg
HMV-No.	10.46.02.3029				

* Frame not foldable



You can find more information on accessories in the price list.

