

momo tricycle.

Let's go for a little ride.

GMFCS
Levels I-V


kg
max. 120




30-98 cm


16-23 cm

I am your cool tricycle.

I'm ready for adventure whenever you want to take a ride around your neighbourhood or by the lake. It's great fun, and we also train and mobilise your joints, muscles and circulation at the same time. You can choose me in lots of different versions, saddle sizes, handlebar types and colours. And your age doesn't matter, I'm available in many sizes for children, teenagers and adults.





momo tricycle.



Diverse range of accessories
for optimal support and
individual equipment, visit.



**You have the choice! Between lots of
handlebars, saddles, foot pans,
head supports and different drives!**



Frame colours
12", 16" & 18" / 12"

-  Window grey / Pure white
-  Telemagenta (pink) / Pure white
-  Turquoise blue / Pure white

Frame colours
20" & 24" / 20"

-  Window grey / Pure white
-  Pastel blue / Pure white
-  Kiwi sorbet green / Pure white

Frame colours
24", 26" & 26" (XL)

-  Window grey / Pure white
-  Ivory / Pure white
-  / Pure white



data	12"	16/ 12"	16"	20"	24/ 20"	24"	26"	26" (XL)
Pedal to top edge of saddle*	30 - 43 cm	35 - 48 cm	42 - 62 cm	49 - 71 cm	53 - 75 cm	55 - 77 cm	65 - 87 cm	72 - 94 cm
Pedal to top edge of saddle**	33 - 40 cm	38 - 45 cm	45 - 60 cm	52 - 67 cm	56 - 71 cm	58 - 81 cm	68 - 91 cm	75 - 98 cm
Seat post tube up to stem tube	33 cm	40 cm	39 cm	45 cm	54 cm	52 cm	58 cm	
wheelbase	68 cm	75 cm	80 - cm	91 cm	99 cm	104 cm	112 cm	
entry level	16 cm		18 cm	20 cm		22 cm	23 cm	
Overall length	102 cm	110 cm	120 cm	140 cm	153 cm	165 cm	178 cm	
Overall width	58 cm	65 cm		69 cm		75 cm		
Max. load	40 kg	50 kg	60 kg	80 kg	90 kg	100 kg	120 kg	
unladen weight	12 kg	13.5 kg	14.5 kg	17 kg	18.5 kg	19.5 kg	21.5 kg	21.7 kg
HMV-No.	22.51.02.0053	22.51.02.0054	22.51.02.0054	22.51.02.0055	22.51.02.0051	22.51.02.0051	22.51.02.0052	22.51.02.0052

* Measurement with standard seat post / ** Measurement with T-seat post



Further information on
accessories and harnesses
can be found in the price list.

