

Stand up for yourself.

GMFCS
Level III-V

kg
max. 80



75–160 cm

(-)15° – 30°



I am your personal standing frame.

Using me, you can practice standing on your own two feet. Together, we will get your muscles fit, strengthen your bones and get your circulation going. Don't be afraid, I will be holding you all the time. Even if you lean forwards or backwards – my tilt angle can be adjusted. All my components can be adjusted in height, depth or angle precisely to suit your requirements.





Chest support pad with lateral guide
stabilises the upper body.

Pelvic frame
for an optimum upright position and easy transfer.

Pelvic positioner
mediates a feeling of security.

Flexible knee positioners (FKP)
adjustable in three dimensions, with adjustable contact bands.

Split footplate
can be adapted to individual requirements.

Castors with total locking device
prevents unwanted rolling away.

Headrest
provides additional support.

Therapy tray
suitable for therapeutic and play purposes.

Hand crank
permits the user to be brought slowly into upright position and to stretch their knees and pelvic area.

Central column
can be adjusted continuously from 15° reclined position to a forward tilt of up to 30°.

Frame geometry
guarantees that you can play together with other children at eye level.



You can order my wooden tray with or without adjustable arm supports, a bowl or a tray cover for therapy or play purposes.

Frame colour
Size 1

 Melon yellow

Frame colour
Size 2

 Dormant Red

Frame colour
Size 3

 Lazurite blue metallic

The pelvic frame helps you into an optimum upright position. Using a crank, your parents or carers can adjust the frame so that you can slowly stretch out your legs.



Flexible knee positioners (FKP) are individually three-dimensionally adjustable. Length-adjustable contact straps adapt optimally to the shape of the knee.

data	Size 1	Size 2	Size 3
Body size	75 – 110 cm	100 – 140 cm	120 cm – 160 cm
Tray height from footplate	48,5 – 80,5 cm	71 – 95 cm	80 – 112 cm
Pelvis depth	12 – 22 cm	16 – 28 cm	20 – 35 cm
Depth adjustment of the pelvis frame (via hand crank)	6 cm		
Total length	84 cm	90 cm	100 cm
Total width	57 cm	65 cm	71 cm
Abduction split footrests	0° – 27°		
Pelvis height	45 – 60 cm	57 – 80 cm	65 – 90 cm
Sternum height	55 – 85 cm	75 – 110 cm	95 – 125 cm
Max. load	40 kg	60 kg	80 kg
Weight	27 kg	29 kg	32 kg
HMV-No.	28.29.01.1059		



You can find more information on accessories in the price list.

