

The clever way to wash.

GMFCS
Level III-V



I am your versatile bathing and shower cradle.

I will adapt completely to your requirements, whether in the bathroom, in the sea or on the beach. My cover is washable to 30 °C, is particularly skin-friendly and leaves no pressure marks on your body. Because you can sit very low in the bathtub with me, I can save a lot of water. Aren't I practical? What's also practical is that my back angle can be adjusted with one hand during use - leaving the other hand free to support you!



Thanks to the adjustable headrest at
no water will run into your eyes.



Head section

can be easily folded.

Cover

can be removed for cleaning purposes.

Chest strap

provides reliable retention.

Optional pelvic mould depression

stabilises and positions the pelvis, with adjustable seat support for ischiatic bone.

Head guide pad

can be adjusted in height and width.

Backrest

can be adjusted from 90° to 180°.

Knee angle adjustment

provides a comfortable sitting position.

Lower leg support

can be adjusted between 0° and 90°.

Sub-frame

for the shower.

Pelvic mould depression

optionally, a cover is available with an integrated pelvic mould depression.

Back angle

can be adjusted even during use with one hand.

Size 4 has two stainless steel joints.





data	Size 1	Size 2	Size 3	Size 4
Seat width	30 cm	35 cm		40 cm
Seat depth	30 cm	34 cm	38cm	44cm
Back height	45 cm	61 cm	71 cm	82 cm
Standing space length	69 cm	74 cm	88 cm	93 cm
Lower leg length	19cm	20cm	28 cm	
Back angle	90° - 180°			
Knee angle	0° - 90°			
Headrest angle	0° - 90°			
Headrest angle	0° / 5cm			
Max. load	35 kg	45 kg	55 kg	75 kg
Weight	4 kg	4.1 kg	4.3 kg	4.8 kg
HMV-No.	04.40.04.0006			

* in combination with the cover with basin recess



You can find more information on accessories & harnesses
in the Price List.

