

GMFCS
Level III-V

I am your shower and toilet chair, I'm a real water fan!

With my many different sub-frames, I can support you in everything you do in the bathroom – whether during transfers, when used as a toilet chair, or in the shower. My back angle can be adjusted precisely, so you can lean forwards or backwards as needed. And because comfort and safety matter just as much as flexibility, my supportive seat unit, padding options and intuitive adjustments ensure you always feel secure and well-positioned. No matter the room layout or routine, I'm here to give you the ideal support you need.



Blue buttons

show the places where you can adjust the height, depth or angle – without the need for tools.

Adjustable seat angle at any height.





mounting bracket

This means that the seat unit can be easily installed on any toilet.

The bar can hold you permanently without restricting use.



Different sub-frames

for each spatial environment.



Mobile sub-frame



Fixed sub-frame



Bath sub-frame



Transfer and care sub-frame

Back angle

Can be set from $(-)\ 6^{\circ}$ to 18° .

Lateral pad

provides increased stability.

Multifunctional seating unit

can be adjusted in seat depth.

Seat pad

with two opening variations. Also available in an enclosed version.

Sub-frame

can be adjusted in height and angle.

Castors

facilitate uncomplicated transfer.

Headrest

can be adjusted in height, depth and angle.

Removable pads

can be cleaned easily.

Armrest

provides the necessary support.

Thigh guide

provides the necessary support.

Abduction wedge

guides the legs.

Footstool

ly adjustable and can be folded up for easy transfer.

Removable table

with padded stomach contact surface. This provides you with arm support and stimulation when going to the toilet.

Travel base hts.

Just take me with you in the bag that comes with the.

Toilet equipment

You can use me as a toilet chair independently of a toilet.

Frame colour  White



data		Size 1	Size 2	Size 3
Width	Mobile sub-frames	53 cm		57 cm
	Stationary sub-frames	51 cm		53 cm
	Transfer and care sub-frame	67 cm		71 cm
Seat height	Mobile sub-frames	39 - 55 cm		44 - 58 cm
	Stationary sub-frames	33 - 43 cm	36 - 44 cm	38 - 48 cm
	Transfer and care sub-frame	51 - 78 cm		
Seat width		20* / 25* / 30 cm	25* / 30* / 36 cm	33* / 38* / 43 cm
Seat depth		20 / 23 / 25 / 28 cm	29 / 32 / 34 / 37 cm	37 / 40 / 42 / 44 cm
Armrest height		13.5 cm	17 cm	20 cm
Back width		29 cm	34 cm	41 cm
Lower leg length, toilet seat		22 - 34 cm	27 - 39 cm	34 - 44 cm
Back height		42 cm	50 cm	62 cm
max. toilet height over which the product can be positioned	Mobile sub-frames	52 cm		52.5 cm
	Stationary sub-frames	38 cm		42 cm
Length	Mobile sub-frames	74 cm		81 cm
	Stationary sub-frames	71 cm		74 cm
	Transfer and care sub-frame	84 cm		87 cm
Front edge of hygiene opening to rear edge of mobile sub-frame		48 cm		
Back angle		(-) 6°, 6°, 18°		
Seat angle		(-) 5° - 25°		
		(-) 15° - 15°		
Max. load		34 kg	68 kg	113 kg
HMV-No.		33.40.01.4005 (use as a toilet attachment) / 33.40.04.1008 (use as a toilet chair)		

* adjustable using optional thigh guide



You can find more information on accessories & harnesses in the price list.

